

Product Information

Luciferase from *Photinus pyralis* (firefly)

Recombinant, expressed in *E. coli*, lyophilized powder, $\geq 10 \times 10^{10}$ units/mg protein

SRE0045

Product Description

CAS Registry Number: 61970-00-1

Enzyme Commission (EC) Number: 1.13.12.7

Synonyms: Luciferin 4-monooxygenase, Firefly Luciferase

Storage Temperature: $-20\text{ }^{\circ}\text{C}$

Firefly luciferase is a 62 kDa protein that catalyzes the production of light. The enzyme requires ATP, molecular oxygen, and luciferin, a heterocyclic compound, to generate light in a two-step process.¹ The light-producing reaction is initiated by luciferin activation (adenylation of its carboxylate group) and proceeds in the presence of molecular oxygen to yield a photon of yellow-green light.^{1,2}

Firefly luciferase is used extensively in molecular and cell biology, in particular for the efficient detection and quantitation of ATP and as a reporter for genetic function.^{3,4}

This product is a recombinant luciferase from *Photinus pyralis* (American firefly) produced from the *luc* gene expressed in *E. coli*. Several references⁵⁻¹² and dissertations^{13,14} have cited use of this product in their research.

Product

This product is lyophilized from a buffered solution that contains HEPES (pH 7.5), NaCl, MgCl₂, EDTA, DTT, and a carbohydrate stabilizer.

Specific Activity: $\geq 10 \times 10^{10}$ light units/mg protein

Unit definition: One luciferase enzyme unit will produce one Relative Light Unit (RLU) at 20-25 °C over a 10-second period, measured in a 100 µL assay mixture that contains 40 pmole ATP and 15 nmole luciferin in Tris-glycine buffer (pH 7.6), using a GloMax™ 20/20 Luminometer.

Precautions and Disclaimer

This product is for R&D use only. Not for drug, household, or other uses. Please consult the Safety Data Sheet for information regarding hazards and safe handling practices.

Storage/Stability

Store the product at $-20\text{ }^{\circ}\text{C}$.

Preparation Instructions

The enzyme can be prepared at a concentration of up to 5 mg protein/mL. To obtain maximal solubility, it is important to reconstitute the enzyme at a high salt concentration, such as 1 M Tris buffer with any counterion at pH 7-8. **Do not vortex. Avoid agitation.**

After reconstitution, the enzyme solutions can be kept at 4-8 °C for up to 2 days or frozen in working aliquots at $-20\text{ }^{\circ}\text{C}$ for at least one month. Repeated freezing and thawing is not recommended. Do not store in a frost-free freezer.

References

1. DeWet, J. R. et al., *Mol. Cell. Biol.*, **7**(2), 725-737 (1987).
2. Stanley, P. E., *J. Biolumin. Chemilumin.*, **4**(1), 375-380 (1989).
3. Kricka, L. J., *Anal. Biochem.*, **175**(1), 14-21 (1988).
4. Chappelle, E. W. et al., *Meth. Enzymol.*, **57**, 65-72 (1978).
5. Chen, D. et al., *Nat. Commun.*, **8**(1), 2265 (2017).
6. Tariq, A. et al., *Cell Rep.*, **28**(8), 2080-2095 (2019).

7. Adelöf, J. *et al.*, *Aging Cell*, **20(4)**, e13336 (2021).
8. Harada, L. *et al.*, *Biosensors (Basel)*, **11(4)**, 124 (2021).
9. Yevtodiyenko, A. *et al.*, *Nat. Commun.*, **12**, 2680 (2021)
10. Bhadra, A. K. *et al.*, *Nat. Commun.*, **13(1)**, 4570 (2022).
11. Vyas, P. *et al.*, *J. Am. Chem. Soc.*, **145(15)**, 8344-8354 (2023).
12. Greiser, M. *et al.*, *eLife*, **12**, e84204 (2023).
13. Godinat, Aurélien, "The Split Luciferin Reaction: From Bioorthogonal Chemistry to Bioluminescence Imaging". École Polytechnique Fédérale de Lausanne, Ph.D. dissertation, pp. 38, 61 (2017).
14. Coustets, Mathilde, "La lectine de *Xerocomellus Chrysenteron*, un nano-objet théranostique pour l'imagerie et le traitement des cancers épithéliaux: preuve de concept appliquée aux carcinomes péritonéaux d'origine ovarienne" ("*Xerocomellus Chrysenteron* lectin, a theranostic nano-object for imaging and treatment of epithelial cancers: proof of concept applied to peritoneal carcinomatosis of ovarian origin"). Université de Toulouse, Ph.D. dissertation, p. 216 (2020).

Notice

We provide information and advice to our customers on application technologies and regulatory matters to the best of our knowledge and ability, but without obligation or liability. Existing laws and regulations are to be observed in all cases by our customers. This also applies in respect to any rights of third parties. Our information and advice do not relieve our customers of their own responsibility for checking the suitability of our products for the envisaged purpose.

The information in this document is subject to change without notice and should not be construed as a commitment by the manufacturing or selling entity, or an affiliate. We assume no responsibility for any errors that may appear in this document.

Technical Assistance

Visit the tech service page at SigmaAldrich.com/techservice.

Terms and Conditions of Use

Warranty, use restrictions, and other conditions of sale may be found at SigmaAldrich.com/terms.

Contact Information

For the location of the office nearest you, go to SigmaAldrich.com/offices.

The life science business of Merck KGaA, Darmstadt, Germany operates as MilliporeSigma in the U.S. and Canada.

MilliporeSigma, and Sigma-Aldrich are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.

© 2016-2023 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved.

SRE0045dat Rev 06/23 DS,DT,GCY,MAM

**MILLIPORE
SIGMA**