


**RABBIT ANTI-NEUROFASCIN 155
AFFINITY PURIFIED
POLYCLONAL ANTIBODY**

CATALOG NUMBER:	AB15189
LOT NUMBER:	
QUANTITY:	100 µg
CONCENTRATION:	1 mg/mL
SPECIFICITY:	Neurofascin 155. By Western blot the antibody recognizes a band at ~155 kDa in mouse brain lysate. An additional band at ~80 kDa may be seen depending on sample and antibody concentration used.
IMMUNOGEN:	Synthetic peptide.
APPLICATIONS:	Western blot: 0.5-1.0 µg/mL on extracts from mouse brain lysate.



Optimal working dilutions must be determined by the end user.

SPECIES REACTIVITY: Mouse. Other species have not been tested.

FORMAT: Affinity purified immunoglobulin.

PRESENTATION: Liquid in PBS with 0.1% sodium azide.

STORAGE/HANDLING: Maintain at 2-8°C in undiluted aliquots for up to 6 months after date of receipt.



Important Note: *During shipment, small volumes of product will occasionally become entrapped in the seal of the product vial. For products with volumes of 200 μ L or less, we recommend gently tapping the vial on a hard surface or briefly centrifuging the vial in a tabletop centrifuge to dislodge any liquid in the container's cap.*

FOR RESEARCH USE ONLY; NOT FOR USE IN DIAGNOSTIC
PROCEDURES. NOT FOR HUMAN OR ANIMAL CONSUMPTION

Unless otherwise stated in our catalog or other company documentation accompanying the product(s), our products are intended for research use only and are not to be used for any other purpose, which includes but is not limited to, unauthorized commercial uses, in vitro diagnostic uses, ex vivo or in vivo therapeutic uses or any type of consumption or application to humans or animals.

©2002 - 2011: Millipore Corporation. All rights reserved. No part of these works may be reproduced in any form without permission in writing.

28820 Single Oak Drive • Temecula, CA 92590
Technical Support: T: 1-800-MILLIPORE (1-800-645-5476) • F: 1-800-437-7502
www.millipore.com