

Product Information

Monoclonal Anti-Rat IgG2a

Clone R2A-2

Mouse Ascites Fluid

Product Number **R 0761**

Product Description

Monoclonal Anti-Rat IgG2a (mouse IgG1 isotype) is derived from the R2A-2 hybridoma produced by the fusion of mouse myeloma cells and splenocytes from BALB/c mice immunized with purified rat IgG. The isotype is determined using the Sigma ImmunoType Kit (Product Code ISO-1) and by a double diffusion immunoassay using Mouse Monoclonal Antibody Isotyping Reagents (Product Code ISO-2).

Monoclonal Anti-Rat IgG2a recognizes an epitope located on the heavy chain of rat IgG2a. The antibody detects rat IgG2a derived from normal serum or myeloma proteins, but not the other rat immunoglobulins. It localizes the denatured-reduced heavy chain molecule of rat IgG2a (γ 2a), when applied in immunoblotting. Weak cross-reaction is observed with guinea pig immunoglobulins but not with IgG and serum preparations of bovine, cat, chicken, dog, goat, horse, human, mouse, pig, rabbit, or sheep when tested by indirect ELISA and dot blot techniques. The antibody is also applicable as a secondary antibody in immunohistochemistry of human tissue where it does not react against the tissue itself.

Rat immunoglobulins are subdivided into five classes; IgM, IgG, IgA, IgE, IgD, and four IgG subclasses; IgG1, IgG2a, IgG2b, and IgG2c, on the basis of the structural, biological, physiochemical and electrophoretic properties of their heavy chains.¹ The rat has been extensively used as research model in pharmacology, oncology, and the study of the immunology of aging. Rat polyclonal and monoclonal antibodies have come into widespread use as primary antibodies.^{2,3} Secondary antibodies to rat immunoglobulin subtypes may be particularly valuable in double labeling experiments and for isotyping and immunoaffinity purification of rat-derived antibodies. Anti-rat antibodies are commonly produced by xenogeneic immunization of rabbits, goats, or sheep, resulting in antibodies that

cross-react with other immunoglobulin subclasses of rat and of other species, unless extensively adsorbed. Monoclonal anti-rat immunoglobulins, which are devoid of any binding capacity to human and many other species, can serve as an essential tool in many applications, especially when used as a secondary reagent in immunohistochemistry.

Monoclonal Anti-Rat IgG2a may be used for the localization of rat IgG2a using various immunochemical assays such as ELISA, immunoblot, dot blot or immunocytochemistry.

Reagents

The product is provided as ascites fluid with 0.1% sodium azide as a preservative.

Precautions and Disclaimer

Due to the sodium azide content a material safety data sheet (MSDS) for this product has been sent to the attention of the safety officer of your institution. Consult the MSDS for information regarding hazardous and safe handling practices.

Storage/Stability

For continuous use, store at 0-5 °C. For extended storage freeze in working aliquots. Repeated freezing and thawing is not recommended. Storage in "frost-free" freezers is not recommended. If slight turbidity occurs upon prolonged storage, clarify the solution by centrifugation before use.

Product Profile

The antibody titer (1:5,000) was determined by indirect ELISA, using 2 μ g/mL freshly prepared rat myeloma protein for coating.

Note: Second antibodies against mouse immunoglobulins may cross-react with the rat proteins coated on the microtiter plate, unless properly adsorbed with rat immunoglobulins.

In order to obtain best results in different techniques and preparations, it is recommended that each individual user determine their optimum working dilutions by titration assay.

References

1. Bazin, H., et al., *Eur. J. Immunol.*, **4**, 44 (1974).
2. Springer, T. A., et al., *Hybridoma*, **1**, 257 (1982).
3. Bazin, H., ed., *Rat Hybridomas and Rat Monoclonal Antibodies*, (CRC Press, Boca Raton, Florida, 1990).

KAA 02/06

Sigma brand products are sold through Sigma-Aldrich, Inc.

Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.