

**MOUSE ANTI-HUMAN TENASCIN
MONOCLONAL ANTIBODY**

CATALOG NUMBER: MAB1911

LOT NUMBER:

QUANTITY: 100 µg

CONCENTRATION: 1 mg/mL

SPECIFICITY: Binds to a recombinant fragment of human tenascin containing the fibronectin type III repeats, numbers 1-3.

CLONE NAME: 9F8

IMMUNOGEN: Purified human tenascin

ISOTYPE: IgG_{1κ}

APPLICATIONS: Immunohistochemistry at 2 µg/mL on tissues fixed with cold acetone
Immunocytochemistry at 2 µg/mL on U251 cells
Immunoblotting (Western) at 2 µg/mL to detect whole tenascin under reducing conditions
ELISA: 2×10^{-10} M endpoint dilution
Affinity constant $K_d = 6.37 \times 10^{-9}$
Final working dilutions must be determined by end user.

FORMAT: Immunoglobulin purified by Protein A chromatography.

PRESENTATION: Liquid in 0.02 M phosphate, 0.25 M Sodium Chloride pH 7.6, with 0.1% Sodium Azide.

STORAGE/HANDLING: Store at 2 and 8°C for up to 6 months.

Important Note: *During shipment, small volumes of antibody will occasionally become entrapped in the seal of the product vial. For antibodies with volumes of 200 µl or less, we recommend gently tapping the vial on a hard surface or briefly centrifuging the vial in a tabletop centrifuge to dislodge any liquid in the container's cap.*

For research use only; not for use as a diagnostic.

© 2003-2005: CHEMICON International, Inc. - By CHEMICON International, Inc. All rights reserved. No part of these works may be reproduced in any form without permissions in writing.

USA & Canada • Phone: +1(800) 437-7500 • Fax: +1 (951) 676-9209 • Europe +44 (0) 23 8026 2233
Australia +61 3 9839 2000 • Germany +49-6192-207300 • ISO Registered worldwide
www.chemicon.com • custserv@chemicon.com • techserv@chemicon.com