

## Product Information

# 15895 *Listeria* mono Enrichment Supplement II

For selective differentiation of *Listeria monocytogenes* from other *Listeria* species

### Composition

(per vial, sufficient for 500 ml medium)

L-phosphatidylinositol 0.5g

Distilled water 15 ml

### Directions

Thaw the contents of 1 vial L. mono Enrichment Supplement II at room temperature. Aseptically add the sterile contents to 470 ml of L. mono confirmatory Agar Base (Cat. No. 92302) along with sterile rehydrated contents of L. mono Selective Supplement I (Cat. No. 92301) and L. mono Selective Supplement II (Cat. No. 91603). Mix well and pour into sterile petri plates.