

17194 Park and Sanders Selective Supplement II

An antibiotic Supplement recommended for the selective isolation and enumeration of thermotolerant *Campylobacter* species.

Composition
(per vial, sufficient for 1000 ml medium)

Cefoperazone	32 mg
Cycloheximide	100 mg

Directions

Rehydrate the contents of 1 vial aseptically with 5 ml of 50% ethanol. Mix well and aseptically add it to 1000 ml of sterile *Campylobacter* Selective Agar, Base (Prod No 21378) containing 50- 70 ml sterile defibrinated blood and Park and Sanders Selective Supplement I (Prod No 17191). Mix well and dispense.

Precautions and Disclaimer

This product is for R&D use only, not for drug, household, or other uses. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

