



**MILLIPLEx[®] MAP Skin Map Magnetic Panel
SKINMAG-50K**

Standard Catalog # L-8050

STANDARD LOT#	ANALYTE	STANDARD VALUE	UNITS
850-14K	Keratin 6	2440	ng/mL
850-14K	Keratin-1,10	37958	ng/mL

Control Catalog # L-6050

<u>Analyte</u>	<u>QC Level</u>	<u>Lot #</u>	<u>Expected Range</u>	<u>Units</u>
Cortisol	Control I	SM-114	6.0 – 13	ng/mL
	Control II	SM-214	37 – 76	ng/mL
HSA	Control I	SM-114	27 – 56	ng/mL
	Control II	SM-214	161 – 359	ng/mL
Keratin-6	Control I	SM-114	38 – 90	ng/mL
	Control II	SM-214	213 – 442	ng/mL
Keratin-1,10	Control I	SM-114	545 – 1132	ng/mL
	Control II	SM-214	2614 – 5430	ng/mL
Involucrin	Control I	SM-114	0.51 – 1.1	ng/mL
	Control II	SM-214	3.2 – 6.7	ng/mL

L-6050: SM-114.214 (07 August 2019)

EMD MILLIPORE • 1-800-645-5476 • 1-800-645-5439 (FAX) • www.emdmillipore.com