

Product Information

SILu™Prot AKT1, RAC-alpha serine/threonine-protein kinase, human recombinant, expressed in HEK 293 cells
SIL MS Protein Standard, ¹³C and ¹⁵N -labeled

Catalog Number **MSST0049**
Storage Temperature **-20 °C**

Synonyms: Protein kinase B (PKB), Protein kinase B alpha (PKB alpha), Proto-oncogene c-Akt, RAC-PK-alpha

Product Description

SILu™Prot AKT1 is a recombinant, stable isotope-labeled human AKT1 which incorporates [¹³C₆, ¹⁵N₄]-Arginine and [¹³C₆, ¹⁵N₂]-Lysine. Expressed in human 293 cells, it is designed to be used as an internal standard for bioanalysis of AKT1 in mass-spectrometry. SILu™Prot AKT1 is a protein of 501 amino acids (including N-terminal polyhistidine and FLAG® tags), with a calculated molecular mass of 58.8 kDa.

AKT1 is a serine/threonine kinase that is a member of the AKT family. AKT1, like the other AKT proteins, is activated in cells in response to diverse stimuli such as hormones, growth factors, and extracellular matrix components.¹ Once activated by phosphorylation at Ser⁴⁷³ and Thr³⁰⁸, AKT1 promotes proliferation, cell survival, motility, and angiogenesis processes, interfering with the apoptotic functions of the cell.²

Overexpressed phospho-AKT is frequently observed in human lung, gastric, hepatocellular, pancreatic, renal, prostate, and endometrial cancer as well as multiple myeloma.³⁻⁴ The aggressiveness of several types of solid tumors and hematologic malignancies is linked to the deregulation of AKT and its upstream signaling partners.⁵ Members of the AKT pathway are therefore potential targets for novel anti-cancer therapeutics.⁵

Each vial contains 10 µg of SILu™Prot AKT1 standard in a solution of phosphate buffered saline with 1mM EDTA and 25% glycerol. Vial content was determined by the Bradford method using BSA as a calibrator.

Purity: ≥95% (SDS-PAGE)

Heavy amino acids incorporation efficiency: ≥98% (MS)

UniProt: P31749

Sequence Information:

The N-terminal polyhistidine and FLAG tags are italicized.

MDYKDDDDKGGHHHHHHHGGQMSDVAIVKEGWLH
KRGEYIKTWRPRYFLLKNDGTFIGYKERPDVDQRE
APLNNFSVAQCQLMKTERPRPNTFIIRCLQWTTVIERT
FHVETPEEREETTAIQTVADGLKKQEEEEEMDFRSG
SPSDNSGAEEMEVS LAKPKHRVTMNEFEYLKLLGKG
TFGKVLVKEKATGRYYAMKILKKEVIVAKDEVAHTLT
ENRVLQNSRHPFLTALKYSFQTHDRLCFVMEYANGG
ELFFHLSRERVFSEDRARFYGAEIVSALDYHSEKNV
VYRDLKLENMLDKDGHKIDDFGLCKEGIKDGATMKT
FCGTPEYLAPEVLEDNDYGRAVDWWGLGVVYEMM
CGRLPFYNQDHEKLFELILMEEIRFPRTLGPPEAKSLLS
GLLKKDPKQRLGGGSEDAKEIMQHRFFAGIVWQHVV
EKKLSPPFKPQVTSETDTRYFDEEFTAQMITITPPDQD
DSMECVDSERRPHFPQFSYSASGTA

Precautions and Disclaimer

This product is for R&D use only, not for drug, household, or other uses. Please consult the Safety Data Sheet for information regarding hazards and safe handling practices.

Storage/Stability

Store the product at **-20 °C**. The product retains its concentration for at least 2 years as supplied. After initial thawing it is recommended to store the protein in working aliquots at **-20 °C**.

References

1. Alessi, D.R., and Cohen, P., Mechanism of activation and function of protein kinase B. *Curr. Opin. Genet. Dev.*, **8(1)**, 55-62 (1998).
2. Coffey, P.G. et al., Protein kinase B (c-Akt): a multifunctional mediator of phosphatidylinositol 3-kinase activation. *Biochem. J.*, **335(1)**, 1-13 (1998).
3. Altomare, D.A., and Testa, J.R., Perturbations of the AKT signaling pathway in human cancer. *Oncogene*, **24(50)**, 7455-7464 (2005).
4. Cicens, J., The potential role of Akt phosphorylation in human cancers. *Int. J. Biol. Markers*, **23(1)**, 1-9 (2008).
5. Garcia-Echeverria, C., and Sellers, W.R., Drug discovery approaches targeting the PI3K/Akt pathway in cancer. *Oncogene*, **27(41)**, 5511-5526 (2008).

Legal Information

Sold under license from DuPont, U.S. Patent No. 7,396,688.

This product is licensed under U.S. Patent No. 7,396,688 and foreign counterparts from E. I. du Pont de Nemours and Company. The purchase of this product conveys to the buyer the nontransferable right to use the purchased amount of the product for research and development only, including services for a third party for consideration. The buyer cannot sell or otherwise transfer this product, its components or materials made using this product or its components to a third party. Information about licenses for excluded uses is available from: E. I. du Pont de Nemours and Company; Attn: Associate Director, Commercial Development; DuPont Experimental Station E268; 200 Powdermill Rd.; Wilmington, DE 19803; 1-877-881-9787 (voice), 1-302-695-1437 (fax), licensing@dupont.com.

SILu is a trademark of Sigma-Aldrich Co. LLC.
FLAG is a registered trademark of Sigma-Aldrich Co. LLC.

JK,MAM 05/17-1