

3050 Spruce Street Saint Louis, Missouri 63103 USA Telephone (800) 325-5832 (314) 771-5765 Fax (314) 286-7828 email: techserv@sial.com sigma-aldrich.com

# **ProductInformation**

# WORT BROTH

Product Number W2011

## **Product Description**

Wort Broth is used for the cultivation of yeasts, especially in syrups and butter.<sup>1</sup> Yeasts grow well in culture media containing dextrose or maltose in an acidic environment. The peptic digest of animal tissue and malt extract provide nitrogenous compounds for the growth of yeast. Dextrin and maltose are the fermentable carbohydrates. The highly acidic pH inhibits many bacteria.

Components	
ltem	<u>g/L</u>
Malt Extract	15.00
Peptic Digest of Animal Tissue	0.78
Maltose	12.75
Dextrin	2.75
Dipotassium Phosphate	1.00
Ammonium Chloride	1.00

Final pH (at 25 °C) 4.8 ± 0.2

## **Precautions and Disclaimer**

For laboratory use only. Not for drug, household or other uses.

#### **Preparation Instructions**

Suspend 33.28 grams of Wort Broth in 1000 mls of distilled water containing 2.35 grams of (G6279) glycerol. Boil to dissolve the medium completely. Sterilize by autoclaving at 15 lbs. pressure (121 °C) for 15 minutes.

## Storage

Store the dehydrated medium at 24  $^\circ C$  and the prepared medium at 2-8  $^\circ C.$ 

### **Product Profile**

Appearance	Yellow colored, homogeneous, free flowing powder.
Color and Clarity.	Yellow colored, opalescent solution forms in tubes, and may contain a flocculent precipitate
Cultural Response	Cultural characteristics are observed after 40-48 hours at 30 °C.

Organisms	(ATCC)	Growth
Aspergillus niger	(1604)	luxuriant
Saccharomyces cerevisiae	(9763)	luxuriant
Saccharomyces uvarum	(9080)	luxuriant
Candida albicans	(10231)	luxuriant

#### References

- 1. Parfitt (1933), J. Dairy Sci., 19: 141.
- MacFaddin, J., (1985). Media for Isolation-Cultivation-Maintenance of Medical Bacteria. Vol.1. Williams and Wilkins. Baltimore, Maryland.
- 3. American Type Culture Collection, Manassas, Va., U.S.A..

Sigma brand products are sold through Sigma-Aldrich, Inc.

Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.