



**SIGMA QUALITY CONTROL TEST
PROCEDURE****Enzymatic Assay of PYRUVATE CARBOXYLASE
(EC 6.4.1.1)****PRINCIPLE:**

Abbreviations used:

ATP = Adenosine 5'-Triphosphate

PC = Pyruvate Carboxylase

AcCOA = Acetyl Coenzyme A

ADP = Adenosine 5'-Diphosphate

P_i = Inorganic Phosphate

β-NADH = β-Nicotinamide Adenine Dinucleotide, Reduced Form

MDH = Malic Dehydrogenase

β-NAD = β-Nicotinamide Adenine Dinucleotide, Oxidized Form

CONDITIONS: T = 30°C, pH = 7.8, A_{340nm}, Light path = 1 cm

METHOD: Continuous Spectrophotometric Rate Determination

REAGENTS:

- A. 135 mM Triethanolamine Buffer with 7 mM Magnesium Sulfate, 9 mM Pyruvic Acid, and 0.15% (w/v) Bovine Serum Albumin, pH 8.0 at 30°C (Substrate)
(Prepare 100 ml in deionized water using Triethanolamine Hydrochloride, Sigma Prod. No. T-1502, Magnesium Sulfate, Anhydrous, Sigma Prod. No. M-7506, Pyruvic Acid, Sodium Salt, Sigma Prod. No. P-2256, and Albumin, Bovine, Sigma Prod. No. A-6003. Adjust to pH 8.0 at 30°C with 1 M HCl.)
- B. 0.3 mM Acetyl Coenzyme A Solution
(Prepare 10 ml in deionized water using Acetyl Coenzyme A, Sodium Salt, Sigma Prod. No. A-2056. **PREPARE FRESH.**)

**Enzymatic Assay of PYRUVATE CARBOXYLASE
(EC 6.4.1.1)**

REAGENTS: (continued)

- C. Malic Dehydrogenase Enzyme Solution (AcCoA/MDH)
(Immediately before use, add 150 units of Malic Dehydrogenase, Sigma Prod. No. M-9004 to 5 ml of Reagent B. Bring the solution to a total volume of 10.0 ml with deionized water.)
- D. 100 mM Triethanolamine Buffer with 30 mM Adenosine 5'-Triphosphate and 450 mM Potassium Bicarbonate, pH 8.0 at 30°C (ATP/KHCO₃)
(Prepare 10 ml in deionized water using Triethanolamine Hydrochloride, Sigma Prod. No. T-1502, Adenosine 5'-Triphosphate, Disodium Salt, Sigma Prod. No. A-5394, and Potassium Bicarbonate, Sigma Prod. No. P-9144. Adjust to pH 8.0 at 30°C with 1 M KOH.)
- E. 2.6 mM β-Nicotinamide Adenine Dinucleotide, Reduced Form (β-NADH)
(Prepare 10 ml in deionized water using β-Nicotinamide Adenine Dinucleotide, Reduced Form, Disodium Salt, Sigma Prod. No. N-8129.)
- F. 50 mM Tris HCl Buffer with 50% (v/v) Glycerol, 2 mM Magnesium Acetate, and 1 mM Ethylenediaminetetraacetic Acid, pH 7.4 at 30°C (Enz Dil)
(Prepare 50 ml in deionized water using Trizma Base, Sigma Prod. No. T-1503, Glycerol, Sigma Prod. No. G-7893, Magnesium Acetate, Tetrahydrate, Sigma Prod. No. M-0631, and Ethylenediaminetetraacetic Acid, Tetrasodium Salt, Hydrate, Sigma Stock No. ED4SS. Adjust to pH 7.4 at 30°C with 1 M HCl.)
- G. Pyruvate Carboxylase Enzyme Solution (PC)
(Immediately before use, prepare a solution containing approximately 30 - 90 units/ml of Pyruvate Carboxylase in cold Reagent F.)

PROCEDURE:

Prepare a reaction cocktail by pipetting (in milliliters) the following reagents into a suitable container:

Reagent A (Substrate)	20.00
Reagent C (AcCoA/MDH)	5.00
Reagent E (β-NADH)	2.50

Mix by swirling. Adjust to pH 7.8 at 30°C if necessary, with either 1 M HCl or 1 M KOH.

**Enzymatic Assay of PYRUVATE CARBOXYLASE
(EC 6.4.1.1)**

PROCEDURE:

Pipette (in milliliters) the following reagents into a suitable container:

	<u>Test</u>	<u>Blank</u>
Reaction Cocktail	2.90	2.90
Reagent G (PC)	0.005	-----
Reagent F (Enz Dil)	-----	0.005

Mix by inversion and equilibrate to 30°C. Monitor the $A_{340\text{nm}}$ until constant using a suitably thermostatted spectrophotometer. Record the $\Delta A_{340\text{nm}}/\text{minute}$ of the Test.¹ Then add:

Reagent D (ATP/KHCO ₃)	0.10	0.10
------------------------------------	------	------

Immediately mix by inversion and record the decrease in $A_{340\text{nm}}$ for approximately 5 minutes. Obtain the $\Delta A_{340\text{nm}}/\text{minute}$ using the maximum linear rate for both the Test and Blank.

CALCULATIONS:

$$\text{Units/ml enzyme} = \frac{(\Delta A_{340\text{nm}}/\text{min Test} - \Delta A_{340\text{nm}}/\text{min Blank})(3.005)(\text{df})}{(6.22)(0.005)}$$

3.005 = Total volume (in milliliters) of the assay

df = Dilution factor

6.22 = Millimolar extinction coefficient of β -NADH at 340 nm

0.005 = Volume (in milliliter) of pyruvate carboxylase used in the assay

$$\text{Units/mg protein} = \frac{\text{units/ml enzyme}}{\text{mg protein/ml enzyme}}$$

UNIT DEFINITION:

One unit will convert 1.0 μmole of pyruvate and CO₂ to oxalacetate per minute at pH 7.8 at 30°C.

**Enzymatic Assay of PYRUVATE CARBOXYLASE
(EC 6.4.1.1)**

FINAL CONCENTRATIONS:

In a 3.005 ml reaction mix, the final concentrations are 134 mM triethanolamine, 5 mM magnesium sulfate, 7 mM pyruvic acid, 0.12% (w/v) bovine serum albumin, 0.23 mM β -nicotinamide adenine dinucleotide, reduced form, 0.05 mM acetyl coenzyme A, 2.63 units malic dehydrogenase, 1 mM adenosine 5'-triphosphate, 15 mM potassium bicarbonate, 0.05% (v/v) glycerol, 0.002 mM magnesium acetate, 0.001 mM ethylenediaminetetraacetic acid, 0.05 mM Tris, 0.15 - 0.45 unit pyruvate carboxylase.

REFERENCE:

Warren, G.B. and Tipton, K.F. (1974) *Biochemical Journal* **139**, 297-310

NOTES:

1. Lactic Dehydrogenase is the principle contaminant that may interfere with the assay for pyruvate carboxylase. If the $\Delta A_{340\text{nm}}/\text{min}$ is not zero, it must be subtracted from the $\Delta A_{340\text{nm}}/\text{min}$ for the Test after the addition of Reagent D (ATP/ KHCO_3).
2. Malic Dehydrogenase Unit Definition: One unit will convert 1.0 μmole of oxalacetate and β -NADH to L-malate and β -NAD per minute at pH 7.5 at 25°C.
3. This assay is based on the cited reference.
4. Where Sigma Product or Stock numbers are specified, equivalent reagents may be substituted.

Sigma warrants that the above procedure information is currently utilized at Sigma and that all Sigma-Aldrich, Inc. products conform to the information in this and other Sigma-Aldrich, Inc. publications. Purchaser must determine the suitability of the information and product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packaging slip.