

## Product Information

## Anti-c-Myc-Peroxidase Antibody

Mouse Monoclonal, Clone 9E10, Purified from Hybridoma Cell Culture

**SAB4200742**

### Product Description

Monoclonal Anti-c-Myc (mouse IgG1 isotype) is derived from the 9E10 hybridoma, produced by the fusion of mouse myeloma cells and splenocytes from a BALB/c mouse immunized with synthetic peptide corresponding to residues 408-439 of the human p62<sup>c-myc</sup> protein, conjugated to KLH<sup>1</sup>. This clone is referenced also as Myc1-9E10.<sup>1</sup> The isotype is determined by ELISA using Mouse Monoclonal Antibody Isotyping Reagents (ISO2). The antibody is purified from culture supernatant of hybridoma cells and is conjugated to horseradish peroxidase.

Monoclonal Anti-c-Myc specifically recognizes an epitope located within the sequence EQKLISEEDL (residues 410-419) of the product of the human oncogene c-Myc, known as the c-Myc tag.<sup>2</sup> The antibody recognizes the c-Myc tag sequence when it is expressed at either the amino or the carboxyl terminus of the fusion protein. The antibody is recommended to use in various immunological techniques, including immunoblotting and immunohistology<sup>3-4</sup>. Reaction of the antibody in immunoblotting is inhibited by the c-Myc peptide (M2435).<sup>5</sup>

c-Myc is a proto-oncogene encoding a transcription factor which is involved in the regulation of cell proliferation, growth and cell death. The human c-Myc proto-oncogene is the cellular homologue of the avian v-Myc gene found in several leukemogenic retroviruses.<sup>1,6</sup> Increased expression of the cellular oncogene c-Myc has been described in a variety of human tumors and may occur by several mechanisms, including gene amplification and chromosomal translocation.<sup>7</sup>

The c-Myc gene encodes a polypeptide with predicted molecular weight of 49 kDa but show aberrant electrophoretic mobility on polyacrylamide gel electrophoresis with an apparent molecular weight of around 62 kDa (p62<sup>c-myc</sup>).<sup>8</sup> p62<sup>c-myc</sup> is associated mainly with the cell nucleus, where it exerts its normal and oncogenic functions.<sup>1,7</sup>

c-Myc sequence is widely used as an epitope tag which provides an opportunity to create a stable fusion product that does not seem to interfere with the protein bioactivity or the biodistribution of the c-Myc tagged product.<sup>2,9-11</sup>

Monoclonal anti-c-Myc-Peroxidase antibody can serve a useful tool for study the prevalence of the protein product of the c-Myc oncogene in cells and tissues or for detection of c-Myc tagged proteins in bacterial and cellular lysates.<sup>11-10</sup>

### Reagent

Supplied as a lyophilized powder.

### Preparation Instructions

Reconstitute the content of the vial with 0.1/0.5 mL of distilled water to a final antibody concentration of ~ 2 mg/mL. After reconstitution, the solution contains 1% BSA, 2.5% trehalose, 0.05% MIT in 0.01 M sodium phosphate buffered saline.

### Precautions and Disclaimer

For R&D use only. Not for drug, household, or other uses. Please consult the Safety Data Sheet for information regarding hazards and safe handling practices.

## Storage/Stability

Store the lyophilized product at 2–8 °C. For extended storage after reconstitution, keep at –20 °C in working aliquots. Avoid repeated freeze-thaw cycles. For continuous use after reconstitution, keep at 2–8 °C for up to 1 month. Solutions at working dilution should be discarded if not used within 12 hours.

## Product Profile

### Immunoblotting

A working dilution of 1:250-1:500 is recommended using lysate of HEK-293T cells over expressing c-Myc fusion protein.

**Note:** In order to obtain best results in different techniques and preparations we recommend determining optimal working concentration by titration test.

## References

1. Evan, G. I., Lewis, G. K., Ramsay, G., & Bishop, J. M. (1985). Isolation of monoclonal antibodies specific for human c-myc proto-oncogene product. *Molecular and cellular biology*, 5(12), 3610-3616.
2. Campbell, A. M., Kessler, P. D., & Fambrough, D. M. (1992). The alternative carboxyl termini of avian cardiac and brain sarcoplasmic reticulum/endoplasmic reticulum Ca (2+)-ATPases are on opposite sides of the membrane. *Journal of Biological Chemistry*, 267(13), 9321-9325.
3. Baker, A. M., Van Noorden, S., Rodriguez-Justo, M., Cohen, P., Wright, N. A., & Lampert, I. A. (2016). Distribution of the c-MYC gene product in colorectal neoplasia. *Histopathology*, 69(2), 222-229.
4. Blancato, J., Singh, B., Liu, A., Liao, D. J., & Dickson, R. B. (2004). Correlation of amplification and overexpression of the c-myc oncogene in high-grade breast cancer: FISH, in situ hybridisation and immunohistochemical analyses. *British journal of cancer*, 90(8), 1612-1619.
5. Jarvik, J. W., & Telmer, C. A. (1998). Epitope tagging. *Annual review of genetics*, 32(1), 601-618.
6. Alitalo, K., Ramsay, G., Bishop, J. M., Pfeifer, S. O., Colby, W. W., & Levinson, A. D. (1983). Identification of nuclear proteins encoded by viral and cellular myc oncogenes. *Nature*, 306(5940), 274-277.
7. Pelengaris, S., & Khan, M. (2001). Oncogenic co-operation in beta-cell tumorigenesis. *Endocrine-Related Cancer*, 8(4), 307-314.
8. Robertson, D., Paterson, H. F., Adamson, P., Hall, A., & Monaghan, P. (1995). Ultrastructural localization of ras-related proteins using epitope-tagged plasmids. *Journal of Histochemistry & Cytochemistry*, 43(5), 471-480.
9. Hu, J., Xiao, X. D., Ogletree, D. F., & Salmeron, M. (1995). Imaging the condensation and evaporation of molecularly thin films of water with nanometer resolution. *Science*, 268(5208), 267-269.
10. Kobayashi, R., Kubota, T., & Hidaka, H. (1994). Purification, characterization, and partial sequence-analysis of a new 25-kDa actin-binding protein from bovine aorta: a SM22 homolog. *Biochemical and biophysical research communications*, 198(3), 1275-1280.
11. Lie, E., Ko, J. S., Choi, S. Y., Roh, J. D., Cho, Y. S., Noh, R., . & Kim, E. (2016). SALM4 suppresses excitatory synapse development by cis-inhibiting trans-synaptic SALM3-LAR adhesion. *Nature communications*, 7(1), 12328. (2016).
12. He, T. C., Sparks, A. B., Rago, C., Hermeking, H., Zawel, L., Da Costa, L. T., . & Kinzler, K. W. (1998). Identification of c-MYC as a target of the APC pathway. *Science*, 281(5382), 1509-1512.

---

## Notice

We provide information and advice to our customers on application technologies and regulatory matters to the best of our knowledge and ability, but without obligation or liability. Existing laws and regulations are to be observed in all cases by our customers. This also applies in respect to any rights of third parties. Our information and advice do not relieve our customers of their own responsibility for checking the suitability of our products for the envisaged purpose.

The information in this document is subject to change without notice and should not be construed as a commitment by the manufacturing or selling entity, or an affiliate. We assume no responsibility for any errors that may appear in this document.

### Technical Assistance

Visit the tech service page at [SigmaAldrich.com/techservice](https://SigmaAldrich.com/techservice).

### Terms and Conditions of Sale

Warranty, use restrictions, and other conditions of sale may be found at [SigmaAldrich.com/terms](https://SigmaAldrich.com/terms).

### Contact Information

For the location of the office nearest you, go to [SigmaAldrich.com/offices](https://SigmaAldrich.com/offices).

The life science business of Merck operates  
as MilliporeSigma in the U.S. and Canada.

Merck, Millipore and Sigma-Aldrich are trademarks of Merck or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.

© 2022-2025 Merck and/or its affiliates. All Rights Reserved.

23267057 Rev 09/25

