

**GOAT ANTI HUMAN IgG  
AFFINITY PURIFIED  
POLYCLONAL ANTIBODY**

<b>CATALOG NUMBER:</b>	AP112
<b>LOT NUMBER:</b>	23040144
<b>QUANTITY:</b>	2 mg
<b>CONCENTRATION:</b>	2 mg/mL
<b>SPECIFICITY:</b>	Reacts with Human IgG (H&L), as well as the light chains from all human immunoglobulin classes. Monospecific on IEP against whole human serum.
<b>IMMUNOGEN:</b>	Highly purified Human IgG isolated from pooled normal human serum.
<b>PURITY:</b>	Single band on IEP at a concentration of 10 mg/mL when run against rabbit anti goat whole serum.
<b>FORMAT:</b>	Liquid in 0.02M Phosphate Buffer, 0.25M NaCl and 0.1% sodium azide. Filtered through a 0.45 micron membrane filter.
<b>STORAGE:</b>	Store at +2-8°C for up to 6 months.

**Important Note:** *During shipment, small volumes of product will occasionally become entrapped in the seal of the product vial. For products with volumes of 200  $\mu$ L or less, we recommend gently tapping the vial on a hard surface or briefly centrifuging the vial in a tabletop centrifuge to dislodge any liquid in the container's cap.*

*For research use only; not for use as a diagnostic.*

Unless otherwise stated in our catalog or other company documentation accompanying the product(s), our products are intended for research use only and are not to be used for any other purpose, which includes but is not limited to, unauthorized commercial uses, in vitro diagnostic uses, ex vivo or in vivo therapeutic uses or any type of consumption or application to humans or animals.

©2007: Millipore Corporation. All rights reserved. No part of these works may be reproduced in any form without permission in writing.

USA & Canada • Phone: +1(800) 437-7500 • Fax: +1 (951) 676-9209  
[www.millipore.com](http://www.millipore.com)