

# QUALITY CONTROL RANGES

## Human Pancreatic Polypeptide ELISA Kit Catalog No. EZHPP40K

### Human PP Standard, Lyophilized Catalog No. E8040K

STANDARD	HUMAN PP ELISA VALUE
Lot Number: EHPP-14K	3066 pg/mL

The exact value of lyophilized standard is shown above.  
Reconstitute with **0.5mL** distilled or deionized water.

### Quality Controls 1 and 2, Lyophilized Control Catalog No. E6040K

CONTROL	LOT NUMBER	EXPECTED RANGE	UNITS
QC 1	HPP-114	125 - 261	pg/mL
QC 2	HPP-214	647 - 1344	pg/mL

Reconstitute each vial with **0.5mL** distilled or deionized water.

The life science business of Merck KGaA, Darmstadt, Germany operates as  
MilliporeSigma in the U.S. and Canada.

The vibrant M, Millipore and Sigma-Aldrich are trademarks of Merck KGaA, Darmstadt,  
Germany or its affiliates.

© 2025 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved.

Cat#: E6040K, Lot#: HPP-114.214 (Rev 17 October 2025)

