

QUALITY CONTROL RANGES



**MILLIPLEX[®] MAP Mouse Kidney Injury
Magnetic Bead Panel 1
Kit Catalog # MKI1MAG-94K
Control Catalog # MKI1-6094-1
Lot# MKI1-105 and MKI1-205**

<u>Analyte</u>	<u>QC Level</u>	<u>Expected Range</u>	<u>Units</u>
b-2-Microglobulin	Control 1	303 - 629	pg/mL
	Control 2	3332 - 6921	pg/mL
Renin	Control 1	308 - 640	pg/mL
	Control 2	3345 - 6947	pg/mL
KIM-1	Control 1	153 - 317	pg/mL
	Control 2	1687 - 3504	pg/mL
TIMP-1	Control 1	136 - 283	pg/mL
	Control 2	1568 - 3256	pg/mL
IP-10	Control 1	33 - 69	pg/mL
	Control 2	318 - 661	pg/mL
VEGF	Control 1	6.2 - 13	pg/mL
	Control 2	69 - 144	pg/mL

Note: The Quality Control Ranges are generated with overnight assay format using serum matrix provided in the kit. Quality Control Values in culture media are not tested.

MKI1-6094-1: MKI1-105.205 (21 July 2021)
EMD Millipore • 1-800-645-5476 • 1-800-645-5439 (FAX) • emdmillipore.com

The M logo is a trademark of Merck KGaA, Darmstadt, Germany.
© 2013 EMD Millipore Corporation, Billerica, MA 01821 USA.