

**MOUSE ANTI-ALPHA TUBULIN, DETYROSINATED
MONOCLONAL ANTIBODY**

CATALOG NUMBER:	MAB5566
LOT NUMBER:	
QUANTITY:	100 µL
SPECIFICITY:	Alpha tubulin, detyrosinated.
IMMUNOGEN:	Synthetic peptide GEEEGE.
ISOTYPE:	IgG _{2a}
APPLICATIONS:	Western blot. Immunohistochemistry: 1:1,000-1:3,000 on rat brain sections. Optimal working dilutions must be determined by end user.
SPECIES REACTIVITY:	All mammalian species. Other species have not yet been tested.
FORMAT:	Ascites fluid.
PRESENTATION:	Liquid. Contains no preservative.
STORAGE/HANDLING:	Maintain at -20°C in undiluted aliquots for up to 6 months after date of receipt. Avoid repeated freeze/thaw cycles.
REFERENCES:	Idriss, Haitham, <i>Cell Motility and Cytoskeleton</i> (2000) 45 :173. Mialhe, A, et al., <i>Cancer Research</i> (2001) 61 :5024. Yoshiyama, Y. et al., <i>Journal of Neuroscience</i> (2003) 23 :10662.

Important Note: During shipment, small volumes of product will occasionally become entrapped in the seal of the product vial. For products with volumes of 200 µL or less, we recommend gently tapping the vial on a hard surface or briefly centrifuging the vial in a tabletop centrifuge to dislodge any liquid in the container's cap.

For research use only; not for use as a diagnostic.

© 2005: CHEMICON International, Inc. - By CHEMICON International, Inc. All rights reserved. No part of these works may be reproduced in any form without permissions in writing.

USA & Canada • Phone: +1(800) 437-7500 • Fax: +1 (951) 676-9209 • Europe +44 (0) 23 8026 2233
Australia +61 3 9839 2000 • Germany +49-6192-207300 • ISO Registered worldwide
www.chemicon.com • custserv@chemicon.com • techserv@chemicon.com