



QUALITY CONTROL RANGES

RAT/MOUSE GROWTH HORMONE ELISA CATALOG NO. EZRMGH-45K

STANDARD, LYOPHILIZED CATALOG NO. E8045-K

STANDARD LOT NUMBER: 845-21K	RMNPY VALUE 53 ng/mL
---------------------------------	-------------------------

The exact value of lyophilized standard is shown above.

Reconstitute with **0.5 mL** distilled or deionized water.

QUALITY CONTROLS 1 AND 2, LYOPHILIZED CATALOG NO. E6045-K

CONTROL	LOT NUMBER	EXPECTED RANGE	UNITS
QC 1	RMGH-121	1.4 – 2.9	ng/mL
QC 2	RMGH-221	10.3 – 21.4	ng/mL

E6045-K: RMGH-121.221 (23 SEPTEMBER 2025)

[SigmaAldrich.com/support](https://www.sigmaaldrich.com/support) • [SigmaAldrich.com](https://www.sigmaaldrich.com)

© 2025 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved.

The vibrant M logo and MILLIPIXEL are trademarks of Merck KGaA, Darmstadt Germany or its affiliates.