



VCX 130

130 Watts

Processing Volume: 200µl – 150mL*

The **VCX 130** is a 130 watt, programmable ultrasonic processor and our most popular small volume system.

This model offers the same programmability as larger **VCX** versions.

The 1/4" probe (10 – 50ml volumes) is included and a variety of accessories are available for various applications and sample volumes.

FEATURES:

- Programmable
- Pulse mode
- 10 hour timer
- Variable power output
- Wattage and Energy display



(stand is optional)

Ordering information: Model no. **VCX 130**
130 watt ultrasonic processor with 1/4" diameter probe and tool kit.

SPECIFICATIONS

Power Output	130 watt maximum	Converter Cable	Permanently attached
Frequency	20kHz	Length	5 ft. (1.5m)
Dimensions	H x W x D: 4.5 x 9.75 x 12.5" (115 x 250 x 320mm)		
Ultrasonic Converter	Part #CV18	Tool Kit	Part #888-00026 open end wrench (2 included)
Diameter	1.25" (32mm)		
Length	5.75" (146mm)		
Standard Probe	Part #630-0435	Optional:	
Tip Diameter	1/4" (6mm)	Stand with Clamp	Part #830-00460
Length	5.4" (137mm)	Sound Enclosure	Part #830-00451
Material	Ti6Al4V titanium alloy		

*The actual processing volume of a specific probe or accessory is application specific.