



MOUSE ANTI-COXSACKIE VIRUS A9 MONOCLONAL ANTIBODY

CATALOG NUMBER:	MAB947
LOT NUMBER:	
QUANTITY:	100 μ L
SPECIFICITY:	Reacts with the Cocksackie virus A9. Neutralizes coxsackie virus A9 Griggs strain at a titer of 1:112.
ISOTYPE:	IgG _{2b}
APPLICATIONS:	IFA at 1:1,000. Dilute with buffer pH 7.5-8.0 containing protein or other stabilizing medium. Final working dilutions must be determined by end user.
FORMAT:	Ascites
STORAGE/HANDLING:	Maintain at -20°C for up to 12 months in convenient aliquots. Avoid repeated freeze/thaw cycles.

Important Note: *During shipment, small volumes of product will occasionally become entrapped in the seal of the product vial. For products with volumes of 200 μ L or less, we recommend gently tapping the vial on a hard surface or briefly centrifuging the vial in a tabletop centrifuge to dislodge any liquid in the container's cap.*

**FOR RESEARCH USE ONLY; NOT FOR USE IN DIAGNOSTIC
PROCEDURES. NOT FOR HUMAN OR ANIMAL CONSUMPTION**

Unless otherwise stated in our catalog or other company documentation accompanying the product(s), our products are intended for research use only and are not to be used for any other purpose, which includes but is not limited to, unauthorized commercial uses, in vitro diagnostic uses, ex vivo or in vivo therapeutic uses or any type of consumption or application to humans or animals.

©2003 - 2011: Millipore Corporation. All rights reserved. No part of these works may be reproduced in any form without permission in writing.