

shelf life information data sheet

Product Family: Lynx® S2S Non-Sterile

shelf life information

Expiration Date: 48 Months

Storage Conditions: Ambient Temperature

Applicable Catalog Number(s): SSCFHB1A24, SSCFHB3A24, SSCFHB2A24, SSCMHB1A24, SSCMHB3A24, SSCMHB2A24, SSCSHB1A01, SSCSHB3A01, SSCSHB2A01

Applicable Lot Number(s): Starting with 201707xxxx and later

quality attributes studied

Bacterial Challenge

Extractable TOC Test

Burst Test

The stated expiration date designates the time prior to which a lot of the product is expected to remain within the approved shelf life specification if stored under defined conditions, and after which must not be put into use.

Not all products listed have undergone stability testing to determine the products' shelf lives. The design of the stability schedule is such that only select product(s) within a product family are tested. The design assumes that the products selected are representative of the family as they utilize the same materials of construction, production processes, and packaging components as the products within the family.

The date of manufacture can be determined from the lot number on the product label as in the following examples. Four digits of the year followed by the month number and lastly the sequence number tied to production order (resets each month).

Year code	Month code	Sequence number	Complete Lot number Example	Mfg date
2014	05	5052	2014055052	May. 2014
2015	12	5054	2015125054	Dec. 2015



24-Apr-19

Thomas Brackett
Supplier Quality Engineering Manager
Third Party Finished Goods Quality

The life science business of Merck KGaA, Darmstadt, Germany operates as MilliporeSigma in the US and Canada.

The vibrant M, Millipore, and Sigma-Aldrich are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. Detailed information on trademarks is available via publicly accessible resources.
© 2018 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved.

