



QUALITY CONTROL RANGES

GLP-1 Active RIA Kit Catalog# GLP1A-35HK

GLP-1 Total RIA Kit Catalog# GLP1T-36HK

Lispro Insulin RIA Kit Catalog# LPI-16K

Control Catalog # 6016-K

<u>ANALYTE</u>	<u>CONTROL</u>	<u>LOT NUMBER</u>	<u>EXPECTED RANGE</u>	<u>UNITS</u>
GLP-1 (Active)	QC 1	123	37 – 77	pM
	QC 2	223	157 - 326	pM
GLP-1 (Total)	QC 1	123	39 – 82	pM
	QC 2	223	154 - 320	pM
Lispro Insulin	QC 1	123	13 – 26	μU/mL
	QC 2	223	58 – 121	μU/mL

6016-K: 123.223 (05 April 2022)

EMD Millipore • 1-800-645-5476 • 1-800-645-5439 (FAX) • emdmillipore.com

The M logo is a trademark of Merck KGaA, Darmstadt, Germany.

© 2013 EMD Millipore Corporation, Billerica, MA 01821 USA