

QUALITY CONTROL RANGES

RAT / MOUSE GHRELIN (TOTAL) ELISA KIT CATALOG # EZRGRT-91K

STANDARD, LYOPHILIZED CATALOG # E8091-K

STANDARD	TOTAL GHRELIN VALUE
Lot Number: ERGA-19K	10 ng/mL

Hydrate standard with **2.0mL** deionized water.

QUALITY CONTROLS 1 AND 2, LYOPHILIZED CATALOG # E6091-K

Reconstitute each vial with **0.5 mL** distilled or deionized water.

CONTROL	LOT NUMBER	EXPECTED RANGE	UNITS
QC 1	RAGR-119	0.33 – 0.68	ng/mL
QC 2	RAGR-219	3.2 – 6.7	ng/mL

E6091-K: RAGR-119.219 (17 December 2025)

SigmaAldrich.com/support • SigmaAldrich.com

© 2025 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved.