



## QUALITY CONTROL RANGES

### SARS-COV2 ANTI-S1 IgM ELISA CATALOG # EZS1M-140K

#### STANDARD, LYOPHILIZED CATALOG # E8140-K

STANDARD	SARS-COV2 ANTI-S1 IgM VALUE
LOT NUMBER: 8140-2K	100 ng/mL

The exact value of lyophilized standard is shown above.

Reconstitute with **1.0mL** distilled or deionized water.

#### QUALITY CONTROLS 1 AND 2, LYOPHILIZED CATALOG #E6140-K

Reconstitute each vial with **1.0mL** distilled or deionized water.

CONTROL	LOT NUMBER	EXPECTED RANGE	UNITS
QC 1	6140-102	5.2 - 11	ng/mL
QC 2	6140-202	20 - 41	ng/mL

E6140-K: 6140-102.202 (23 January 2023)

EMD MILLIPORE • 1-800-645-5476 • 1-800-645-5439 (FAX) [www.emdmillipore.com](http://www.emdmillipore.com)

The M logo is a trademark of Merck KGaA, Darmstadt, Germany.

© 2013 EMD Millipore Corporation, Billerica, MA 01821 USA