

3050 Spruce Street, St. Louis, MO 63103 USA
Tel: (800) 521-8956 (314) 771-5765 Fax: (800) 325-5052 (314) 771-5757
email: techservice@sial.com sigma-aldrich.com

Product Information

FabRICATOR® from *Streptococcus pyogenes* recombinant, expressed in *E. coli*

Catalog Numbers **07298** (2,000 units for cleaving 2 mg IgG) **77661** (5,000 units for cleaving 5 mg IgG) Storage Temperature –20 °C

Synonyms: IdeS, Immunoglobulin degrading enzyme

Product Description

FabRICATOR® is a unique proteolytic enzyme that cleaves IgG just below the hinge region, thereby, generating an intact F(ab')2 fragment and a Fc fragment. This enzyme is a modified cysteine protease first isolated from *Streptococcus pyogenes*. The scientific name of the enzyme is IdeS (immunoglobulin degrading enzyme), which has a biological role of circumventing the host defense.

This product is supplied as a lyophilized powder containing sodium phosphate, pH 6.6, and sodium chloride.

Purity: ≥95% (SDS-PAGE)

One unit is defined as the amount of enzyme required to fragment 95% of 1 μg of human IgG in 30 minutes at 37 °C, pH 6.6, as monitored by SDS-PAGE.

Precautions and Disclaimer

This product is for R&D use only, not for drug, household, or other uses. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

Preparation Instructions

Reconstitute 07298 in 30 μ l of ultrapure water and 77661 in 75 μ l of ultrapure water to prepare a solution with a concentration of 67 units/ μ l. To prevent microbial contamination, sodium azide can be added to the solution to a final concentration of 0.02–0.05% (w/v). After reconstitution, a FabRICATOR solution retains activity for 1 month at 2–8 °C.

Storage/Stability

The product ships at ambient temperature and storage at –20 °C is recommended. When stored at –20 °C, the protein retains activity for at least 1 year.

Procedure

Add 1 unit of FabRICATOR per 1 μ g of IgG for digestion in the recommended cleavage buffer of 50 mM sodium phosphate, pH 6.6, with 150 mM NaCl at 37 °C for 30 minutes. The recommended antibody concentration range is 0.5–10 mg/ml. The Fc fragments can be removed with protein A or protein G.

While optimal activity is obtained at pH 6.6 and 37 °C, it is possible to use a buffer with a higher pH and to increase the reaction time. Digestion can also be done at room temperature with prolonged incubation time.

References

- Björck, L., et al., Bacterial growth blocked by a synthetic peptide based on the structure of a human proteinase inhibitor. Nature, 337, 385-386 (1989).
- von Pawel-Rammingen, U., et al., IdeS, a novel streptococcal cysteine proteinase with unique specificity for immunoglobulin G. EMBO J., 21, 1607-1615 (2002).
- 3. Agniswamy, J., et al., Insight of host immune evasion mediated by two variants of group a *Streptococcus* Mac protein. J. Biol. Chem., **279**, 52789-52796 (2004).
- Vincents, B., et al., Enzymatic characterization of the streptococcal endopeptidase, IdeS, reveals that it is a cysteine protease with strict specificity for IgG cleavage due to exosite binding. Biochemistry, 43, 15540-15549 (2004).
- Wenig, K., et al., Structure of the streptococcal endopeptidase IdeS, a cysteine proteinase with strict specificity for IgG. Proc. Natl. Acad. Sci. U.S.A., 101, 17371-17376 (2004).
- Hess, J.L., et al., Immunoglobulin cleavage by the streptococcal cysteine protease IdeS can be detected using protein G capture and mass spectrometry. J. Microbiol. Methods, 70, 284-291 (2007).

- 7. Ryan, M.H., et al., Proteolysis of purified IgGs by human and bacterial enzymes in vitro and the detection of specific proteolytic fragments of endogenous IgG in rheumatoid synovial fluid. Mol. Immunol., **45**, 1837-1846 (2008).
- 8. Rispens, T., et al., Mechanism of immunoglobulin G4 Fab-arm exchange. J. Am. Chem. Soc., **133**, 10302-10311 (2011).
- Goetze, A.M., et al., Rapid LC-MS screening for IgG Fc modifications and allelic variants in blood. Mol. Immunol., 49, 338-352 (2011).
- Chevreux, G., et al., Fast analysis of recombinant monoclonal antibodies using IdeS proteolytic digestion and electrospray mass spectrometry. Anal. Biochem., 415, 212-214 (2011).
- Lišková, M., et al., Conjugation reactions in the preparations of quantum dot-based immunoluminescent probes for analysis of proteins by capillary electrophoresis. Anal. Bioanal. Chem., 400, 369-379 (2011).
- Leontyev, D., et al., Sialylation-independent mechanism involved in the amelioration of murine immune thrombocytopenia using intravenous gammaglobulin. Transfusion, 52(8), 1799-1805 (2012).

FabRICATOR is a registered trademark of Genovis AB.

RBG,GCY,MAM 03/14-1