


Anti-Arp1/Centractin (A 5601)



Immunoblot

Whole cell extract of A431 cells was separated on SDS-PAGE and blotted with Rabbit Anti-Arp1/Centractin (Product No. A 5601).

The antibody was developed with Goat Anti-Rabbit IgG Peroxidase (Product No. A 0545) and chemiluminescent substrate.

Lane A: Antibody dilution 1:1,000

Lane B: Antibody dilution 1:1,000 + 5 μ g/ml
Arp1/Centractin immunizing peptide