



QUALITY CONTROL RANGES

GLUCAGON-LIKE PEPTIDE 2 ELISA CATALOG #: EZGLP2-37K

GLUCAGON-LIKE PEPTIDE 2 STANDARD, LYOPHILIZED CATALOG #: E8037-K

STANDARD	GLP2 VALUE
LOT NUMBER: 837-18K	64 ng/mL

The exact value of lyophilized standard is shown above.

Reconstitute with **0.5mL** distilled or deionized water.

QUALITY CONTROLS 1 AND 2, LYOPHILIZED CATALOG #: E6037-K

Reconstitute each vial with **0.5mL** distilled or deionized water.

CONTROL	LOT NUMBER	EXPECTED RANGE	UNITS
QC 1	GLP2-118	2.3 – 4.8	ng/mL
QC 2	GLP2-218	10 – 21	ng/mL

E6037-K: GLP2-118.218 (07 October 2024)

EMD MILLIPORE • 1-800-645-5476 • 1-800-645-5439 (FAX) www.emdmillipore.com

The M logo is a trademark of Merck KGaA, Darmstadt, Germany.

© 2013 EMD Millipore Corporation, Billerica, MA 01821 USA