



Quality Control Ranges
MILLIPLEX[®] MAP Human Neurodegenerative Disease
Magnetic Bead Panel 2 Kit
Catalog # HNDG2MAG-36K
Control Catalog # HNDG2-6036-2
Lot: HNDG2-112.212

<u>Hormone</u>	<u>QC Level</u>	<u>Lot #</u>	<u>Expected Range</u>	<u>Units</u>
HCRP	Control I	112	0.17 - 0.43	ng/mL
	Control II	212	1.1 - 2.3	ng/mL
HA1AT	Control I	112	3.9 - 8.2	ng/mL
	Control II	212	20 - 41	ng/mL
HPEDF	Control I	112	0.88 - 1.8	ng/mL
	Control II	212	4.5 - 9.2	ng/mL
HSAP	Control I	112	2.2 - 4.6	ng/mL
	Control II	212	11 - 22	ng/mL
HMIP-4	Control I	112	0.09 - 0.18	ng/mL
	Control II	212	0.45 - 0.93	ng/mL
HCC4	Control I	112	4.4 - 9.1	ng/mL
	Control II	212	22 - 45	ng/mL

HNDG2-6036-2: HNDG2-112.212 (21 APRIL 2025)

EMD Millipore • 1-800-645-5476 • 1-800-645-5439 (FAX) • emdmillipore.com

The M logo is a trademark of Merck KGaA, Darmstadt, Germany.

© 2013 EMD Millipore Corporation, Billerica, MA 01821 USA