



QUALITY CONTROL RANGES

MILLIPLEX® MAP Human Kidney Injury Magnetic Bead Panel 4
Kit Catalog # HKI4MAG-99K
Control Catalog # HKI4-6099-4

<u>Analyte</u>	<u>QC Level</u>	<u>Lot #</u>	<u>Expected Range</u>	<u>Units</u>
EGF	Control 1	HKI4-107	0.07 - 0.14	ng/mL
	Control 2	HKI4-207	0.8 - 1.6	ng/mL
KIM-1	Control 1	HKI4-107	1.3 - 2.6	ng/mL
	Control 2	HKI4-207	12 - 26	ng/mL
OPN	Control 1	HKI4-107	8.5 - 18	ng/mL
	Control 2	HKI4-207	74 - 155	ng/mL
IP-10	Control 1	HKI4-107	0.03 - 0.07	ng/mL
	Control 2	HKI4-207	0.34 - 0.70	ng/mL
Renin	Control 1	HKI4-107	4.7 - 9.8	ng/mL
	Control 2	HKI4-207	44 - 92	ng/mL
FABP1	Control 1	HKI4-107	1.11 - 2.3	ng/mL
	Control 2	HKI4-207	9.9 - 20	ng/mL
PTH	Control 1	HKI4-107	0.15 - 0.30	ng/mL
	Control 2	HKI4-207	1.3 - 2.7	ng/mL

The Quality Control Ranges are generated using serum matrix provided in the kit.
HKI4-6099-4: HKI4-107/-207 (02 September 2025)

[SigmaAldrich.com/support](https://www.sigmaaldrich.com/support) • [SigmaAldrich.com](https://www.sigmaaldrich.com)

© 2025 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved. The vibrant M logo and MILLIPLEX are trademarks of Merck KGaA, Darmstadt Germany or its affiliates.