

Don't use printout for evaluation!!!

Supelco.

www.sigmaaldrich.com

1.14418

MQuant®

Copper Test

Kupfer-Test

Test Cuivre

Test Cobre

Cu



© 2024


Color card

Farbkarte


Carte colorimétrique

Tarjeta colorimétrica


1.



2.



3.



7.69144.1808-XXXXXX 04/2024

Merck KGaA, 64271 Darmstadt, Germany,  
Tel. +49(0)6151 72-2440  
[www.sigmaaldrich.com/colorimetry](http://www.sigmaaldrich.com/colorimetry)

EMD Millipore Corporation, 400 Summit Drive  
Burlington MA 01803, USA, Tel. +1-978-715-4321

MilliporeSigma Canada Ltd., 2149 Winston Park Dr,  
Oakville, Ontario, L6H 6J8, Canada,  
Phone: +1 800-565-1400

0

0.3

0.6

1.0

1.5

2.0


2.5

3

4

5

mg/l  
Cu



5 ml



1 x Cu-1



pH 7.0 - 9.5



5 x Cu-2



5 min



B  
A





