



CHEMICON International, Inc. 28835 Single Oak Drive • Temecula, CA 92590
Phone: 909-676-8080 • 800-437-7500 • Fax: 909-676-9209
custserv@chemicon.com • techserv@chemicon.com • www.chemicon.com

MOUSE ANTI-RAT P450 11 β -HYDROXYLASE MONOCLONAL ANTIBODY

CATALOG NUMBER: MAB6020

LOT NUMBER: 18090298

QUANTITY: 500 μ L

SPECIFICITY: Rat p450 11 β -Hydroxylase

IMMUNOGEN: Aa 272-285 (KNVYRELAE GRQQS-C)

ISOTYPE: IgG₁

CLONE NAME: 11B2-1A4

APPLICATIONS: Western blot: 1:300
Immunohistochemistry (STF on paraffin embedded tissue): 1:150
Optimal working dilutions must be determined by end user.

SPECIES REACTIVITIES: Rat

FORMAT: Supernatant

PRESENTATION: Liquid in 0.1% sodium azide.

STORAGE/HANDLING: Maintain refrigerated at 2-8°C in undiluted aliquots for up to 6/12 months.

For research use only; not for use as a diagnostic.

Important Note: *During shipment, small volumes of product will occasionally become entrapped in the seal of the product vial. For products with volumes of 200 μ L or less, we recommend gently tapping the vial on a hard surface or briefly centrifuging the vial in a tabletop centrifuge to dislodge any liquid in the container's cap.*