

QUALITY CONTROL RANGES
Mouse CVD Magnetic Bead Panel 2
Kit Catalog # MCVD2MAG-77K
Control Catalog # MCVDMG2-6077-2, Lot # MCVD2-104 and MCVD2-204



<u>Analyte</u>	<u>QC Level</u>	<u>Expected Range</u>	<u>Units</u>	<u>Analyte</u>	<u>QC Level</u>	<u>Expected Range</u>	<u>Units</u>
CXCL16	Control 1	178 – 370	pg/mL	sCD40L	Control 1	525 – 1091	pg/mL
	Control 2	1292 – 2683	pg/mL		Control 2	3801 – 7895	pg/mL
Endocan-1	Control 1	27 – 57	pg/mL	PIGF-2	Control 1	14 – 30	pg/mL
	Control 2	244 – 506	pg/mL		Control 2	111 – 230	pg/mL
LIGHT	Control 1	2629 – 5459	pg/mL	Troponin I	Control 1	2806 – 5829	pg/mL
	Control 2	19760 – 41040	pg/mL		Control 2	23007 – 47784	pg/mL
Follistatin	Control 1	1022 – 2123	pg/mL	Troponin T	Control 1	668 – 1387	pg/mL
	Control 2	7854 – 16313	pg/mL		Control 2	6093 – 12654	pg/mL
Oncostatin M	Control 1	151 – 314	pg/mL				
	Control 2	1085 – 2253	pg/mL				

MCVDMG2-6077-2: MCVD2-104.204 (17 Nov 2020)
 EMD Millipore • 1-800-645-5476 • 1-800-645-5439 (FAX) • emdmillipore.com