

75258 Yersinia selective Supplement

An antibiotic supplement recommended for the selective isolation of *Yersinia enterocolitica*.

Composition

(per vial, sufficient for 500 ml medium)

Cefsulodin 7.5 mg

Triclosan 2 mg

Novobiocin 1.25 mg

Directions

Rehydrate the contents of 1 vial aseptically with 1 ml of sterile distilled water and 1 ml ethanol. Mix gently to dissolve the contents completely and aseptically add to 500 ml of sterile, molten *Yersinia* Selective Agar base (Cat. No. 95760). Mix well and pour into sterile petri plates.

Precautions and Disclaimer

This product is for R&D use only, not for drug, household, or other uses. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

