

51962 m-CP selective Supplement I

An antibiotic supplement recommended by the Directive of the Council of the European Union 98/83/EC for the selective isolation of *Clostridium perfringens*.

Composition

(per vial, is sufficient for 500 ml medium)

D-Cycloserine 200.0mg

Polymyxin B sulphate 12.5mg

Directions

Rehydrate the contents of 1 vial in 5 ml sterile distilled water, mix well, and aseptically add to 485 ml sterile molten M-CP Agar Basis (Cat. No. 75605; cooled to < 50°C) alongwith one vial of m-CP-Selektiv-Supplement II (Cat. No. 82265). Mix well before pouring into sterile petri plates.

Precautions and Disclaimer

This product is for R&D use only, not for drug, household, or other uses. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

