


MERCK

pharmaceutical
secondary standards

Impurities & Related Compounds

Reference Materials to meet your needs in Pharma QC



Supelco[®]
Analytical Products



Merck's pharmaceutical secondary standards provide pharmaceutical laboratories and manufacturers with a convenient and cost effective alternative to making in-house working standards.

The US Food and Drug Administration (FDA) and European Pharmacopoeia (EP) recognize the use of secondary standards, or working standards, which are established with reference to the corresponding primary standard (references available on request). Using these commercially available secondary standards allows you to focus resources on core business activities, rather than consuming time and resources preparing internal working standards. In addition, our secondary standards provide larger quantities at less expense.

The most important product features are:

- **Traceability** to United States Pharmacopeia (USP) and also to European Pharmacopoeia (EP) and British Pharmacopoeia (BP) standards where available.
- Analysis is performed using **instruments validated according to GMP** using pharmacopeia monograph methods
- **Certified purity value** according to ISO 17034 and ISO/IEC 17025 using mass balance approach
- **Comprehensive certificate** according to ISO Guide 31

Our secondary standards are offered as certified reference materials (CRMs) produced in an ISO/IEC 17025 and ISO 17034 accredited facility. Their comprehensive certificates supply values for the shipped material demonstrating traceability to USP, EP and BP primary standards (the latter two occurring when a corresponding primary CRS (certified reference substance) is available).

In addition, an independent certified purity value is given. This allows the secondary standards to be used as reference materials for quantitative uses.

Comprehensive Certificate

Our certificate for these pharmaceutical secondary standards consists of up to 8 pages and contains traceability assay results, certified purity as well as information about analytical methods, chromatograms and additional analytical data. Example certificates are available for free download on the web (SigmaAldrich.com/pharmastandards). An example how the pharmacopeia traceable values are reported in the certificate is shown in **Figure 1** on page 3.

The values on the certificate are always traceable to the current Pharmacopeial lots. If a valid Pharmacopeial lot changes, then our corresponding Pharmaceutical Secondary Standard will be recertified with traceability to the new lot and a new certificate will be made available online. For this reason, the valid certificate always needs to be downloaded from the website prior to use of the material.

Please note the following: The Reference Standards of the Pharmacopeial authorities are for the specified use described in the monograph. The use of the standard beyond that limited scope is not supported by the Pharmacopeia. Other data supplied in the certificate of secondary standards is for informational purposes only and does not support Pharmacopeial use beyond that stated in the corresponding monograph.

The comprehensive certificate contains the following data and information:

- Handling and storage instructions
- Traceability results versus Pharmacopeial primary standards
- Certified purity value by mass balance (according to ISO/IEC 17025)
- Comprehensive additional analytical data

Traceability Assay Results

Comparative assay demonstrates direct traceability to Pharmacopeial Standards

Specification: 97.0–103.0% (anhydrous, USP/NF) 
98.5–101.0% (dried substance, EP)

Determination Method: HPLC
(ref., Current Compendial Monograph: Ibuprofen)

Column: Ascentis® C18, 4.6 x 250mm, 5µm 

Mobile Phase: Acetonitrile/Water/Chloroacetic acid
(600:400:4)

Column Temperature: 30° C

Flow Rate: 2 mL/min

Injection: 5 µl 

Detector Wavelength: 254 nm

Assay vs. USP Reference Standard (as is basis)

Assay Value	vs. USP LOT 
99.5% 	R024X0
	Labeled Content = 0.998 mg/mg

-  Analytical specification
-  Reference to Pharmacopeial monograph (if available) used for the certification
-  Analytical conditions of the method used
-  Reference to the Pharmacopeial lot of the primary standard
-  Assay value (traceable to the given Pharmacopeial lot)

Figure 1: Certificate excerpt of traceability assay results (vs the USP reference standard) for Ibuprofen (PHR1004).

The following list contains a partial overview of our catalog of over 2,000 products. Our portfolio of pharmaceutical secondary standards is continually growing:

Please visit our website at [SigmaAldrich.com/pharmastandards](https://www.sigmaaldrich.com/pharmastandards) for an up-to-date product listing.

Please also have a look at our offering of pharmaceutical impurity standards [SigmaAldrich.com/pharmaimpurities](https://www.sigmaaldrich.com/pharmaimpurities) and pharmaceutical proficiency testing materials [SigmaAldrich.com/ptpharma](https://www.sigmaaldrich.com/ptpharma)

Custom Services

Let us partner with your organization to provide manufacturing and packaging services customized to your specifications. Our scopes of accreditation include a broad range of technologies and competencies, ensuring that custom projects can be accommodated and with the highest quality. Your projects will be processed in a dedicated clean room production facility and under cGMP conditions appropriate to the manufacturing of reference materials.

- We can supply custom reference materials based on your specific needs and following your required specifications. These materials can be qualified and characterized at our Laramie Wyoming USA facility and can be supplied as Certified Reference Materials in accordance with ISO 17034 and ISO/IEC 17025.
- We can supply packaging services customized to your specifications. We currently offer packaging services for both solids and liquids and have blister packaging capabilities.

Capabilities

Powder Dispensing

- Serum vials: multiple sizes, crimp cap, fluoropolymer stoppers
- Glass/plastic bottles
- HEPA filtered clean-room dispensing cabinets
- Low actinic light control
- Automated powder filling
- Screw cap vials
- Clean room dispensing suites
- Humidity control
- Anaerobic Glovebox

Liquid Dispensing

- Automated flame-sealed ampuling
- Inert gas overlay
- Amber or clear ampules: 2, 5, 10, or 20 mL

Packaging

- Solids/Liquids
- Flammable materials
- Oxygen sensitive
- Toxic materials
- Light sensitive
- Hygroscopic materials

Monthly COA Updates



Customers are notified in advance of any COA changes that may occur due to re-qualification etc. This service helps users ensure that they always have the most current version of a product certificate.

Sign up for these email updates on our website at <https://www.sigmaaldrich.com/US/en/campaigns/coa-update-form>



Now Available

Primary Pharmaceutical Standards

A primary standard in metrology is a standard that is accurate enough that it is not calibrated by or subordinate to other standards. In a pharmaceutical setting, primary standards are produced and supplied by scientific organizations such as the EDQM (European Pharmacopoeia) and the US Pharmacopoeia.

Now, we provide a stress-free, time-saving and efficient route to procuring USP, EP and BP standards. Following full compliance checks, these are available for easy ordering either through the web or following your usual B2B purchasing routes.

A large number of these reference materials are currently held in stock in our warehouses, and Merck has a major program to increase this week on week. Our web pages provide live stock levels and those not in stock will generally be delivered in 3–4 weeks.

Use of Primary Pharmaceutical Standards

The use of primary reference standards is prescribed by the issuing pharmacopoeia.

USP Reference Standards are intended only for use in analytical or laboratory applications as specified in USP compendia.

EP Materials are supplied exclusively as European Pharmacopoeia Chemical Reference Substances, Biological Reference Preparations or Reference Spectra (Ph. Eur. CRS, BRP or RS) for use as reference standards in tests and assays carried out in accordance with the official methods of the European Pharmacopoeia (Ph. Eur.) and for no other purpose.

View our complete offering of primary pharmaceutical standards at [SigmaAldrich.com/pharmaceuticalstandards](https://www.sigmaaldrich.com/pharmaceuticalstandards)

Cat. No.	Description	Pack Size	CAS No.
ABACAVIR			
PHR1256-500MG	ABACAVIR SULFATE	500MG	188062-50-2
PHR2073-50MG	ABACAVIR IMPURITY E	50MG	208762-35-0
PHR2074-50MG	ABACAVIR IMPURITY F	50MG	
PHR2064-50MG	ABACAVIR RELATED COMPOUND A	50MG	124752-25-6
PHR2065-50MG	ABACAVIR RELATED COMPOUND B	50MG	
PHR2066-50MG	ABACAVIR RELATED COMPOUND C	50MG	172015-79-1
ACAMPROSATE			
PHR2821-400MG	ACAMPROSATE CALCIUM	400MG	77337-73-6
PHR2822-50MG	ACAMPROSATE RELATED COMPOUND A	50MG	3687-18-1
ACEBUTOLOL			
PHR1907-250MG	ACEBUTOLOL HYDROCHLORIDE	250MG	34381-68-5
PHR2824-50MG	ACEBUTOLOL RELATED COMPOUND A	50MG	40188-45-2
PHR2825-50MG	ACEBUTOLOL RELATED COMPOUND B	50MG	22568-64-5
ACETAMINOPHEN			
PHR1005-1G	ACETAMINOPHEN	1G	103-90-2
PHR1005-10G	ACETAMINOPHEN	10G	103-90-2
PHR1612-50MG	ACETAMINOPHEN RELATED COMPOUND A	50MG	2623-33-8
PHR1642-50MG	ACETAMINOPHEN RELATED COMPOUND B	50MG	1693-37-4
PHR1531-200MG	ACETAMINOPHEN RELATED COMPOUND C	200MG	614-80-2
PHR1530-200MG	ACETAMINOPHEN RELATED COMPOUND D	200MG	103-84-4
PHR1150-1G	ACETAMINOPHEN RELATED COMPOUND F	1G	100-02-7
PHR1149-1G	ACETAMINOPHEN RELATED COMPOUND J	1G	539-03-7
PHR1148-1G	ACETAMINOPHEN RELATED COMPOUND K	1G	123-30-8
ACETAZOLAMIDE			
PHR1908-5G	ACETAZOLAMIDE	5G	59-66-5
PHR2826-50MG	ACETAZOLAMIDE RELATED COMPOUND D	50MG	14949-00-9
ACITRETIN			
PHR1523-1G	ACITRETIN	1G	55079-83-9
PHR2919-100MG	ACITRETIN RELATED COMPOUND B	100MG	54350-48-0
ACYCLOVIR			
PHR1254-1G	ACYCLOVIR	1G	59277-89-3
ADAPALENE			
PHR1388-500MG	ADAPALENE	500MG	106685-40-9
PHR2827-100MG	ADAPALENE RELATED COMPOUND A	100MG	33626-98-1
PHR2828-100MG	ADAPALENE RELATED COMPOUND B	100MG	106685-41-0
PHR2829-50MG	ADAPALENE RELATED COMPOUND C	50MG	43109-77-9
PHR2830-50MG	ADAPALENE RELATED COMPOUND D	50MG	932033-57-3
PHR2935-50MG	ADAPALENE RELATED COMPOUND E	50MG	
PHR3168-1G	ADEMETIONINE DISULFATE TOSYLATE	1G	97540-22-2
ADENOSINE			
PHR1138-1G	ADENOSINE	1G	58-61-7
PHR2244-100MG	ADENOSINE IMPURITY F	100MG	58-96-8
PHR2412-100MG	ADENOSINE IMPURITY G	100MG	58-63-9
ALBUTEROL			
PHR1053-1G	ALBUTEROL SULFATE	1G	51022-70-9
PHR2411-100MG	ADENOSINE IMPURITY H, (GUANOSINE)	100MG	118-00-3
PHR2040-50MG	ALBUTEROL RELATED COMPOUND B	50MG	41489-89-8
ALLOPURINOL			
PHR1377-1G	ALLOPURINOL	1G	315-30-0
PHR2005-50MG	ALLOPURINOL RELATED COMPOUND A	50MG	27511-79-1
PHR2006-50MG	ALLOPURINOL RELATED COMPOUND B	50MG	22407-20-1
PHR2007-50MG	ALLOPURINOL RELATED COMPOUND C	50MG	1346604-13-4
PHR2008-50MG	ALLOPURINOL RELATED COMPOUND D	50MG	6994-25-8
PHR2009-50MG	ALLOPURINOL RELATED COMPOUND E	50MG	31055-19-3

Cat. No.	Description	Pack Size	CAS No.
ALOSETRON			
PHR2990-500MG	ALOSETRON HYDROCHLORIDE	500MG	122852-69-1
PHR2991-50MG	ALOSETRON RELATED COMPOUND A	50MG	
AMANTADINE			
PHR1711-500MG	AMANTADINE HYDROCHLORIDE	500MG	665-66-7
PHR2831-50MG	AMANTADINE RELATED COMPOUND A	50MG	935-56-8
PHR2832-50MG	AMANTADINE RELATED COMPOUND B	50MG	880-52-4
AMIFOSTINE			
PHR2940-300MG	AMIFOSTINE	300MG	
PHR2789-200MG	AMIFOSTINE DISULFIDE	200MG	10027-65-3
PHR2992-100MG	AMIFOSTINE THIOL	100MG	
AMILORIDE			
PHR1839-1G	AMILORIDE HYDROCHLORIDE	1G	17440-83-4
PHR2061-50MG	AMILORIDE RELATED COMPOUND A	50MG	1458-01-1
AMINOLEVULINIC ACID			
PHR2993-500MG	AMINOLEVULINIC ACID HYDROCHLORIDE	500MG	
PHR2995-50MG	AMINOLEVULINIC ACID RELATED COMPOUND A	50MG	
PHR2994-50MG	AMINOLEVULINIC ACID RELATED COMPOUND B	50MG	
AMIODARONE			
PHR1164-1G	AMIODARONE HYDROCHLORIDE	1G	19774-82-4
PHR2415-100MG	AMIODARONE RELATED COMPOUND D	100MG	1951-26-4
PHR2416-100MG	AMIODARONE RELATED COMPOUND E	100MG	52490-15-0
PHR1646-200MG	AMIODARONE RELATED COMPOUND H	200MG	869-24-9
AMITRAZ			
PHR2996-100MG	AMITRAZ RELATED COMPOUND A	100MG	
PHR2997-100MG	AMITRAZ RELATED COMPOUND B	100MG	33089-74-6
PHR2998-100MG	AMITRAZ RELATED COMPOUND C	100MG	
AMITRIPTYLINE			
PHR1384-1G	AMITRIPTYLINE HYDROCHLORIDE	1G	549-18-8
PHR2836-50MG	AMITRIPTYLINE RELATED COMPOUND A	50MG	1210-35-1
AMLODIPINE			
PHR1185-1G	AMLODIPINE BESYLATE	1G	111470-99-6
PHR2051-50MG	AMLODIPINE IMPURITY B	50MG	721958-72-1
PHR2057-50MG	AMLODIPINE IMPURITY D	50MG	
PHR2054-50MG	AMLODIPINE IMPURITY E	50MG	
PHR2055-50MG	AMLODIPINE IMPURITY F	50MG	
PHR2058-50MG	AMLODIPINE IMPURITY H	50MG	
PHR2050-50MG	AMLODIPINE RELATED COMPOUND A	50MG	2138811-33-1
PHR2052-50MG	AMLODIPINE RELATED COMPOUND C	50MG	43067-01-2
PHR2053-50MG	AMLODIPINE RELATED COMPOUND D	50MG	88150-62-3
PHR2056-50MG	AMLODIPINE RELATED COMPOUND K	50MG	1621516-91-3
PHR2941-50MG	AMLODIPINE RELATED COMPOUND K	50MG	
AMOXAPINE			
PHR1841-500MG	AMOXAPINE	500MG	14028-44-5
PHR2837-100MG	AMOXAPINE RELATED COMPOUND D	100MG	3158-91-6
AMOXICILLIN			
PHR1127-1G	AMOXICILLIN TRIHYDRATE	1G	61336-70-7
PHR2191-50MG	AMOXICILLIN IMPURITY D	50MG	
PHR2195-100MG	AMOXICILLIN IMPURITY H	100MG	205826-86-4
PHR2188-50MG	AMOXICILLIN RELATED COMPOUND A	50MG	
PHR2189-50MG	AMOXICILLIN RELATED COMPOUND B	50MG	26889-93-0
PHR2190-50MG	AMOXICILLIN RELATED COMPOUND C	50MG	
PHR2192-50MG	AMOXICILLIN RELATED COMPOUND E	50MG	178738-48-2
PHR2193-100MG	AMOXICILLIN RELATED COMPOUND F	100MG	126247-63-0
PHR2196-50MG	AMOXICILLIN RELATED COMPOUND I	50MG	

Cat. No.	Description	Pack Size	CAS No.
AMPICILLIN			
PHR2838-500MG	AMPICILLIN	500MG	69-53-4
PHR1424-1G	AMPICILLIN SODIUM	1G	69-52-3
ANASTROZOLE			
PHR1783-500MG	ANASTROZOLE	500MG	120511-73-1
PHR1776-50MG	AMITRYPTYLENE RELATED COMPOUND B	50MG	
ARIPIRAZOLE			
PHR1784-500MG	ARIPIRAZOLE	500MG	129722-12-9
PHR2320-100MG	ARIPIRAZOLE RELATED COMPOUND B	100MG	889443-20-3
PHR2321-100MG	ARIPIRAZOLE RELATED COMPOUND C	100MG	119532-26-2
PHR2322-100MG	ARIPIRAZOLE RELATED COMPOUND F	100MG	573691-09-5
PHR2323-100MG	ARIPIRAZOLE RELATED COMPOUND G	100MG	129722-25-4
PHR2324-50MG	ARIPIRAZOLE RELATED COMPOUND H	50MG	
PHR3000-40MG	ARTEMETHER RELATED COMPOUND B	40MG	
ARTICAINE			
PHR1882-100MG	ARTICAINE	100MG	23964-58-1
PHR1881-500MG	ARTICAINE HYDROCHLORIDE	500MG	23964-57-0
PHR2114-100MG	ARTICAINE RELATED COMPOUND A	100MG	1712677-79-6
PHR2115-50MG	ARTICAINE RELATED COMPOUND B	50MG	114176-52-2
PHR2116-100MG	ARTICAINE RELATED COMPOUND E	100MG	
ASCORBIC ACID			
PHR1008-10G	ASCORBIC ACID	10G	50-81-7
PHR1008-2G	ASCORBIC ACID	2G	50-81-7
PHR1708-30MG	ASCORBIC ACID IMP C	30MG	
PHR1709-30MG	ASCORBIC ACID IMP D	30MG	
PHR2418-100MG	ASCORBIC ACID IMPURITY F	100MG	89-65-6
PHR2942-50MG	ASCORBIC ACID IMPURITY G	50MG	
PHR2943-50MG	ASCORBIC ACID IMPURITY H	50MG	
PHR1433-1G	ASCORBYL GLUCOSIDE	1G	129499-78-1
PHR1455-2G	ASCORBYL PALMITATE	2G	137-66-6
PHR2349-300MG	ASPARAGINE ANHYDROUS	300MG	
PHR2350-300MG	ASPARAGINE MONOHYDRATE	300MG	
ASPARTAME			
PHR1381-1G	ASPARTAME	1G	22839-47-0
PHR2419-300MG	ASPARTAME ACESULFAME	300MG	106372-55-8
PHR2420-50MG	ASPARTAME RELATED COMPOUND A	50MG	5262-10-2
ATAZANAVIR			
PHR3002-500MG	ATAZANAVIR SULFATE	500MG	
PHR3001-50MG	ATAZANAVIR RELATED COMPOUND A	50MG	
A-TOCOPHEROL			
PHR1031-500MG	A-TOCOPHEROL	500MG	10191-41-0
PHR1030-500MG	A-TOCOPHERYL ACETATE	500MG	7695-91-2
PHR1029-500MG	A-TOCOPHERYL ACID SUCCINATE	500MG	4345-03-3
PHR1692-100MG	PHENYTOIN RELATED COMPOUND B	100MG	
ATORVASTATIN			
PHR1422-1G	ATORVASTATIN CALCIUM	1G	344423-98-9
PHR1868-30MG	ATORVASTATIN RELATED COMPOUND A	30MG	
PHR1869-30MG	ATORVASTATIN RELATED COMPOUND B	30MG	887196-25-0
PHR1870-30MG	ATORVASTATIN RELATED COMPOUND C	30MG	693793-53-2
PHR1871-30MG	ATORVASTATIN RELATED COMPOUND D	30MG	148146-51-4
PHR1872-30MG	ATORVASTATIN RELATED COMPOUND E	30MG	1105067-88-6
PHR1873-30MG	ATORVASTATIN RELATED COMPOUND H	30MG	125995-03-1
PHR1874-30MG	ATORVASTATIN RELATED COMPOUND I	30MG	125971-95-1

Cat. No.	Description	Pack Size	CAS No.
AZATHIOPRINE			
PHR1282-1G	AZATHIOPRINE	1G	446-86-6
PHR2841-30MG	AZATHIOPRINE RELATED COMPOUND A	30MG	4531-54-8
PHR2918-50MG	AZATHIOPRINE RELATED COMPOUND G	50MG	
PHR8758-50MG	AZATHIOPRINE RELATED COMPOUND G	50MG	
AZELASTINE			
PHR1636-1G	AZELASTINE HYDROCHLORIDE	1G	79307-93-0
PHR2842-50MG	AZELASTINE RELATED COMPOUND D	50MG	53242-88-9
PHR2843-50MG	AZELASTINE RELATED COMPOUND E	50MG	20526-97-0
AZITHROMYCIN			
PHR1088-1G	AZITHROMYCIN	1G	117772-70-0
PHR2782-30MG	AZITHROMYCIN N-OXIDE	30MG	90503-06-3
PHR2781-50MG	AZITHROMYCIN RELATED COMPOUND F	50MG	612069-28-0
BACLOFEN			
PHR1682-1G	BACLOFEN	1G	1134-47-0
PHR8947-100MG	BACLOFEN RELATED COMPOUND A	100MG	22518-27-0
BENAZEPRIL			
PHR1912-250MG	BENAZEPRIL HYDROCHLORIDE	250MG	86541-74-4
PHR8936-50MG	BENAZEPRIL RELATED COMPOUND A	50MG	215447-89-5
PHR8937-30MG	BENAZEPRIL RELATED COMPOUND B	30MG	86541-77-7
PHR2922-50MG	BENAZEPRIL RELATED COMPOUND D	50MG	
PHR2422-50MG	BENAZEPRIL RELATED COMPOUND E	50MG	88372-47-8
PHR2423-100MG	BENAZEPRIL RELATED COMPOUND F	100MG	109010-60-8
PHR2930-75MG	BENAZEPRIL RELATED COMPOUND G	75MG	103129-58-4
BENDAMUSTINE			
PHR2333-500MG	BENDAMUSTINE HYDROCHLORIDE MONOHYDRATE	500MG	
PHR2944-50MG	BENDAMUSTINE RELATED COMPOUND A	50MG	
PHR2945-40MG	BENDAMUSTINE RELATED COMPOUND B	40MG	
PHR2946-40MG	BENDAMUSTINE RELATED COMPOUND D	40MG	
PHR2947-40MG	BENDAMUSTINE RELATED COMPOUND E	40MG	
PHR2948-40MG	BENDAMUSTINE RELATED COMPOUND F	40MG	
PHR2949-40MG	BENDAMUSTINE RELATED COMPOUND G	40MG	
PHR2950-50MG	BENDAMUSTINE RELATED COMPOUND H	50MG	
PHR2951-50MG	BENDAMUSTINE RELATED COMPOUND I	50MG	
BETAMETHASONE			
PHR1398-1G	BETAMETHASONE	1G	378-44-9
PHR1779-1G	BETAMETHASONE ACETATE	1G	987-24-6
PHR2707-500MG	BETAMETHASONE BENZOATE	500MG	22298-29-9
PHR1399-1G	BETAMETHASONE DIPROPIONATE	1G	5593-20-4
PHR2424-1G	BETAMETHASONE SODIUM PHOSPHATE	1G	
PHR1780-1G	BETAMETHASONE VALERATE	1G	2152-44-5
BUDESONIDE			
PHR1178-500MG	BUDESONIDE	500MG	51333-22-3
PHR1822-50MG	BUDESONIDE RELATED COMPOUND G	50MG	137174-25-5
PHR1823-50MG	BUDESONIDE RELATED COMPOUND L	50MG	
PHR2574-300MG	BUDESONIDE ACETONIDE	300MG	638-94-8
BUPIVACAINE			
PHR1128-1G	BUPIVACAINE HYDROCHLORIDE	1G	
PHR1738-40MG	BUPIVACAINE IMPURITY E	40MG	
PHR1756-20MG	BUPIVACAINE IMPURITY B	20MG	
BUPRENORPHINE			
PHR1729-50MG	BUPRENORPHINE HYDROCHLORIDE	50MG	53152-21-9
PHR8955-50MG	BUPRENORPHINE RELATED COMPOUND A	50MG	

Cat. No.	Description	Pack Size	CAS No.
BUPROPION			
PHR3178-500MG	BUPROPION HYDROBROMIDE	500MG	
PHR1730-1G	BUPROPION HYDROCHLORIDE	1G	31677-93-7
PHR1752-30MG	BUPROPION HYDROCHLORIDE RELATED COMPOUND A	30MG	
PHR1753-30MG	BUPROPION HYDROCHLORIDE RELATED COMPOUND B	30MG	
PHR1754-30MG	BUPROPION RELATED COMPOUND C	30MG	152943-33-4
PHR1755-30MG	BUPROPION RELATED COMPOUND F	30MG	857233-13-7
BUSIPRONE			
PHR1917-500MG	BUSPIRONE HYDROCHLORIDE	500MG	33386-08-2
PHR2845-2X100MG	BUSPIRONE RELATED COMPOUND A	2X100MG	20980-22-7
PHR2846-100MG	BUSPIRONE RELATED COMPOUND G	100MG	84746-24-7
PHR2847-100MG	BUSPIRONE RELATED COMPOUND K	100MG	1075-89-4
PHR2848-100MG	BUSPIRONE RELATED COMPOUND N	100MG	257877-44-4
CAFFEINE			
PHR1009-1G	CAFFEINE	1G	58-08-2
PHR2038-50MG	CAFFEINE IMPURITY B	50MG	7597-60-6
PHR1999-50MG	CAFFEINE IMPURITY C	50MG	519-32-4
PHR1951-30MG	CAFFEINE IMPURITY E NITRATE	30MG	
PHR2000-50MG	CAFFEINE IMPURITY F	50MG	611-59-6
CANDESARTAN			
PHR1854-200MG	CANDESARTAN CILEXETIL	200MG	145040-37-5
PHR2035-20MG	CANDESARTAN CILEXETIL IMPURITY C	20MG	
PHR2036-20MG	CANDESARTAN CILEXETIL IMPURITY E	20MG	914613-35-7
PHR2022-20MG	CANDESARTAN CILEXETIL IMPURITY I	20MG	139481-69-9
PHR2017-20MG	CANDESARTAN CILEXETIL RELATED COMPOUND A	20MG	139481-58-6
PHR2018-20MG	CANDESARTAN CILEXETIL RELATED COMPOUND B	20MG	869631-11-8
PHR2019-20MG	CANDESARTAN CILEXETIL RELATED COMPOUND D	20MG	
PHR2020-20MG	CANDESARTAN CILEXETIL RELATED COMPOUND F	20MG	
PHR2021-20MG	CANDESARTAN CILEXETIL RELATED COMPOUND G	20MG	139481-59-7
CAPECITABINE			
PHR1405-1G	CAPECITABINE	1G	154361-50-9
PHR2849-50MG	CAPECITABINE RELATED COMPOUND A	50MG	66335-38-4
PHR2850-50MG	CAPECITABINE RELATED COMPOUND B	50MG	3094-09-5
CARBAMAZEPINE			
PHR1067-1G	CARBAMAZEPINE	1G	298-46-4
PHR2237-50MG	CARBAMAZEPINE IMPURITY B	50MG	
PHR2330-50MG	CARBAMAZEPINE IMPURITY C	50MG	1219170-51-0
PHR2433-100MG	CARBAMAZEPINE IMPURITY F	100MG	33948-22-0
PHR1918-50MG	CARBAMAZEPINE RELATED COMPOUND A	50MG	3564-73-6
PHR1197-500MG	CARBAMAZEPINE RELATED COMPOUND B	500MG	
CARPROFEN			
PHR1452-1G	CARPROFEN	1G	53716-49-7
PHR1564-200MG	CARPROFEN RELATED COMPOUND A	200MG	86-74-8
CARVEDILOL			
PHR1265-1G	CARVEDILOL	1G	72956-09-3
PHR2952-40MG	CARVEDILOL RELATED COMPOUND B	40MG	
PHR2953-40MG	CARVEDILOL RELATED COMPOUND C	40MG	
PHR2434-50MG	CARVEDILOL RELATED COMPOUND D	50MG	51997-51-4
PHR2435-50MG	CARVEDILOL RELATED COMPOUND E	50MG	
CEFAZOLIN			
PHR1291-500MG	CEFAZOLIN	500MG	25953-19-9
PHR8907-100MG	CEFAZOLIN EP IMPURITY A	100MG	
PHR2437-100MG	CEFAZOLIN IMPURITY E	100MG	29490-19-5
PHR2436-100MG	CEFAZOLIN IMPURITY H	100MG	957-68-6
PHR2681-100MG	CEFAZOLIN SODIUM SALT	100MG	27164-46-1

Cat. No.	Description	Pack Size	CAS No.
CEFDINIR			
PHR2665-25MG	CEFDINIR	25MG	91832-40-5
PHR2665-500MG	CEFDINIR	500MG	91832-40-5
PHR2666-25MG	CEFDINIR RELATED COMPOUND A	25MG	
PHR2667-20MG	CEFDINIR RELATED COMPOUND B	20MG	79350-10-0
CEFEPIME			
PHR1763-1G	CEFEPIME HYDROCHLORIDE	1G	123171-59-5
PHR2851-50MG	CEFEPIME RELATED COMPOUND D	50MG	65872-41-5
CELECOXIB			
PHR1683-1G	CELECOXIB	1G	169590-42-5
PHR1706-20MG	CELECOXIB RELATED COMPOUND A	20MG	170570-01-1
PHR1707-20MG	CELECOXIB RELATED COMPOUND B	20MG	331943-04-5
CETRIZINE			
PHR1656-1G	CETRIZINE HYDROCHLORIDE	1G	83881-52-1
PHR1710-40MG	CETRIZINE HYDROCHLORIDE RELATED COMPOUND A	40MG	246870-46-2
CHLORHEXIDINE			
PHR1421-1G	CHLORHEXIDINE	1G	55-56-1
PHR1222-500MG	CHLORHEXIDINE ACETATE	500MG	56-95-1
PHR1294-1G	CHLORHEXIDINE DIGLUCONATE SOLUTION	1G	18472-51-0
PHR1459-500MG	CHLORHEXIDINE HYDROCHLORIDE	500MG	3697-42-5
PHR1200-1G	CHLORHEXIDINE IMP P	1G	106-47-8
CHLOROQUINE PHOSPHATE			
PHR1258-1G	CHLOROQUINE PHOSPHATE	1G	50-63-5
PHR2438-50MG	CHLOROQUINE RELATED COMPOUND A	50MG	86-98-6
PHR2933-40MG	CHLOROQUINE RELATED COMPOUND G	40MG	
CHLORPHENIRAMINE			
PHR1016-500MG	CHLORPHENIRAMINE MALEATE	500MG	113-92-8
PHR2105-100MG	CHLORPHENAMINE N-OXIDE MALEATE	100MG	120244-82-8
PHR2004-50MG	CHLORPHENIRAMINE IMPURITY A	50MG	1246816-57-8
PHR2103-100MG	CHLORPHENAMINE IMPURITY D	100MG	
PHR1647-200MG	CHLORPHENIRAMINE RELATED COMPOUND B	200MG	1202-34-2
PHR1688-50MG	CHLORPHENIRAMINE RELATED COMPOUND C	50MG	22630-25-7
CHOLECALCIFEROL			
PHR1237-500MG	CHOLECALCIFEROL	500MG	67-97-0
PHR2098-5X100MG	CHOLECALCIFEROL	5X100MG	
PHR2439-100MG	CHOLECALCIFEROL IMPURITY B	100MG	434-16-2
CICLOPIROX			
PHR1920-500MG	CICLOPIROX	500MG	29342-05-0
PHR2853-250MG	CICLOPIROX OLAMINE	250MG	41621-49-2
PHR2854-50MG	CICLOPIROX RELATED COMPOUND B	50MG	14818-35-0
CILOSTAZOL			
PHR1503-1G	CILOSTAZOL	1G	73963-72-1
PHR2855-100MG	CILOSTAZOL RELATED COMPOUND A	100MG	54197-66-9
PHR1653-100MG	CILOSTAZOL RELATED COMPOUND B	100MG	73963-62-9
PHR1676-100MG	CILOSTAZOL RELATED COMPOUND C	100MG	865792-18-3
CIMETIDINE			
PHR1075-1G	CIMETIDINE	1G	51481-61-9
PHR2365-100MG	CIMETIDINE GUANIDINE DIHYDROCHLORIDE	100MG	
PHR1089-1G	CIMETIDINE HYDROCHLORIDE	1G	70059-30-2
PHR2220-100MG	CIMETIDINE IMPURITY A	100MG	52378-40-2
PHR2250-50MG	CIMETIDINE IMPURITY B	50MG	73785-35-0
PHR2251-100MG	CIMETIDINE IMPURITY C	100MG	52568-80-6
PHR2242-25MG	CIMETIDINE IMPURITY D	25MG	59660-24-1
PHR2243-25MG	CIMETIDINE IMPURITY E	25MG	54237-72-8

Cat. No.	Description	Pack Size	CAS No.
PHR2258-50MG	CIMETIDINE IMPURITY F	50MG	
PHR2259-50MG	CIMETIDINE IMPURITY G	50MG	31857-31-5
PHR2255-50MG	CIMETIDINE IMPURITY H	50MG	74886-59-2
PHR2260-50MG	CIMETIDINE IMPURITY I	50MG	29636-87-1
PHR2252-100MG	CIMETIDINE IMPURITY J	100MG	38603-72-4
PHR2767-1G	CIMETIDINE, POLYMORPH B	1G	
CIPROFLOXACIN			
PHR1167-1G	CIPROFLOXACIN	1G	85721-33-1
PHR2440-200MG	CIPROFLOXACIN FORMAMIDE	200MG	93594-39-9
PHR1044-1G	CIPROFLOXACIN HYDROCHLORIDE	1G	86393-32-0
PHR1174-500MG	CIPROFLOXACIN IMP A	500MG	86393-33-1
PHR1910-50MG	CIPROFLOXACIN RELATED COMPOUND A	50MG	526204-10-4
CITALOPRAM			
PHR1640-1G	CITALOPRAM HYDROBROMIDE	1G	59729-32-7
PHR1901-50MG	CITALOPRAM RELATED COMPOUND A	50MG	64372-56-1
PHR1902-50MG	CITALOPRAM RELATED COMPOUND C	50MG	
PHR1903-30MG	CITALOPRAM RELATED COMPOUND D	30MG	1188264-72-3
PHR1904-50MG	CITALOPRAM RELATED COMPOUND E	50MG	62498-71-9
PHR1863-200MG	CITALOPRAM RELATED COMPOUND F	200MG	55011-89-7
PHR1905-25MG	CITALOPRAM RELATED COMPOUND G	25MG	1332724-08-9
PHR1905-50MG	CITALOPRAM RELATED COMPOUND G	50MG	1332724-08-9
PHR1906-50MG	CITALOPRAM RELATED COMPOUND H	50MG	479065-02-6
CLOPIDOGREL			
PHR1431-1G	CLOPIDOGREL BISULFATE	1G	120202-66-6
PHR1990-40MG	CLOPIDOGREL RELATED COMPOUND A	40MG	144750-42-5
PHR1991-40MG	CLOPIDOGREL RELATED COMPOUND B	40MG	2137493-39-9
PHR1992-40MG	CLOPIDOGREL RELATED COMPOUND C	40MG	120202-71-3
CLOTRIMAZOLE			
PHR1058-1G	CLOTRIMAZOLE	1G	23593-75-1
PHR3005-50MG	CLOTRIMAZOLE IMPURITY F	50MG	15469-97-3
PHR2700-50MG	CLOTRIMAZOLE RELATED COMPOUND A	50MG	66774-02-5
PHR1183-500MG	CLOTRIMAZOLE IMPURITY E	500MG	
CYPROHEPTADINE			
PHR2684-500MG	CYPROHEPTADINE HYDROCHLORIDE	500MG	41354-29-4
PHR2685-50MG	CYPROHEPTADINE RELATED COMPOUND A	50MG	256-81-5
PHR2686-50MG	CYPROHEPTADINE RELATED COMPOUND C	50MG	3967-32-6
CAMPBOR			
PHR1119-1G	D-CAMPBOR	1G	464-49-3
PHR2427-100MG	D-CAMPBOR IMPURITY A	100MG	80-56-8
PHR2425-100MG	D-CAMPBOR IMPURITY D	100MG	7787-20-4
PHR2430-100MG	D-CAMPBOR IMPURITY F	100MG	2217-02-9
PHR2431-100MG	D-CAMPBOR IMPURITY I	100MG	124-76-5
PHR2429-100MG	D-CAMPBOR IMPURITY J	100MG	464-45-9
DANTROLENE			
PHR2291-500MG	DANTROLENE	500MG	7261-97-4
PHR2819-200MG	DANTROLENE SODIUM	200MG	24868-20-0
PHR2292-100MG	DANTROLENE RELATED COMPOUND A	100MG	
PHR2293-100MG	DANTROLENE RELATED COMPOUND B	100MG	
PHR2294-100MG	DANTROLENE RELATED COMPOUND C	100MG	
DES Loratadine			
PHR1680-1G	DES Loratadine	1G	100643-71-8
PHR2317-100MG	DES Loratadine IMPURITY A	100MG	
PHR2314-50MG	DES Loratadine RELATED COMPOUND A	50MG	117796-50-6
PHR2315-50MG	DES Loratadine RELATED COMPOUND B	50MG	432543-89-0
PHR2316-100MG	DES Loratadine RELATED COMPOUND F	100MG	

Cat. No.	Description	Pack Size	CAS No.
DEXAMETHASONE			
PHR1526-1G	DEXAMETHASONE	1G	50-02-2
PHR1527-1G	DEXAMETHASONE ACETATE	1G	1177-87-3
PHR1628-500MG	DEXAMETHASONE PHOSPHATE	500MG	
PHR1768-1G	DEXAMETHASONE SODIUM PHOSPHATE	1G	2392-39-4
DEXTRONETHORPHAN			
PHR1018-500MG	DEXTROMETHORPHAN HYDROBROMIDE	500MG	6700-34-1
PHR3300-50MG	DEXTROMETHORPHAN IMPURITY A	50MG	
PHR1976-30MG	DEXTROMETHORPHAN IMPURITY D	30MG	1453167-99-1
PHR2206-100MG	DEXTROMETHORPHAN N-OXIDE HYDROCHLORIDE	100MG	1177419-85-0
DEXTROSE			
PHR1000-1G	DEXTROSE	1G	50-99-7
PHR2085-100MG	DEXTROSE IMPURITY B	100MG	499-40-1
PHR2086-1G	DEXTROSE IMPURITY C	1G	1109-28-0
DICLOFENAC			
PHR1144-1G	DICLOFENAC SODIUM	1G	15307-79-6
PHR2441-500MG	DICLOFENAC POTASSIUM	500MG	15307-81-0
PHR2444-100MG	DICLOFENAC IMPURITY B	100MG	22121-58-0
PHR2445-100MG	DICLOFENAC IMPURITY C	100MG	27204-57-5
PHR2447-100MG	DICLOFENAC IMPURITY D	100MG	127792-23-8
PHR2443-100MG	DICLOFENAC IMPURITY E	100MG	59-48-3
PHR2446-100MG	DICLOFENAC IMPURITY F	100MG	560075-65-2
PHR2442-100MG	DICLOFENAC RELATED COMPOUND A	100MG	15362-40-0
DILTIAZEM			
PHR2856-500MG	DILTIAZEM HYDROCHLORIDE	500MG	33286-22-5
PHR8702-50MG	DILTIAZEM RELATED COMPOUND A	50MG	
DIPHENHYDRAMINE			
PHR2268-1G	DIPHENHYDRAMINE CITRATE	1G	88637-37-0
PHR1015-1G	DIPHENHYDRAMINE HYDROCHLORIDE	1G	147-24-0
PHR1897-50MG	DIPHENHYDRAMINE IMPURITY B	50MG	4024-34-4
PHR2071-100MG	DIPHENHYDRAMINE IMPURITY D	100MG	91-01-0
PHR1993-50MG	DIPHENHYDRAMINE RELATED COMPOUND A	50MG	53499-40-4
PHR1899-50MG	DIPHENHYDRAMINE RELATED COMPOUND B	50MG	19574-66-4
PHR1898-50MG	DIPHENHYDRAMINE N-OXIDE HYDROCHLORIDE	50MG	13168-00-8
DONEPEZIL			
PHR1584-1G	DONEPEZIL HYDROCHLORIDE	1G	120011-70-3
PHR2024-50MG	DONEPEZIL-N-OXIDE	50MG	120013-84-5
PHR2023-50MG	DONEPEZIL RELATED COMPOUND A	50MG	145546-80-1
DOXEPIN			
PHR8941-100MG	DOXEPIN RELATED COMPOUND A	100MG	
PHR8924-100MG	DOXEPIN RELATED COMPOUND B	100MG	4504-88-5
PHR8923-50MG	DOXEPIN RELATED COMPOUND C	50MG	2887-91-4
DULOXETINE			
PHR1865-500MG	DULOXETINE HYDROCHLORIDE	500MG	136434-34-9
PHR2226-25MG	DULOXETINE RELATED COMPOUND A	25MG	910138-96-4
PHR2166-50MG	DULOXETINE RELATED COMPOUND F	50MG	
PHR2167-50MG	DULOXETINE RELATED COMPOUND H	50MG	199191-66-7
ECONAZOLE			
PHR9029-100MG	ECONAZOLE NITRATE SALT	100MG	24169-02-6
PHR2232-500MG	ECONAZOLE IMPURITY D	500MG	24168-96-5
PHR2231-100MG	ECONAZOLE RELATED COMPOUND A	100MG	24155-42-8
EDETATE DISODIUM			
PHR1068-1G	EDETATE DISODIUM	1G	6381-92-6
PHR1277-500MG	EDETATE DISODIUM IMP A	500MG	139-13-9

Cat. No.	Description	Pack Size	CAS No.
EFAVIRENZ			
PHR2072-500MG	EFAVIRENZ	500MG	154598-52-4
PHR2068-50MG	EFAVIRENZ RELATED COMPOUND A	50MG	209414-27-7
PHR2069-50MG	EFAVIRENZ RELATED COMPOUND B	50MG	440124-96-9
PHR2070-50MG	EFAVIRENZ RELATED COMPOUND C	50MG	391860-73-4
PHR1039-1G	ERYTHROMYCIN	1G	114-07-8
PHR1610-250MG	ERYTHROMYCIN ETHYLSUCCINATE	250MG	1264-62-6
ESTRADIOL			
PHR1353-1G	ESTRADIOL	1G	50-28-2
PHR2228-50MG	ESTRADIOL 3-METHYL ETHER	50MG	1035-77-4
PHR2858-500MG	ESTRADIOL BENZOATE	500MG	50-50-0
PHR2839-500MG	ESTRADIOL CYPIONATE	500MG	313-06-4
PHR2227-50MG	ESTRADIOL IMPURITY L	50MG	57-91-0
PHR2859-25MG	ESTRADIOL RELATED COMPOUND C	25MG	571-92-6
PHR2857-200MG	ESTRADIOL VALERATE	200MG	979-32-8
ESTRIOL			
PHR2800-200MG	ESTRIOL	200MG	50-27-1
PHR8707-50MG	ESTRIOL RELATED COMPOUND A	50MG	
EXEMESTANE			
PHR1634-1G	EXEMESTANE	1G	107868-30-4
PHR1894-30MG	EXEMESTANE RELATED COMPOUND A	30MG	19457-55-7
PHR1895-30MG	EXEMESTANE RELATED COMPOUND B	30MG	184972-09-6
PHR1896-30MG	EXEMESTANE RELATED COMPOUND D	30MG	
FAMCICLOVIR			
PHR1696-1G	FAMCICLOVIR	1G	104227-87-4
PHR2860-50MG	FAMCICLOVIR RELATED COMPOUND E	50MG	452-06-2
PHR2861-50MG	FAMCICLOVIR RELATED COMPOUND F	50MG	10310-21-1
FAMOTIDINE			
PHR1055-1G	FAMOTIDINE	1G	76824-35-6
PHR2909-50MG	FAMOTIDINE RELATED COMPOUND B	50MG	89268-62-2
PHR2910-50MG	FAMOTIDINE RELATED COMPOUND C	50MG	76824-17-4
PHR2911-50MG	FAMOTIDINE RELATED COMPOUND D	50MG	76824-16-3
PHR2912-50MG	FAMOTIDINE RELATED COMPOUND E	50MG	129083-44-9
PHR2908-50MG	FAMOTIDINE RELATED COMPOUND F	50MG	107880-74-0
FENBENDAZOLE			
PHR1832-200MG	FENBENDAZOLE	200MG	43210-67-9
PHR2862-50MG	FENBENDAZOLE RELATED COMPOUND A	50MG	10605-21-7
PHR2863-50MG	FENBENDAZOLE RELATED COMPOUND B	50MG	20367-38-8
FENOFIBRATE			
PHR1246-500MG	FENOFIBRATE	500MG	49562-28-9
PHR2864-50MG	FENOFIBRATE RELATED COMPOUND A	50MG	42019-78-3
FEXOFENADINE			
PHR2709-500MG	FEXOFENADINE	500MG	885946-90-7
PHR1685-1G	FEXOFENADINE HYDROCHLORIDE	1G	153439-40-8
PHR9103-50MG	FEXOFENADINE METHYL ESTER	50MG	154825-96-4
PHR2454-50MG	FEXOFENADINE RELATED COMPOUND A	50MG	76811-98-8
FLUCONAZOLE			
PHR1160-1G	FLUCONAZOLE	1G	86386-73-4
PHR1900-50MG	FLUCONAZOLE IMPURITY A	50MG	89429-59-4
PHR2046-30MG	FLUCONAZOLE IMPURITY E	30MG	86404-63-9
PHR2047-30MG	FLUCONAZOLE IMPURITY F	30MG	118689-07-9
PHR2048-30MG	FLUCONAZOLE IMPURITY G	30MG	86386-76-7
PHR2037-30MG	FLUCONAZOLE IMPURITY H	30MG	
PHR2049-30MG	FLUCONAZOLE IMPURITY I	30MG	

Cat. No.	Description	Pack Size	CAS No.
PHR1651-30MG	FLUCONAZOLE RELATED COMPOUND A	30MG	871550-15-1
PHR1862-20MG	FLUCONAZOLE RELATED COMPOUND B	20MG	81886-51-3
PHR1826-20MG	FLUCONAZOLE RELATED COMPOUND C	20MG	514222-44-7
FLUNIXIN			
PHR1442-1G	FLUNIXIN MEGLUMINE	1G	42461-84-7
PHR2305-50MG	FLUNIXIN RELATED COMPOUND A	50MG	
PHR2823-50MG	FLUNIXIN RELATED COMPOUND B	50MG	54396-44-0
PHR2413-50MG	FLUNIXIN RELATED COMPOUND C	50MG	1452-94-4
FLUOROURACIL			
PHR1227-500MG	FLUOROURACIL	500MG	51-21-8
PHR2865-50MG	FLUOROURACIL RELATED COMPOUND A	50MG	67-52-7
PHR2866-50MG	FLUOROURACIL RELATED COMPOUND B	50MG	496-76-4
PHR2867-50MG	FLUOROURACIL RELATED COMPOUND E	50MG	1820-81-1
PHR2868-50MG	FLUOROURACIL RELATED COMPOUND F	50MG	56177-80-1
FLUOXETINE			
PHR1394-1G	FLUOXETINE HYDROCHLORIDE	1G	56296-78-7
PHR1952-30MG	FLUOXETINE RELATED COMPOUND A	30MG	79088-29-2
PHR2041-5ML	FLUOXETINE RELATED COMPOUND B	5ML	23580-89-4
PHR1953-30MG	FLUOXETINE RELATED COMPOUND C	30MG	1026723-45-4
FLURBIPROFEN			
PHR1499-500MG	FLURBIPROFEN	500MG	5104-49-4
PHR2710-100MG	FLURBIPROFEN RELATED COMPOUND A	100MG	6341-72-6
PHR2711-500MG	FLURBIPROFEN SODIUM	500MG	56767-76-1
FLUTICASONE			
PHR1702-200MG	FLUTICASONE PROPIONATE	200MG	80474-14-2
PHR8918-100MG	FLUTICASONE EP IMPURITY I	100MG	
PHR2746-25MG	FLUTICASONE HYDROXY ACID/THIO ACID PROPI	25MG	
PHR1842-50MG	FLUTICASONE RELATED COMPOUND B	50MG	
PHR1843-10MG	FLUTICASONE RELATED COMPOUND C	10MG	
PHR1844-30MG	FLUTICASONE RELATED COMPOUND D	30MG	73205-13-7
PHR1845-50MG	FLUTICASONE RELATED COMPOUND F	50MG	1219174-94-3
PHR1846-30MG	FLUTICASONE RELATED COMPOUND J	30MG	28416-82-2
FOLIC ACID			
PHR1035-1G	FOLIC ACID	1G	59-30-3
PHR1334-500MG	FOLIC ACID IMPURITY A	500MG	4271-30-1
PHR2001-20MG	FOLIC ACID IMPURITY B	20MG	35011-47-3
PHR2409-50MG	FOLIC ACID IMPURITY C	50MG	
PHR2002-100MG	FOLIC ACID IMPURITY D	100MG	119-24-4
PHR2003-20MG	FOLIC ACID IMPURITY F	20MG	1391194-56-1
FOSCARNET			
PHR1436-1G	FOSCARNET SODIUM	1G	34156-56-4
PHR2869-50MG	FOSCARNET RELATED COMPOUND B	50MG	55920-24-6
PHR2870-0.5ML	FOSCARNET RELATED COMPOUND D	0.5ML	
FOSINOPRIL			
PHR8945-30MG	FOSINOPRIL RELATED COMPOUND C	30MG	
PHR8710-30MG	FOSINOPRIL RELATED COMPOUND G	30MG	
PHR8711-50MG	FOSINOPRIL RELATED COMPOUND H	50MG	
FUROSEMIDE			
PHR1057-1G	FUROSEMIDE	1G	54-31-9
PHR2457-100MG	FUROSEMIDE IMPURITY E	100MG	50-84-0
GABAPENTIN			
PHR1049-1G	GABAPENTIN	1G	60142-96-3
PHR1322-100MG	GABAPENTIN RELATED COMPOUND A	100MG	64744-50-9
PHR2012-20MG	GABAPENTIN RELATED COMPOUND D	20MG	1076198-17-8

Cat. No.	Description	Pack Size	CAS No.
GLIMEPIRIDE			
PHR1617-1G	GLIMEPIRIDE	1G	93479-97-1
PHR1726-40MG	GLIMEPIRIDE SULFONAMIDE	40MG	
PHR1728-40MG	GLIMEPIRIDE 3-ISOMER	40MG	
PHR1725-40MG	GLIMEPIRIDE RELATED COMPOUND A	40MG	684286-46-2
PHR1727-40MG	GLIMEPIRIDE RELATED COMPOUND C	40MG	119018-30-3
GLIPIZIDE			
PHR1818-500MG	GLIPIZIDE	500MG	29094-61-9
PHR1813-50MG	GLIPIZIDE RELATED COMPOUND A	50MG	33288-71-0
PHR1814-50MG	GLIPIZIDE RELATED COMPOUND B	50MG	1346600-54-1
PHR1815-50MG	GLIPIZIDE RELATED COMPOUND C	50MG	
GLUCOSAMINE			
PHR1199-500MG	GLUCOSAMINE HYDROCHLORIDE	500MG	66-84-2
PHR2290-50MG	GLUCOSAMINE IMPURITY B	50MG	13185-73-4
PHR2460-100MG	GLUCOSAMINE IMPURITY E	100MG	67-47-0
GLUTATHIONE			
PHR1359-500MG	GLUTATHIONE	500MG	70-18-8
PHR2462-100MG	GLUTATHIONE IMPURITY A	100MG	19246-18-5
PHR2461-100MG	GLUTATHIONE IMPURITY C	100MG	27025-41-8
PHR2463-100MG	GLUTATHIONE IMPURITY D	100MG	636-58-8
GLYCOPYRROLATE			
PHR2464-500MG	GLYCOPYRROLATE	500MG	51186-83-5
PHR2712-50MG	GLYCOPYRROLATE ERYTHRO ISOMER	50MG	
PHR2465-100MG	GLYCOPYRROLATE RELATED COMPOUND A	100MG	
PHR2466-100MG	GLYCOPYRROLATE RELATED COMPOUND B	100MG	13118-11-1
PHR2467-50MG	GLYCOPYRROLATE RELATED COMPOUND C	50MG	427-49-6
PHR2713-30MG	GLYCOPYRROLATE RELATED COMPOUND I	30MG	
PHR2468-100MG	GLYCOPYRROLATE RELATED COMPOUND L	100MG	19833-96-6
GUAIACOL			
PHR1136-1.5G	GUAIACOL	1.5G	90-05-1
PHR2474-100MG	GUAIACOL IMPURITY A	100MG	120-80-9
PHR2469-100MG	GUAIACOL IMPURITY C	100MG	91-16-7
PHR2475-100MG	GUAIACOL IMPURITY E	100MG	93-58-3
PHR2471-100MG	GUAIACOL IMPURITY F	100MG	
PHR2473-100MG	GUAIACOL IMPURITY G	100MG	150-76-5
PHR2472-100MG	GUAIACOL IMPURITY H	100MG	150-19-6
GUAIFENESIN			
PHR1027-1G	GUAIFENESIN	1G	93-14-1
PHR2325-50MG	GUAIFENESIN IMPURITY C	50MG	
PHR2326-50MG	GUAIFENESIN IMPURITY D	50MG	16929-60-5
PHR1751-500MG	GUANFACINE HYDROCHLORIDE	500MG	29110-48-3
PHR1543-50MG	GUAIFENESIN IMPURITY B	50MG	14007-09-1
HALOPERIDOL			
PHR1724-400MG	HALOPERIDOL	400MG	52-86-8
PHR8951-25MG	HALOPERIDOL RELATED COMPOUND A	25MG	67987-08-0
PHR8714-25MG	HALOPERIDOL RELATED COMPOUND B	25MG	
HYDROCODONE			
PHR8980-125MG	HYDROCODONE	125MG	
PHR8981-100MG	HYDROCODONE BITARTRATE	100MG	
PHR8982-100MG	HYDROCODONE RELATED COMPOUND A	100MG	
HYDROCORTISONE			
PHR1014-500MG	HYDROCORTISONE	500MG	50-23-7
PHR1606-1G	HYDROCORTISONE ACETATE	1G	50-03-3
PHR3010-50MG	HYDROCORTISONE EP IMPURITY O	50MG	

Cat. No.	Description	Pack Size	CAS No.
PHR1926-500MG	HYDROCORTISONE HEMISUCCINATE	500MG	83784-20-7
PHR2156-100MG	HYDROCORTISONE IMPURITY D	100MG	53-35-0
PHR2157-100MG	HYDROCORTISONE IMPURITY E	100MG	
PHR2158-100MG	HYDROCORTISONE IMPURITY F	100MG	152-58-9
PHR2159-50MG	HYDROCORTISONE IMPURITY G	50MG	14760-49-7
PHR2274-100MG	HYDROCORTISONE IMPURITY H	100MG	
PHR2275-100MG	HYDROCORTISONE IMPURITY I	100MG	
PHR2160-50MG	HYDROCORTISONE IMPURITY J	50MG	
PHR2161-50MG	HYDROCORTISONE IMPURITY K	50MG	640-87-9
PHR2162-50MG	HYDROCORTISONE IMPURITY L	50MG	641-77-0
PHR2163-50MG	HYDROCORTISONE IMPURITY M	50MG	566-35-8
PHR2276-100MG	HYDROCORTISONE IMPURITY N	100MG	
PHR2276-25MG	HYDROCORTISONE IMPURITY N	25MG	
PHR3185-500MG	HYDROCORTISONE VALERATE	500MG	
HYDROXYZINE			
PHR2714-500MG	HYDROXYZINE PAMOATE	500MG	10246-75-0
PHR2872-50MG	HYDROXYZINE RELATED COMPOUND A	50MG	
IBUPROFEN			
PHR1004-1G	IBUPROFEN	1G	15687-27-1
PHR1146-500MG	IBUPROFEN RELATED COMPOUND C	500MG	38861-78-8
PHR1933-20MG	IBUPROFEN IMPURITY A	20MG	66622-47-7
PHR1934-20MG	IBUPROFEN IMPURITY B	20MG	
PHR1964-30MG	IBUPROFEN IMPURITY F	30MG	65322-85-2
PHR1978-50MG	IBUPROFEN IMPURITY J	50MG	65813-55-0
PHR1965-50MG	IBUPROFEN IMPURITY G	50MG	1391054-15-1
PHR2028-50MG	IBUPROFEN IMPURITY H	50MG	
PHR1935-20MG	IBUPROFEN IMPURITY N	20MG	3585-52-2
PHR2145-50MG	IBUPROFEN 1,2,3-PROPANETRIOL ESTERS	50MG	
PHR2144-50MG	IBUPROFEN 1,2-PROPYLENE GLYCOL ESTERS	50MG	
PHR2142-100MG	IBUPROFEN 1,3-BUTYLENE GLYCOL ESTERS	100MG	
PHR2143-50MG	IBUPROFEN 2,3-BUTYLENE GLYCOL ESTER	50MG	
PHR2138-50MG	IBUPROFEN ETHYL ESTER	50MG	41283-72-1
PHR2481-100MG	IBUPROFEN INTERMEDIATE	100MG	1368714-46-8
PHR2137-1ML	IBUPROFEN METHYL ESTER	1ML	61566-34-5
PHR2137-50MG	IBUPROFEN METHYL ESTER	50MG	61566-34-5
IMIQUIMOD			
PHR2716-500MG	IMIQUIMOD	500MG	99011-02-6
PHR2489-50MG	IMIQUIMOD RELATED COMPOUND A	50MG	99010-24-9
PHR3014-50MG	IMIQUIMOD RELATED COMPOUND B	50MG	
PHR2490-50MG	IMIQUIMOD RELATED COMPOUND C	50MG	99010-64-7
PHR3013-50MG	IMIQUIMOD RELATED COMPOUND D	50MG	
INDOMETHACIN			
PHR1247-1G	INDOMETHACIN	1G	53-86-1
PHR1247-500MG	INDOMETHACIN	500MG	53-86-1
PHR2770-50MG	INDOMETHACIN RELATED COMPOUND A	50MG	2882-15-7
PHR2581-50MG	INDOMETHACIN RELATED COMPOUND B	50MG	74-11-3
IPRATROPIUM			
PHR2796-400MG	IPRATROPIUM BROMIDE	400MG	66985-17-9
PHR3017-25MG	IPRATROPIUM BROMIDE RELATED COMPOUND A	25MG	
PHR2873-100MG	IPRATROPIUM BROMIDE RELATED COMPOUND C	100MG	
PHR2368-20MG	IPRATROPIUM IMPURITY B	20MG	
PHR2369-20MG	IPRATROPIUM IMPURITY C	20MG	

Cat. No.	Description	Pack Size	CAS No.
IRBESARTAN			
PHR1443-1G	IRBESARTAN	1G	138402-11-6
PHR1958-50MG	IRBESARTAN RELATED COMPOUND A	50MG	748812-53-5
IRINOTECAN			
PHR2717-200MG	IRINOTECAN HYDROCHLORIDE	200MG	136572-09-3
PHR2491-25MG	IRINOTECAN RELATED COMPOUND A	25MG	19685-09-7
PHR2718-20MG	IRINOTECAN RELATED COMPOUND B	20MG	86639-52-3
PHR2492-30MG	IRINOTECAN RELATED COMPOUND E	30MG	78287-27-1
ISOFLURANE			
PHR2874-1.2ML	ISOFLURANE	1.2ML	
PHR2875-0.5ML	ISOFLURANE RELATED COMPOUND A	0.5ML	
PHR2876-0.5ML	ISOFLURANE RELATED COMPOUND B	0.5ML	
KETOPROFEN			
PHR1375-500MG	KETOPROFEN	500MG	22071-15-4
PHR3018-50MG	KETOPROFEN RELATED COMPOUND C	50MG	
PHR2955-40MG	KETOPROFEN RELATED COMPOUND D	40MG	
PHR2877-50MG	KETOPROFEN RELATED COMPOUND F	50MG	
KETOROLAC			
PHR1140-500MG	KETOROLAC TROMETHAMINE	500MG	74103-07-4
PHR2493-50MG	KETOROLAC RELATED COMPOUND A	50MG	167105-80-8
PHR8716-50MG	KETOROLAC RELATED COMPOUND B	50MG	
PHR2494-50MG	KETOROLAC RELATED COMPOUND C	50MG	113502-52-6
LAMOTRIGINE			
PHR1392-1G	LAMOTRIGINE	1G	84057-84-1
PHR8729-20MG	LAMOTRIGINE RELATED COMPOUND C	20MG	
PHR8730-20MG	LAMOTRIGINE RELATED COMPOUND D	20MG	
PHR1368-100MG	LAMOTRIGINE RELATED COMPOUND B	100MG	50-45-3
LANSOPRAZOLE			
PHR1390-1G	LANSOPRAZOLE	1G	103577-45-3
PHR3019-50MG	LANSOPRAZOLE RELATED COMPOUND A	50MG	
PHR2956-100MG	LANSOPRAZOLE RELATED COMPOUND B	100MG	
LATANOPROST			
PHR1884-50MG	LATANOPROST	50MG	130209-82-4
PHR3020-100MG	LATANOPROST RELATED COMPOUND A	100MG	
PHR3022-50MG	LATANOPROST RELATED COMPOUND B	50MG	
PHR3021-50MG	LATANOPROST RELATED COMPOUND E	50MG	
LEFLUNOMIDE			
PHR1378-1G	LEFLUNOMIDE	1G	75706-12-6
PHR1568-200MG	LEFLUNOMIDE RELATED COMPOUND A	200MG	455-14-1
LEVETIRACETAM			
PHR1447-1G	LEVETIRACETAM	1G	102767-28-2
PHR2417-30MG	LEVETIRACETAM RACEMIC MIXTURE	30MG	
PHR2428-50MG	LEVETIRACETAM RELATED COMPOUND A	50MG	
PHR2432-50MG	LEVETIRACETAM RELATED COMPOUND B	50MG	7682-20-4
LEVODOPA			
PHR1271-500MG	LEVODOPA	500MG	59-92-7
PHR2095-100MG	LEVODOPA RELATED COMPOUND A	100MG	27244-64-0
PHR2096-100MG	LEVODOPA RELATED COMPOUND B	100MG	7636-26-2
LEVOFLOXACIN			
PHR1697-1G	LEVOFLOXACIN	1G	138199-71-0
PHR2229-100MG	LEVOFLOXACIN RELATED COMPOUND B	100MG	100986-89-8
PHR2230-50MG	LEVOFLOXACIN RELATED COMPOUND C	50MG	177472-30-9
PHR3023-50MG	LEVOFLOXACIN RELATED COMPOUND E	50MG	
LIDOCAINE			
PHR1034-1G	LIDOCAINE	1G	137-58-6

Cat. No.	Description	Pack Size	CAS No.
PHR1257-500MG	LIDOCAINE HYDROCHLORIDE	500MG	6108-05-0
PHR1669-100MG	LIDOCAINE RELATED COMPOUND A	100MG	201-758-7
PHR1670-100MG	LIDOCAINE RELATED COMPOUND B	100MG	
PHR1671-100MG	LIDOCAINE RELATED COMPOUND C	100MG	2198-53-0
PHR1672-100MG	LIDOCAINE RELATED COMPOUND E	100MG	
PHR1673-100MG	LIDOCAINE RELATED COMPOUND H	100MG	1131-01-7
PHR1674-100MG	LIDOCAINE RELATED COMPOUND I	100MG	
PHR1675-100MG	LIDOCAINE RELATED COMPOUND J	100MG	
PHR1977-50MG	LIDOCAINE IMPURITY D	50MG	
PHR1994-50MG	LIDOCAINE IMPURITY F	50MG	857170-72-0
PHR1995-50MG	LIDOCAINE IMPURITY G	50MG	
LINEZOLID			
PHR1885-500MG	LINEZOLID	500MG	165800-03-3
PHR3027-50MG	LINEZOLID RELATED COMPOUND A	50MG	
PHR3028-50MG	LINEZOLID RELATED COMPOUND B	50MG	
PHR3024-50MG	LINEZOLID RELATED COMPOUND C	50MG	
PHR3025-50MG	LINEZOLID RELATED COMPOUND D	50MG	
PHR3026-50MG	LINEZOLID R-ISOMER	50MG	
LISINOPRIL			
PHR1143-1G	LISINOPRIL	1G	83915-83-7
PHR2239-50MG	LISINOPRIL IMPURITY A	50MG	1012-05-1
PHR2495-100MG	LISINOPRIL IMPURITY B	100MG	6192-52-5
PHR8894-100MG	LISINOPRIL R,S,S-DIKETOPIPERAZINE	100MG	
PHR8931-100MG	LISINOPRIL S,S,S-DIKETOPIPERAZINE	100MG	
LOMUSTINE			
PHR3029-500MG	LOMUSTINE	500MG	
PHR3030-50MG	LOMUSTINE RELATED COMPOUND B	50MG	
PHR2924-500MG	LOMUSTINE RELATED COMPOUND C	500MG	2387-23-7
LORATADINE			
PHR1376-1G	LORATADINE	1G	79794-75-5
PHR2310-50MG	LORATADINE IMPURITY A	50MG	133284-74-9
PHR2311-50MG	LORATADINE IMPURITY E	50MG	170727-59-0
PHR2312-50MG	LORATADINE IMPURITY F	50MG	125743-80-8
PHR2313-100MG	LORATADINE IMPURITY H	100MG	29976-53-2
PHR2307-50MG	LORATADINE RELATED COMPOUND A	50MG	100643-71-8
PHR2308-50MG	LORATADINE RELATED COMPOUND B	50MG	38092-89-6
PHR2309-50MG	LORATADINE RELATED COMPOUND C	50MG	31251-41-9
LOVASTATIN			
PHR1285-1G	LOVASTATIN	1G	75330-75-5
PHR2240-50MG	LOVASTATIN IMPURITY A	50MG	73573-88-3
PHR2278-50MG	LOVASTATIN IMPURITY C	50MG	109273-98-5
PHR2279-20MG	LOVASTATIN IMPURITY D	20MG	149756-62-7
PHR2277-100MG	LOVASTATIN RELATED COMPOUND A	100MG	77517-29-4
MELOXICAM			
PHR1799-500MG	MELOXICAM	500MG	71125-38-7
PHR2327-50MG	MELOXICAM IMPURITY C	50MG	
PHR2131-50MG	MELOXICAM IMPURITY D	50MG	
PHR2503-50MG	MELOXICAM IMPURITY E	50MG	35511-15-0
PHR2500-50MG	MELOXICAM RELATED COMPOUND A	50MG	24683-26-9
PHR2501-50MG	MELOXICAM RELATED COMPOUND B	50MG	7305-71-7
PHR2502-50MG	MELOXICAM RELATED COMPOUND C	50MG	118854-48-1
MEMANTINE			
PHR1886-1G	MEMANTINE HYDROCHLORIDE	1G	41100-52-1
PHR2449-100MG	MEMANTINE RELATED COMPOUND A	100MG	702-79-4

Cat. No.	Description	Pack Size	CAS No.
PHR2450-50MG	MEMANTINE RELATED COMPOUND B	50MG	707-37-9
PHR2451-100MG	MEMANTINE RELATED COMPOUND C	100MG	707-36-8
PHR2583-100MG	MEMANTINE RELATED COMPOUND D	100MG	941-37-7
PHR2931-30MG	MEMANTINE RELATED COMPOUND E	30MG	351329-88-9
PHR2759-100MG	MEMANTINE RELATED COMPOUND F	100MG	356572-08-2
PHR2761-50MG	MEMANTINE RELATED COMPOUND G	50MG	33103-93-4
PHR2764-50MG	MEMANTINE RELATED COMPOUND H	50MG	15210-60-3
PHR2765-50MG	MEMANTINE RELATED COMPOUND I	50MG	80121-67-1
MEROPENEM			
PHR1772-500MG	MEROPENEM	500MG	119478-56-7
PHR2505-50MG	MEROPENEM IMPURITY A	50MG	1053703-36-8
PHR2504-50MG	MEROPENEM IMPURITY B	50MG	
MESALAMINE			
PHR1060-1G	MESALAMINE (MESALAZINE)	1G	89-57-6
PHR1225-500MG	M-AMINOPHENOL(MESALAZINE IMPURITY B)	500MG	591-27-5
PHR1739-200MG	2-AMINOPHENOL (MESALAZINE IMPURITY C)	200MG	95-55-6
MESNA			
PHR1570-500MG	MESNA	500MG	19767-45-4
PHR2732-100MG	MESNA RELATED COMPOUND A	100MG	
PHR2818-100MG	MESNA RELATED COMPOUND A	100MG	
PHR2121-100MG	MESNA RELATED COMPOUND B	100MG	
METAMIZOLE			
PHR1800-100MG	METAMIZOLE SODIUM	100MG	5907-38-0
PHR2783-100MG	METAMIZOLE IMPURITY A	100MG	1672-58-8
PHR2784-100MG	METAMIZOLE IMPURITY B	100MG	83-07-8
PHR2785-100MG	METAMIZOLE IMPURITY C	100MG	856307-27-2
PHR2786-100MG	METAMIZOLE IMPURITY D	100MG	58-15-1
METFORMINE			
PHR1084-500MG	METFORMIN HYDROCHLORIDE	500MG	1115-70-4
PHR1331-500MG	METFORMIN RELATED COMPOUND A	500MG	461-58-5
PHR1505-50MG	METFORMIN RELATED COMPOUND B	50MG	1674-62-0
PHR1506-50MG	METFORMIN RELATED COMPOUND C	50MG	1985-46-2
PHR1532-200MG	METFORMIN RELATED COMPOUND F	200MG	506-59-2
PHR1274-1G	METFORMIN IMPURITY D	1G	108-78-1
METHIMAZOLE			
PHR3037-500MG	METHIMAZOLE	500MG	
PHR3038-100MG	METHIMAZOLE RELATED COMPOUND A	100MG	
PHR3039-100MG	METHIMAZOLE RELATED COMPOUND C	100MG	
METHYLDOPA			
PHR2672-1G	METHYLDOPA	1G	41372-08-1
PHR2673-30MG	METHYLDOPA RELATED COMPOUND B	30MG	35026-10-9
PHR2674-30MG	METHYLDOPA RELATED COMPOUND C	30MG	
METHYLPREDNISOLONE			
PHR1717-500MG	METHYLPREDNISOLONE	500MG	83-43-2
PHR1718-500MG	METHYLPREDNISOLONE ACETATE	500MG	53-36-1
PHR3044-500MG	METHYLPREDNISOLONE HEMISUCCINATE	500MG	
PHR3045-100MG	METHYLPREDNISOLONE HEMISUCCINATE RELATED	100MG	
PHR3043-50MG	METHYLPREDNISOLONE RELATED COMPOUND C	50MG	
METOCLOPRAMIDE			
PHR1132-1G	METOCLOPRAMIDE HYDROCHLORIDE	1G	
PHR2507-100MG	METOCLOPRAMIDE IMPURITY C	100MG	7206-70-4
PHR2286-100MG	METOCLOPRAMIDE IMPURITY E	100MG	100-36-7
PHR2506-100MG	METOCLOPRAMIDE IMPURITY H	100MG	
PHR2283-100MG	METOCLOPRAMIDE RELATED COMPOUND A	100MG	5608-13-9

Cat. No.	Description	Pack Size	CAS No.
PHR2284-100MG	METOCLOPRAMIDE RELATED COMPOUND B	100MG	4093-31-6
PHR2285-100MG	METOCLOPRAMIDE RELATED COMPOUND D	100MG	4093-29-2
METOPROLOL			
PHR3046-500MG	METOPROLOL FUMARATE	500MG	
PHR2509-50MG	METOPROLOL RELATED COMPOUND A	50MG	109632-08-8
PHR2510-30MG	METOPROLOL RELATED COMPOUND C	30MG	
PHR2511-500MG	METOPROLOL SUCCINATE	500MG	98418-47-4
PHR1076-1G	METOPROLOL TARTRATE	1G	56392-17-7
METRONIDAZOLE			
PHR1052-1G	METRONIDAZOLE	1G	443-48-1
PHR2514-100MG	METRONIDAZOLE IMPURITY B	100MG	3034-38-6
PHR2512-100MG	METRONIDAZOLE IMPURITY E	100MG	705-19-1
PHR2165-100MG	METRONIDAZOLE IMPURITY F	100MG	16156-94-8
MINOXIDIL			
PHR2123-200MG	MINOXIDIL	200MG	38304-91-5
PHR3050-50MG	MINOXIDIL RELATED COMPOUND E	50MG	
MIRTAZAPINE			
PHR8952-500MG	MIRTAZAPINE	500MG	
PHR8736-50MG	MIRTAZAPINE RELATED COMPOUND B	50MG	
MOXIDECTIN			
PHR1827-500MG	MOXIDECTIN	500MG	113507-06-5
MOXIFLOXACIN			
PHR1542-1G	MOXIFLOXACIN HYDROCHLORIDE	1G	186826-86-8
PHR3054-50MG	MOXIFLOXACIN RELATED COMPOUND A	50MG	
PHR3055-30MG	MOXIFLOXACIN RELATED COMPOUND G	30MG	
NAPROXEN			
PHR1040-500MG	NAPROXEN	500MG	22204-53-1
PHR2219-100MG	NAPROXEN 1,2 PROPYLENE GLYCOL ESTERS	100MG	
PHR2254-50MG	NAPROXEN 1,2,3 PROPANETRIOL ESTERS	50MG	
PHR2256-50MG	NAPROXEN 1,3-BUTYLENE GLYCOL ESTERS	50MG	
PHR2248-50MG	NAPROXEN 2,3-BUTYLENE GLYCOL ESTER	50MG	
PHR2521-100MG	NAPROXEN IMPURITY A	100MG	52079-10-4
PHR2520-100MG	NAPROXEN IMPURITY B	100MG	89617-86-7
PHR2519-100MG	NAPROXEN IMPURITY C	100MG	84236-26-0
PHR2529-100MG	NAPROXEN IMPURITY E	100MG	26159-35-3
PHR2528-100MG	NAPROXEN IMPURITY F	100MG	31220-35-6
PHR2526-100MG	NAPROXEN IMPURITY H	100MG	5111-66-0
PHR2525-100MG	NAPROXEN IMPURITY M	100MG	93-04-9
PHR2524-100MG	NAPROXEN IMPURITY N	100MG	5111-65-9
PHR2257-50MG	NAPROXEN ISOPROPYL ESTER	50MG	156967-24-7
PHR2527-100MG	NAPROXEN RELATED COMPOUND A	100MG	
PHR2516-100MG	NAPROXEN RELATED COMPOUND G	100MG	23979-41-1
PHR2518-100MG	NAPROXEN RELATED COMPOUND G	100MG	23979-41-1
PHR2517-100MG	NAPROXEN RELATED COMPOUND K	100MG	77301-42-9
PHR2523-100MG	NAPROXEN RELATED COMPOUND L	100MG	3900-45-6
PHR1165-1G	NAPROXEN SODIUM	1G	26159-34-2
NEVIRAPINE			
PHR1757-1G	NEVIRAPINE	1G	129618-40-2
PHR1760-50MG	NEVIRAPINE RELATED COMPOUND A	50MG	133627-17-5
PHR1761-50MG	NEVIRAPINE RELATED COMPOUND B	50MG	287980-84-1
PHR1762-50MG	NEVIRAPINE RELATED COMPOUND C	50MG	287980-85-2
NICOTINAMIDE			
PHR2087-100MG	NICOTINAMIDE IMPURITY B	100MG	100-54-9
PHR2088-100MG	NICOTINAMIDE IMPURITY C	100MG	1452-77-3

Cat. No.	Description	Pack Size	CAS No.
PHR2089-100MG	NICOTINAMIDE IMPURITY D	100MG	1453-82-3
PHR2090-100MG	NICOTINAMIDE IMPURITY E	100MG	1986-81-8
PHR2091-100MG	NICOTINAMIDE IMPURITY F	100MG	1005-24-9
NICOTINE			
PHR2768-100MG	NICOTINE BASE	100MG	
PHR2532-500MG	NICOTINE BITARTRATE DIHYDRATE	500MG	6019-06-3
PHR3060-50MG	NICOTINE RELATED COMPOUND A/EP IMPURITY	50MG	
PHR3061-50MG	NICOTINE RELATED COMPOUND B	50MG	
PHR2533-50MG	NICOTINE RELATED COMPOUND C	50MG	486-56-6
PHR3062-50MG	NICOTINE RELATED COMPOUND D	50MG	
PHR3063-50MG	NICOTINE RELATED COMPOUND F	50MG	
PHR3064-50MG	NICOTINE RELATED COMPOUND G	50MG	
NILOTINIB			
PHR2833-500MG	NILOTINIB	500MG	923288-90-8
PHR2834-50MG	NILOTINIB RELATED COMPOUND A	50MG	641571-11-1
PHR2835-50MG	NILOTINIB RELATED COMPOUND C	50MG	2458-12-0
NIMODIPINE			
PHR1293-1G	NIMODIPINE	1G	66085-59-4
PHR3065-100MG	NIMODIPINE RELATED COMPOUND A	100MG	
NITAZOXANIDE			
PHR2888-400MG	NITAZOXANIDE	400MG	
PHR2889-50MG	NITAZOXANIDE RELATED COMPOUND A	50MG	
NORFLOXACIN			
PHR1614-1G	NORFLOXACIN	1G	70458-96-7
PHR2233-100MG	NORFLOXACIN RELATED COMPOUND A	100MG	68077-26-9
OFLOXACIN			
PHR1168-1G	OFLOXACIN	1G	82419-36-1
PHR2234-100MG	OFLOXACIN IMPURITY A	100MG	82419-35-0
PHR2235-50MG	OFLOXACIN RELATED COMPOUND A	50MG	82419-52-1
OLANZAPINE			
PHR1825-1G	OLANZAPINE	1G	132539-06-1
PHR2536-30MG	OLANZAPINE RELATED COMPOUND A	30MG	138564-59-7
PHR2537-30MG	OLANZAPINE RELATED COMPOUND B	30MG	221176-49-4
PHR2538-30MG	OLANZAPINE RELATED COMPOUND C	30MG	174794-02-6
OLMESARTAN			
PHR1851-200MG	OLMESARTAN MEDOXOMIL	200MG	144689-63-4
PHR2016-40MG	OLMESARTAN MEDOXOMIL RELATED COMPOUND A	40MG	849206-43-5
OMEPRAZOLE			
PHR1059-1G	OMEPRAZOLE	1G	73590-58-6
PHR1987-30MG	OMEPRAZOLE IMPURITY B	30MG	110374-16-8
PHR1988-30MG	OMEPRAZOLE IMPURITY C	30MG	73590-85-9
PHR1989-30MG	OMEPRAZOLE IMPURITY H	30MG	863029-89-4
PHR1648-30MG	OMEPRAZOLE RELATED COMPOUND A	30MG	88546-55-8
PHR1649-30MG	OMEPRAZOLE RELATED COMPOUND E	30MG	176219-04-8
PHR1650-30MG	OMEPRAZOLE RELATED COMPOUND I	30MG	
PHR2329-50MG	OMEPRAZOLE RELATED COMPOUND F AND G	50MG	
ONDANSETRON			
PHR1141-1G	ONDANSETRON HYDROCHLORIDE	1G	103639-04-9
PHR2223-100MG	ONDANSETRON RELATED COMPOUND A	100MG	119812-29-2
PHR2224-100MG	ONDANSETRON RELATED COMPOUND C	100MG	27387-31-1
PHR2221-100MG	ONDANSETRON RELATED COMPOUND D	100MG	99614-64-9
PHR1181-500MG	ONDANSETRON IMPURITY F	500MG	693-98-1
ORLISTAT			
PHR1445-1G	ORLISTAT	1G	96829-58-2

Cat. No.	Description	Pack Size	CAS No.
PHR2543-50MG	ORLISTAT RELATED COMPOUND B	50MG	19740-72-8
PHR2544-50MG	ORLISTAT RELATED COMPOUND C	50MG	791-28-6
ORPHENADRINE			
PHR1777-1G	ORPHENADRINE CITRATE	1G	4682-36-4
PHR2957-40MG	ORPHENADRINE RELATED COMPOUND B	40MG	
PHR2958-40MG	ORPHENADRINE RELATED COMPOUND C	40MG	
PHR2959-100MG	ORPHENADRINE RELATED COMPOUND E	100MG	
PHR3070-100MG	ORPHENADRINE RELATED COMPOUND F	100MG	
OXALIPLATIN			
PHR1528-200MG	OXALIPLATIN	200MG	61825-94-3
PHR2545-200MG	OXALIPLATIN RELATED COMPOUND A	200MG	6153-56-6
PHR2546-50MG	OXALIPLATIN RELATED COMPOUND F	50MG	66845-32-7
OXCARBAZEPINE			
PHR1635-1G	OXCARBAZEPINE	1G	28721-07-5
PHR3071-25MG	OXCARBAZEPINE RELATED COMPOUND C	25MG	
PHR2238-50MG	OXCARBAZEPINE RELATED COMPOUND E	50MG	21737-58-6
OXYBUTYNIN			
PHR2547-500MG	OXYBUTYNIN CHLORIDE	500MG	1508-65-2
PHR2960-50MG	OXYBUTYNIN IMPURITY A	50MG	
PHR2548-200MG	OXYBUTYNIN RELATED COMPOUND A	200MG	4335-77-7
PHR2549-50MG	OXYBUTYNIN RELATED COMPOUND B	50MG	10399-13-0
OXYCODONE			
PHR8802-200MG	OXYCODONE	200MG	
PHR8800-25MG	OXYCODONE RELATED COMPOUND A	25MG	
PHR8801-30MG	OXYCODONE RELATED COMPOUND B	30MG	
PHR8999-30MG	OXYCODONE RELATED COMPOUND C	30MG	
OXYMETAZOLIN			
PHR1572-500MG	OXYMETAZOLINE HYDROCHLORIDE	500MG	2315-02-8
PHR2215-100MG	OXYMETAZOLINE RELATED SUBSTANCE IV	100MG	1391053-50-1
PHR2215-50MG	OXYMETAZOLINE RELATED SUBSTANCE IV	50MG	1391053-50-1
PACLITAXEL			
PHR1803-200MG	PACLITAXEL	200MG	33069-62-4
PHR2113-50MG	PACLITAXEL IMPURITY C	50MG	153415-45-3
PHR2111-50MG	PACLITAXEL RELATED COMPOUND A	50MG	71610-00-9
PHR2112-50MG	PACLITAXEL RELATED COMPOUND B	50MG	78454-17-8
PANTOPRAZOLE			
PHR1604-1G	PANTOPRAZOLE SODIUM	1G	164579-32-2
PHR2458-50MG	PANTOPRAZOLE RELATED COMPOUND A	50MG	127780-16-9
PHR3074-50MG	PANTOPRAZOLE RELATED COMPOUND B	50MG	
PHR3075-50MG	PANTOPRAZOLE RELATED COMPOUND C	50MG	
PAROXETINE			
PHR1804-500MG	PAROXETINE HYDROCHLORIDE	500MG	110429-35-1
PHR3076-50MG	PAROXETINE RELATED COMPOUND C	50MG	
PHR3077-50MG	PAROXETINE RELATED COMPOUND F	50MG	
PERINDOPRI			
PHR1887-500MG	PERINDOPRIL ERBUMINE	500MG	107133-36-8
PHR3078-50MG	PERINDOPRIL RELATED COMPOUND A	50MG	
PHENYL BUTYRATE			
PHR3080-50MG	PHENYL BUTYRATE RELATED COMPOUND A	50MG	
PHR3081-200MG	PHENYL BUTYRATE RELATED COMPOUND B	200MG	
PHR3082-50MG	PHENYL BUTYRATE RELATED COMPOUND C	50MG	
PHENYLEPHRINE			
PHR1017-500MG	PHENYLEPHRINE HYDROCHLORIDE	500MG	61-76-7
PHR2151-100MG	PHENYLEPHRINE IMPURITY 01	100MG	57196-61-9

Cat. No.	Description	Pack Size	CAS No.
PHR2152-50MG	PHENYLEPHRINE IMPURITY 02	50MG	
PHR2153-100MG	PHENYLEPHRINE IMPURITY 03	100MG	
PHR2154-100MG	PHENYLEPHRINE IMPURITY 04	100MG	
PHR1695-50MG	PHENYLEPHRINE RELATED COMPOUND C	50MG	94240-17-2
PHR2148-100MG	PHENYLEPHRINE RELATED COMPOUND D	100MG	1367567-95-0
PHR1712-50MG	PHENYLEPHRINE RELATED COMPOUND E	50MG	71786-67-9
PHR2149-100MG	PHENYLEPHRINE RELATED COMPOUND F	100MG	
PHR2150-100MG	PHENYLEPHRINE RELATED COMPOUND G	100MG	1094089-46-9
PHENYTOIN			
PHR1139-1G	PHENYTOIN	1G	57-41-0
PHR1574-100MG	PHENYTOIN RELATED COMPOUND A	100MG	3060-50-2
PHR1692-100MG	PHENYTOIN RELATED COMPOUND B	100MG	6802-95-5
PHR2554-100MG	PHENYTOIN IMPURITY D	100MG	
PHR2553-100MG	PHENYTOIN IMPURITY F	100MG	
PHR2246-100MG	PHENYTOIN INTERMEDIATE COMPOUND M01	100MG	6859-11-6
PHR1492-1G	PHENYTOIN SODIUM	1G	630-93-3
PIROCTONE			
PHR2187-500MG	PIROCTONE OLAMINE	500MG	68890-66-4
PHR3086-50MG	PIROCTONE RELATED COMPOUND B	50MG	
PRAMIPEXOLE			
PHR1598-500MG	PRAMIPEXOLE DIHYDROCHLORIDE	500MG	191217-81-9
PHR2010-50MG	PRAMIPEXOLE RELATED COMPOUND A	50MG	106092-09-5
PHR1824-50MG	PRAMIPEXOLE RELATED COMPOUND D	50MG	104632-27-1
PRAZIQUANTEL			
PHR1391-1G	PRAZIQUANTEL	1G	55268-74-1
PHR2961-40MG	PRAZIQUANTEL RELATED COMPOUND A	40MG	
PHR2962-100MG	PRAZIQUANTEL RELATED COMPOUND B	100MG	
PHR2963-100MG	PRAZIQUANTEL RELATED COMPOUND C	100MG	
PREDNISOLONE			
PHR1043-500MG	PREDNISOLONE	500MG	50-24-8
PHR1630-1G	PREDNISOLONE ACETATE	1G	52-21-1
PHR2556-50MG	PREDNISOLONE IMPURITY H	50MG	2427-64-7
PHR2557-50MG	PREDNISOLONE IMPURITY I	50MG	13479-38-4
PHR2737-1G	PREDNISOLONE SODIUM PHOSPHATE	1G	
PHR2816-500MG	PREDNISOLONE SODIUM PHOSPHATE	500MG	125-02-0
PREDNISON			
PHR1042-1G	PREDNISON	1G	53-03-2
PHR2288-20MG	PREDNISON IMPURITY C	20MG	70522-55-3
PHR2289-50MG	PREDNISON IMPURITY J	50MG	1807-14-3
PRIMIDONE			
PHR3109-500MG	PRIMIDONE	500MG	
PHR3110-50MG	PRIMIDONE RELATED COMPOUND A	50MG	
PHR3111-50MG	PRIMIDONE RELATED COMPOUND C	50MG	
PROGESTERONE			
PHR1142-1G	PROGESTERONE	1G	57-83-0
PHR2280-50MG	PROGESTERONE IMPURITY J	50MG	1162-54-5
PHR2904-100MG	PROGESTERONE IMPURITY K	100MG	
PHR2905-100MG	PROGESTERONE RELATED COMPOUND H	100MG	
PROGUANIL			
PHR1713-1G	PROGUANIL HYDROCHLORIDE	1G	637-32-1
PHR3115-50MG	PROGUANIL RELATED COMPOUND A	50MG	
PHR3116-50MG	PROGUANIL RELATED COMPOUND C	50MG	
PHR3117-50MG	PROGUANIL RELATED COMPOUND E	50MG	

Cat. No.	Description	Pack Size	CAS No.
PROMETHAZINE			
PHR1467-1G	PROMETHAZINE HYDROCHLORIDE	1G	58-33-3
PHR2777-50MG	PROMETHAZINE RELATED COMPOUND B	50MG	5568-90-1
PHR3563-200MG	PROMETHAZINE SULPHOXIDE	200MG	
PROPOFOL			
PHR1663-1ML	PROPOFOL	1ML	2078-54-8
PHR2565-50MG	PROPOFOL RELATED COMPOUND A	50MG	2416-95-7
PHR2566-50MG	PROPOFOL RELATED COMPOUND B	50MG	1988-11-0
PHR3118-50MG	PROPOFOL RELATED COMPOUND C/EP IMPURITY&	50MG	
PROPRANOLOL			
PHR1308-1G	PROPRANOLOL HYDROCHLORIDE	1G	318-98-9
PHR3119-50MG	PROPRANOLOL RELATED COMPOUND A	50MG	
PYRANTEL			
PHR2080-1G	PYRANTEL PAMOATE	1G	22204-24-6
PHR2082-50MG	PYRANTEL IMPURITY B	50MG	36700-39-7
PHR2083-50MG	PYRANTEL IMPURITY C	50MG	98-03-3
PHR2084-50MG	PYRANTEL IMPURITY D	50MG	4271-96-9
PHR2081-50MG	PYRANTEL RELATED COMPOUND A	50MG	
PYRIDOXINE			
PHR1036-500MG	PYRIDOXINE HYDROCHLORIDE	500MG	58-56-0
PHR1693-50MG	PYRIDOXINE IMPURITY A	50MG	5196-20-3
PHR1694-50MG	PYRIDOXINE IMPURITY B	50MG	
QUETIAPINE			
PHR1856-1G	QUETIAPINE FUMARATE	1G	111974-72-2
PHR3121-40MG	QUETIAPINE RELATED COMPOUND B	40MG	
PHR3122-40MG	QUETIAPINE RELATED COMPOUND G	40MG	
RABEPRAZOLE			
PHR1525-1G	RABEPRAZOLE SODIUM	1G	117976-90-6
PHR3123-50MG	RABEPRAZOLE RELATED COMPOUND C	50MG	
PHR3124-50MG	RABEPRAZOLE RELATED COMPOUND E	50MG	
PHR3125-30MG	RABEPRAZOLE RELATED COMPOUND F	30MG	
RALOXIFENE			
PHR1852-500MG	RALOXIFENE HYDROCHLORIDE	500MG	82640-04-8
PHR2015-40MG	RALOXIFENE RELATED COMPOUND C	40MG	195454-31-0
RAMIPRIL			
PHR1446-1G	RAMIPRIL	1G	87333-19-5
PHR3128-50MG	RAMIPRIL RELATED COMPOUND A	50MG	
RANITIDINE			
PHR1026-500MG	RANITIDINE HYDROCHLORIDE	500MG	66357-59-3
PHR1773-100MG	RANITIDINE RELATED COMPOUND A	100MG	91224-69-0
PHR1774-100MG	RANITIDINE RELATED COMPOUND B	100MG	72126-78-4
PHR1775-100MG	RANITIDINE RELATED COMPOUND C	100MG	73851-70-4
PHR2175-100MG	RANITIDINE IMPURITY D	100MG	112251-56-6
PHR2176-100MG	RANITIDINE IMPURITY E	100MG	
PHR2177-100MG	RANITIDINE IMPURITY F	100MG	
PHR2178-100MG	RANITIDINE IMPURITY G	100MG	112233-23-5
PHR2180-100MG	RANITIDINE IMPURITY I	100MG	
PHR2181-100MG	RANITIDINE IMPURITY J	100MG	1331637-48-9
PHR2182-100MG	RANITIDINE IMPURITY K	100MG	61832-41-5
RANOLAZINE			
PHR2676-500MG	RANOLAZINE	500MG	95635-55-5
PHR2677-100MG	RANOLAZINE RELATED COMPOUND A	100MG	2210-74-4
PHR2678-100MG	RANOLAZINE RELATED COMPOUND B	100MG	
PHR2679-100MG	RANOLAZINE RELATED COMPOUND C	100MG	5294-61-1
PHR2680-100MG	RANOLAZINE RELATED COMPOUND D	100MG	380204-72-8

Cat. No.	Description	Pack Size	CAS No.
REPAGLINIDE			
PHR3408-400MG	REPAGLINIDE	400MG	
PHR3129-25MG	REPAGLINIDE RELATED COMPOUND A	25MG	
PHR3130-50MG	REPAGLINIDE RELATED COMPOUND C	50MG	
PHR3131-50MG	REPAGLINIDE RELATED COMPOUND E	50MG	
RETINYL			
PHR1236-1G	RETINYL ACETATE	1G	127-47-9
PHR1235-1G	RETINYL PALMITATE	1G	79-81-2
PHR2569-25MG	RETINYL PALMITATE IMPURITY A	25MG	
PHR2567-25MG	RETINYL PALMITATE IMPURITY B	25MG	1224-78-8
PHR2568-25MG	RETINYL PALMITATE IMPURITY C	25MG	16729-22-9
PHR1668-1G	RETINYL PROPIONATE	1G	7069-42-3
RIBAVIRIN			
PHR1700-1G	RIBAVIRIN	1G	36791-04-5
PHR2778-30MG	RIBAVIRIN RELATED COMPOUND A	30MG	39925-19-4
PHR2779-30MG	RIBAVIRIN RELATED COMPOUND D	30MG	3641-08-5
RILUZOLE			
PHR2499-300MG	RILUZOLE	300MG	1744-22-5
PHR3132-50MG	RILUZOLE RELATED COMPOUND A	50MG	
RISEDRONATE			
PHR1888-1G	RISEDRONATE SODIUM	1G	329003-65-8
PHR2078-50MG	RISEDRONATE IMPURITY C	50MG	105462-25-7
PHR2079-50MG	RISEDRONATE IMPURITY D	50MG	6419-36-9
PHR2075-50MG	RISEDRONATE RELATED COMPOUND A	50MG	105462-23-5
PHR2076-50MG	RISEDRONATE RELATED COMPOUND B	50MG	
PHR2077-50MG	RISEDRONATE RELATED COMPOUND C	50MG	75755-10-1
RIVASTIGMINE			
PHR1867-1G	RIVASTIGMINE TARTRATE	1G	129101-54-8
PHR2029-30MG	RIVASTIGMINE RELATED COMPOUND A	30MG	32634-68-7
PHR2030-30MG	RIVASTIGMINE RELATED COMPOUND B	30MG	25081-93-0
PHR2031-30MG	RIVASTIGMINE RELATED COMPOUND C	30MG	139306-10-8
PHR2032-100MG	RIVASTIGMINE RELATED COMPOUND D	100MG	855300-09-3
PHR2033-30MG	RIVASTIGMINE RELATED COMPOUND E	30MG	
PHR2034-50MG	RIVASTIGMINE RELATED COMPOUND F	50MG	1346602-84-3
ROPINIROLE			
PHR1633-200MG	ROPINIROLE HYDROCHLORIDE	200MG	91374-20-8
PHR3133-100MG	ROPINIROLE RELATED COMPOUND B	100MG	
ROPIVACAINE			
PHR2571-500MG	ROPIVACAINE	500MG	84057-95-4
PHR2303-500MG	ROPIVACAINE HYDROCHLORIDE	500MG	132112-35-7
PHR2572-50MG	ROPIVACAINE RELATED COMPOUND A	50MG	21436-98-6
PHR2306-100MG	ROPIVACAINE RELATED COMPOUND B/ROPIVACAI	100MG	
SAFINAMIDE			
PHR2379-50MG	SAFINAMIDE MESYLATE	50MG	202825-46-5
PHR2381-20MG	SAFINAMIDE DEFLUOR-DERIVATIVE	20MG	
PHR2383-20MG	SAFINAMIDE DIALANINE AMIDE ADDUCT	20MG	
PHR2380-20MG	SAFINAMIDE FREE ACID	20MG	
PHR2382-20MG	SAFINAMIDE IMIDAZOLIDINONE	20MG	
SALBUTAMOL			
PHR2738-500MG	SALBUTAMOL BASE	500MG	
PHR2378-10MG	SALBUTAMOL ETHYL ETHER	10MG	
PHR1960-50MG	SALBUTAMOL IMPURITY D	50MG	156339-88-7
PHR1961-50MG	SALBUTAMOL IMPURITY G	50MG	
PHR1962-20MG	SALBUTAMOL IMPURITY I	20MG	56796-66-8

Cat. No.	Description	Pack Size	CAS No.
PHR2939-50MG	SALBUTAMOL IMPURITY K	50MG	
PHR2936-50MG	SALBUTAMOL IMPURITY L	50MG	
PHR2938-50MG	SALBUTAMOL IMPURITY N	50MG	
PHR2395-1MG	SALBUTAMOL SIDE-BY-SIDE DIMER	1MG	
PHR2374-1ML	SALBUTAMOL SULPHATE TEST MIX 1	1ML	
PHR2375-1ML	SALBUTAMOL TEST MIX 2	1ML	
PHR2376-1MG	SALBUTAMOL TEST MIX 3	1MG	
SALICYLIC ACID			
PHR1013-1G	SALICYLIC ACID	1G	69-72-7
PHR1643-200MG	SALICYLIC ACID RELATED COMPOUND B	200MG	636-46-4
SESAME OIL			
PHR2967-1G	SESAME OIL	1G	
PHR2966-50MG	SESAME OIL RELATED COMPOUND A	50MG	
PHR2968-50MG	SESAME OIL RELATED COMPOUND B	50MG	
SPIRONOLACTONE			
PHR2886-200MG	SPIRONOLACTONE	200MG	
PHR2887-150MG	SPIRONOLACTONE RELATED COMPOUND A	150MG	
PHR3135-50MG	SPIRONOLACTONE RELATED COMPOUND C	50MG	
SULFAMETHOXAZOLE			
PHR1126-1G	SULFAMETHOXAZOLE	1G	723-46-6
PHR2578-30MG	SULFAMETHOXAZOLE IMPURITY A	30MG	21312-10-7
PHR2577-100MG	SULFAMETHOXAZOLE RELATED COMPOUND C	100MG	1072-67-9
SUMATRIPTAN			
PHR2579-100MG	SUMATRIPTAN	100MG	103628-46-2
PHR2579-500MG	SUMATRIPTAN	500MG	103628-46-2
PHR2580-100MG	SUMATRIPTAN SUCCINATE	100MG	103628-48-4
PHR2580-500MG	SUMATRIPTAN SUCCINATE	500MG	103628-48-4
PHR2913-40MG	SUMATRIPTAN SUCCINATE RELATED COMPOUND A	40MG	
PHR2914-75MG	SUMATRIPTAN SUCCINATE RELATED COMPOUND C	75MG	
TERAZOSIN			
PHR1721-1G	TERAZOSIN HYDROCHLORIDE	1G	70024-40-7
PHR3137-100MG	TERAZOSIN RELATED COMPOUND A	100MG	
PHR3138-100MG	TERAZOSIN RELATED COMPOUND B	100MG	
PHR3139-50MG	TERAZOSIN RELATED COMPOUND C	50MG	
TERBINAFINE			
PHR1298-1G	TERBINAFINE HYDROCHLORIDE	1G	78628-80-5
PHR2885-50MG	TERBINAFINE RELATED COMPOUND A	50MG	
PHR2927-50MG	TERBINAFINE RELATED COMPOUND B	50MG	176168-78-8
TETRACAINE			
PHR1129-1G	TETRACAINE HYDROCHLORIDE	1G	136-47-0
PHR2236-50MG	TETRACAINE RELATED COMPOUND B	50MG	4740-24-3
PHR2261-100MG	TETRACAINE RELATED COMPOUND C	100MG	71839-12-8
PHR3140-100MG	TETRACAINE	100MG	
THEOPHYLLINE			
PHR1023-1G	THEOPHYLLINE	1G	58-55-9
PHR2587-100MG	THEOPHYLLINE IMPURITY B	100MG	1076-22-8
PHR1151-1G	THEOPHYLLINE MELTING POINT STANDARD	1G	58-55-9
PHR2542-100MG	THEOPHYLLINE RELATED COMPOUND B	100MG	1076-22-8
PHR2550-50MG	THEOPHYLLINE RELATED COMPOUND C	50MG	7597-60-6
PHR2551-50MG	THEOPHYLLINE RELATED COMPOUND D	50MG	116131-08-9
PHR2541-100MG	THEOPHYLLINE RELATED COMPOUND F	100MG	519-37-9
PHR2588-100MG	THEOPHYLLINE RELATED COMPOUND F	100MG	519-37-9
THIAMINE			
PHR1037-1G	THIAMINE HYDROCHLORIDE	1G	67-03-8

Cat. No.	Description	Pack Size	CAS No.
PHR2216-1G	THIAMINE DISULFIDE	1G	67-16-3
PHR2216-50MG	THIAMINE DISULFIDE	50MG	67-16-3
PHR1996-20MG	THIAMINE IMPURITY A	20MG	
PHR1997-20MG	THIAMINE IMPURITY B	20MG	
PHR1998-20MG	THIAMINE IMPURITY C	20MG	
PHR2168-100MG	THIAMINE IMPURITY D	100MG	490-82-4
PHR1759-20MG	THIAMINE IMPURITY E	20MG	299-35-4
PHR2170-100MG	THIAMINE IMPURITY G	100MG	3419-28-1
PHR2171-50MG	THIAMINE IMPURITY H	50MG	89285-03-0
PHR2172-50MG	THIAMINE IMPURITY I	50MG	13100-57-7
PHR1337-1G	THIAMINE NITRATE	1G	532-43-4
PHR1369-500MG	THIAMINE PYROPHOSPHATE	500MG	154-87-0
PHR2207-100MG	THIAMINE THIOCARBAMATE	100MG	89285-03-0
THIMEROSAL			
PHR1587-1G	THIMEROSAL	1G	54-64-8
PHR2590-100MG	THIMEROSAL RELATED COMPOUND A	100MG	119-80-2
TIMOLOL			
PHR2593-500MG	TIMOLOL MALEATE	500MG	26921-17-5
PHR3144-50MG	TIMOLOL RELATED COMPOUND A	50MG	
PHR2594-50MG	TIMOLOL RELATED COMPOUND F	50MG	30165-96-9
TRIAMTERENE			
PHR1722-400MG	TRIAMTERENE	400MG	396-01-1
PHR2014-100MG	TRIAMTERENE RELATED COMPOUND A	100MG	1006-23-1
PHR1819-100MG	TRIAMTERENE RELATED COMPOUND B	100MG	19375-89-4
PHR1820-100MG	TRIAMTERENE RELATED COMPOUND C	100MG	19152-93-3
VALACYCLOVIR			
PHR1601-1G	VALACYCLOVIR HYDROCHLORIDE	1G	124832-27-5
PHR2925-500MG	VALACYCLOVIR RELATED COMPOUND G	500MG	1122-58-3
VALPROIC ACID			
PHR1061-1G	VALPROIC ACID	1G	99-66-1
PHR1652-0.5ML	VALPROIC ACID RELATED COMPOUND A	0.5ML	99-67-2
PHR2281-100MG	VALPROIC ACID RELATED COMPOUND B	100MG	62391-99-5
PHR2614-100MG	VALPROIC ACID IMPURITY H	100MG	110-59-8
PHR2282-100MG	VALPROIC ACID IMPURITY K	100MG	5343-52-2
VALSARTAN			
PHR1315-1G	VALSARTAN	1G	137862-53-4
PHR1875-30MG	VALSARTAN RELATED COMPOUND A	30MG	
PHR1876-30MG	VALSARTAN RELATED COMPOUND B	30MG	952652-79-8
PHR1877-30MG	VALSARTAN RELATED COMPOUND C	30MG	
PHR1878-30MG	VALSARTAN RELATED COMPOUND E	30MG	137863-17-3
VERAPAMIL			
PHR1131-1G	VERAPAMIL HYDROCHLORIDE	1G	152-11-4
PHR2616-100MG	VERAPAMIL IMPURITY K	100MG	20850-49-1
PHR8890-100MG	VERAPAMIL IMPURITY P	100MG	2086275-13-8
PHR2360-100MG	VERAPAMIL RELATED COMPOUND A	100MG	67775-97-7
PHR2361-100MG	VERAPAMIL RELATED COMPOUND B	100MG	1794-55-4
PHR1194-500MG	VERAPAMIL RELATED COMPOUND E	500MG	120-14-9
PHR1193-500MG	VERAPAMIL RELATED COMPOUND F	500MG	
VIGABATRIN			
PHR3158-100MG	VIGABATRIN RELATED COMPOUND A	100MG	
PHR3159-50MG	VIGABATRIN RELATED COMPOUND E	50MG	
ZALEPLON			
PHR8816-500MG	ZALEPLON	500MG	
PHR3162-30MG	ZALEPLON RELATED COMPOUND A	30MG	

Cat. No.	Description	Pack Size	CAS No.
ZIDOVUDINE			
PHR1292-1G	ZIDOVUDINE	1G	30516-87-1
PHR1345-250MG	ZIDOVUDINE RELATED COMPOUND C	250MG	65-71-4
PHR3163-100MG	ZIDOVUDINE RELATED COMPOUND D	100MG	50-89-5
ZILEUTON			
PHR2555-200MG	ZILEUTON	200MG	111406-87-2
PHR3164-100MG	ZILEUTON RELATED COMPOUND B	100MG	
PHR2724-500MG	ZINC OXIDE	500MG	1314-13-2
ZOLEDRONIC ACID			
PHR1893-1G	ZOLEDRONIC ACID	1G	165800-06-6
PHR2459-50MG	ZOLEDRONIC ACID RELATED COMPOUND A	50MG	22884-10-2

Cat. No.	Description	Pack Size	CAS No.
New Products			
PHR1963-50MG	(+)-CATECHIN	50MG	154-23-4
PHR2197-100MG	(2R)-2-AMINO-2-(4-HYDROXYPHENYL)-ACETAMIDE	100MG	54397-23-8
PHR2104-100MG	2-(4-CHLOROPHENYL)-4-DIMETHYLAMINO BUTYRONITRILE HYDROCHLORIDE	100MG	854863-55-1
PHR2146-100MG	2-(4-ISOBUTYLPHENYL)-N-(1-PHENYLETHYL)PROPIONAMIDE	100MG	105959-56-6
PHR2141-50MG	2-(4-ISOBUTYLPHENYL)PROP-2-ENOIC ACID	50MG	
PHR2122-50MG	2-(4-ISOBUTYRYLPHENYL)PROPANAMIDE	50MG	
PHR2147-50MG	2-(4-SEC-BUTYLPHENYL)PROP-2-ENOIC ACID	50MG	1368654-26-5
PHR1915-3X1.2ML	2,2,4-TRIMETHYLPENTANE	3X1.2ML	540-84-1
PHR1719-1G	2-AMINO-6-METHOXYPURINE	1G	20535-83-5
PHR1914-3X1.2ML	2-ETHYLHEXANOIC ACID	3X1.2ML	149-57-5
PHR2136-100MG	4-(1-HYDROXY-2-METHYLPROPYL)-ACETOPHENONE	100MG	
PHR2117-200MG	4,4'-DIFLUOROBENZOPHENONE	200MG	345-92-6
PHR2212-100MG	4-AMINO-5-AMINOMETHYL-2-METHYLPYRIMIDINE DIHYDROCHLORIDE	100MG	874-43-1
PHR1745-500MG	4-AMINOPHENOL HYDROCHLORIDE	500MG	51-78-5
PHR1916-3X1.2ML	4-HYDROXY-4-METHYL-2-PENTANONE	3X1.2ML	123-42-2
PHR2124-100MG	4-ISOBUTYRYLACETOPHENONE	100MG	103931-20-0
PHR1546-500MG	ACETYLCHOLINE CHLORIDE	500MG	60-31-1
PHR1942-50MG	ACETYLSALICYLIC ANHYDRIDE	50MG	1466-82-6
PHR1940-50MG	ACETYLSALICYLSALICYLIC ACID	50MG	530-75-6
PHR1599-1G	ALENDRONATE SODIUM	1G	121268-17-5
PHR2063-300MG	AMBROXOL HYDROCHLORIDE	300MG	23828-92-4
PHR1654-1G	AMIKACIN	1G	37517-28-5
PHR1860-1G	AMIKACIN SULFATE	1G	39831-55-5
PHR1840-500MG	AMINO BENZOIC ACID	500MG	150-13-0
PHR2198-75MG	APIGENIN	75MG	520-36-5
PHR2199-75MG	APIGENIN-7-GLUCOSIDE	75MG	578-74-5
PHR2100-500MG	ARABINOSE	500MG	10323-20-3
PHR1909-500MG	ATENOLOL	500MG	29122-68-7
PHR1911-300MG	AZACITIDINE	300MG	320-67-2
PHR1785-500MG	AZTREONAM	500MG	78110-38-0
PHR2097-500MG	BETAHISTINE HYDROCHLORIDE	500MG	5579-84-0
PHR1913-500MG	BISOCTRIZOLE	500MG	103597-45-1
PHR1831-200MG	BROMHEXINE HYDROCHLORIDE	200MG	611-75-6
PHR1861-50MG	BROMODIPHENHYDRAMINE HYDROCHLORIDE	50MG	1808-12-4
PHR1864-500MG	BROMPHENIRAMINE MALEATE	500MG	980-71-2
PHR1716-500MG	CAPREOMYCIN SULFATE	500MG	1405-37-4
PHR1684-1G	CARBIMAZOLE	1G	22232-54-8
PHR2013-50MG	CARBOXYMETHYL-CYCLOHEXANECARBOXYLIC ACID	50MG	67950-95-2
PHR2208-20MG	CARNOSIC ACID	20MG	3650-09-7
PHR2099-1G	CEFOTAXIME SODIUM	1G	64485-93-4
PHR1847-500MG	CEFTAZIDIME PENTAHYDRATE	500MG	78439-06-2
PHR1848-1G	CEPHALEXIN	1G	15686-71-2
PHR2202-100MG	CHLOROGENIC ACID	100MG	327-97-9
PHR1919-500MG	CHLORZOXAZONE	500MG	95-25-0
PHR1786-500MG	CHONDROITIN SULFATE SODIUM	500MG	9082-07-9
PHR1921-500MG	CLOBETASOL PROPIONATE	500MG	25122-46-7
PHR1720-1G	CLOMIPRAMINE HYDROCHLORIDE	1G	17321-77-6
PHR1922-500MG	CLOXACILLIN SODIUM	500MG	7081-44-9
PHR1764-400MG	COLCHICINE	400MG	64-86-8
PHR1737-500MG	COLISTIMETHATE SODIUM	500MG	8068-28-8
PHR1605-1G	COLISTIN SULFATE	1G	1264-72-8
PHR1461-1G	COPPER CHLOROPHYLLIN	1G	11006-34-1
PHR2155-500MG	CORTISONE	500MG	53-06-5

Cat. No.	Description	Pack Size	CAS No.
PHR1923-300MG	CORTISONE ACETATE	300MG	50-04-4
PHR2209-50MG	CURCUMIN	50MG	458-37-7
PHR1765-1G	CYCLOSERINE	1G	68-41-7
PHR1787-500MG	CYTARABINE	500MG	147-94-4
PHR1924-250MG	DACARBAZINE	250MG	4342-03-4
PHR1858-100MG	DAIDZEIN	100MG	486-66-8
PHR1735-1G	DEMECLOCYCLINE HYDROCHLORIDE	1G	64-73-3
PHR2025-50MG	DESBENZYL DONEPEZIL HYDROCHLORIDE	50MG	120013-39-0
PHR1830-100MG	DESOGESTREL	100MG	54024-22-5
PHR1956-200MG	D-GLUCURONIC ACID	200MG	6556-12-3
PHR2094-5G	DIETHYLAMINO HYDROXYBENZOYL HEXYL BENZOATE	5G	302776-68-7
PHR1771-300MG	DIGOXIN	300MG	20830-75-5
PHR1517-500MG	DIHYDROSTREPTOMYCIN SULFATE	500MG	5490-27-7
PHR2118-500MG	DIMETHYL FUMARATE	500MG	624-49-7
PHR1690-1G	DISULFIRAM	1G	97-77-8
PHR1883-200MG	DOCETAXEL	200MG	148408-66-6
PHR1788-1G	DOCOSAHEXAENOIC ACID ETHYL ESTER	1G	84494-72-4
PHR1925-500MG	DORZOLAMIDE HYDROCHLORIDE	500MG	130693-82-2
PHR1789-200MG	DOXORUBICIN HYDROCHLORIDE	200MG	25316-40-9
PHR1849-500MG	DYCLONINE HYDROCHLORIDE	500MG	536-43-6
PHR1790-1G	EICOSAPENTAENOIC ACID ETHYL ESTER	1G	86227-47-6
PHR2120-500MG	EMTRICITABINE	500MG	143491-57-0
PHR1766-1G	EPRINOMECTIN	1G	123997-26-2
PHR1453-1G	ERYTHRULOSE	1G	533-50-6
PHR1733-1G	ESCITALOPRAM OXALATE	1G	219861-08-2
PHR1930-500MG	ETHAMBUTOL HYDROCHLORIDE	500MG	1070-11-7
PHR1374-1G	ETHYLCELLULOSE	1G	
PHR2169-100MG	ETHYLTHIAMINE HYDROCHLORIDE	100MG	
PHR1828-500MG	EUCALYPTOL	500MG	470-82-6
PHR1866-1G	EZETIMIBE	1G	163222-33-1
PHR1482-1G	FERROUS FUMARATE	1G	141-01-5
PHR1644-1G	FERROUS LACTATE	1G	5905-52-2
PHR1483-1G	FERROUS SULFATE	1G	7782-63-0
PHR1791-100MG	FERULIC ACID	100MG	537-98-4
PHR1792-1G	FLUPHENAZINE HYDROCHLORIDE	1G	146-56-5
PHR1793-500MG	FLUVOXAMINE MALEATE	500MG	61718-82-9
PHR1816-100MG	FUROSEMIDE RELATED COMPOUND B	100MG	3086-91-7
PHR1833-1G	GALACTITOL	1G	608-66-2
PHR2042-200MG	GALACTOSAMINE HYDROCHLORIDE	200MG	1772-03-8
PHR1623-500MG	GALANTAMINE HYDROBROMIDE	500MG	1953-04-4
PHR1593-1G	GANCICLOVIR	1G	82410-32-0
PHR1859-100MG	GENISTEIN	100MG	446-72-0
PHR1567-500MG	GIBBERELIC ACID	500MG	77-06-5
PHR1748-3X1.5ML	GLACIAL ACETIC ACID	3X1.5ML	64-19-7
PHR1818-1G	GLIPIZIDE	1G	29094-61-9
PHR1484-1G	GLYCERYL MONOOLEATE 40%	1G	25496-72-4
PHR1516-200MG	GLYCYRRHIZIC ACID	200MG	1405-86-3
PHR1794-500MG	HESPERIDIN	500MG	520-26-3
PHR1782-1G	HYDROXYCHLOROQUINE SULFATE	1G	747-36-4
PHR1939-500MG	HYDROXYPROLINE	500MG	51-35-4
PHR1440-1G	HYDROXYPROPYL-BETA-CYCLODEXTRIN	1G	128446-35-5
PHR1607-1G	HYDROXYZINE DIHYDROCHLORIDE	1G	2192-20-3
PHR1829-500MG	HYPROMELLOSE	500MG	9004-65-3
PHR1795-100MG	IMINODIBENZYL	100MG	494-19-9

Cat. No.	Description	Pack Size	CAS No.
PHR1796-200MG	IMIPENEM	200MG	74431-23-5
PHR1797-1G	IMIPRAMINE HYDROCHLORIDE	1G	113-52-0
PHR1931-1G	IOHEXOL	1G	66108-95-0
PHR1769-1G	ISOMALT	1G	64519-82-0
PHR1932-500MG	ISONIAZID	500MG	54-85-3
PHR1834-200MG	ITRACONAZOLE	200MG	84625-61-6
PHR1798-250MG	LEVAMISOLE HYDROCHLORIDE	250MG	16595-80-5
PHR1850-500MG	LEVONORGESTREL	500MG	797-63-7
PHR1504-500MG	LIOTHYRONINE	500MG	6893-02-3
PHR1927-1G	LOPINAVIR	1G	192725-17-0
PHR2186-500MG	LUMEFANTRINE	500MG	82186-77-4
PHR1936-10MG	LUMIFLAVINE	10MG	
PHR1699-3X100MG	LUTEIN	3X100MG	127-40-2
PHR1770-500MG	LYCOPENE	500MG	502-65-8
PHR1770-5X100MG	LYCOPENE	5X100MG	502-65-8
PHR1363-2.5G	MAGNESIUM STEARATE	2.5G	
PHR2044-1G	MANNOSE	1G	3458-28-4
PHR1767-500MG	MELATONIN	500MG	73-31-4
PHR1955-500MG	MEPIVACAINE HYDROCHLORIDE	500MG	1722-62-9
PHR1943-1G	METHACHOLINE CHLORIDE	1G	62-51-1
PHR1594-1G	METHACYCLINE HYDROCHLORIDE	1G	3963-95-9
PHR2287-100MG	METHYL 4-(ACETYLAMINO)-2-HYDROXYBENZOATE	100MG	4093-28-1
PHR1835-500MG	METHYL NICOTINATE	500MG	93-60-7
PHR1571-1G	METHYL OLEATE	1G	112-62-9
PHR1414-1G	METHYL STEARATE	1G	112-61-8
PHR1801-500MG	MINOCYCLINE HYDROCHLORIDE	500MG	13614-98-7
PHR1490-1ML	MONOETHANOLAMINE	1ML	141-43-5
PHR1954-30MG	N- METHYL OMEPRAZOLE	30MG	
PHR1802-300MG	NALOXONE	300MG	51481-60-8
PHR1836-500MG	NAPHAZOLINE HYDROCHLORIDE	500MG	550-99-2
PHR1698-1G	NEOTAME	1G	243116-62-3
PHR1731-1G	NITISINONE	1G	104206-65-7
PHR2179-100MG	N-METHYL-2-NITROACETAMIDE	100MG	72078-82-1
PHR1714-1G	NORETHINDRONE	1G	68-22-4
PHR2119-500MG	NORGESTIMATE	500MG	35189-28-7
PHR2210-50MG	NORVALINE	50MG	760-78-1
PHR1880-10MG	OCTREOTIDE ACETATE	10MG	79517-01-4
PHR1781-500MG	OSELTAMIVIR PHOSPHATE	500MG	204255-11-8
PHR1945-300MG	OUABAIN	300MG	11018-89-6
PHR1632-500MG	PIOGLITAZONE HYDROCHLORIDE	500MG	112529-15-4
PHR1805-1G	PIPERACILLIN	1G	66258-76-2
PHR1588-200MG	PIROXICAM	200MG	36322-90-4
PHR2173-1G	PLASTIC ADDITIVE 2	1G	6683-19-8
PHR2174-1G	PLASTIC ADDITIVE 4	1G	2082-79-3
PHR1518-1G	POLYDIMETHYLSILOXANE	1G	9016-00-6
PHR2011-100MG	POTASSIUM 1-CYANOCYCLOHEXANEACETATE	100MG	133481-09-1
PHR2203-1G	POTASSIUM NITRATE	1G	7757-79-1
PHR2253-100MG	PROCAINE N-OXIDE	100MG	1027072-50-9
PHR1575-500MG	PROPANTHELINE BROMIDE	500MG	50-34-0
PHR1538-250MG	PSEUDOEPHEDRINE HYDROCHLORIDE	250MG	345-78-8
PHR1944-1G	QUINAPRIL HYDROCHLORIDE	1G	82586-55-8
PHR1750-1G	RESERPINE	1G	50-55-5
PHR2101-500MG	RIBOSE	500MG	50-69-1
PHR1946-300MG	RIFAXIMIN	300MG	80621-81-4

Cat. No.	Description	Pack Size	CAS No.
PHR1734-1G	RITONAVIR	1G	155213-67-5
PHR1889-500MG	RIZATRIPTAN BENZOATE	500MG	145202-66-0
PHR1928-1G	ROSUVASTATIN CALCIUM	1G	147098-20-2
PHR1947-500MG	SALMETEROL XINAFOATE	500MG	94749-08-3
PHR1941-200MG	SALSALATE	200MG	552-94-3
PHR1807-500MG	SILDENAFIL CITRATE	500MG	171599-83-0
PHR1715-1G	SISOMICIN SULFATE	1G	53179-09-2
PHR1857-1G	SITAGLIPTIN PHOSPHATE	1G	654671-77-9
PHR1471-1G	SODIUM ALGINATE	1G	9005-38-3
PHR1808-100MG	SODIUM AMINOSALICYLATE DIHYDRATE	100MG	6018-19-5
PHR1948-2G	SODIUM CARBONATE ANHYDROUS	2G	497-19-8
PHR1645-1G	SODIUM IRON EDTA	1G	15708-41-5
PHR1949-2G	SODIUM LAURYL SULFATE	2G	151-21-3
PHR1444-1G	SODIUM STEARYL FUMARATE	1G	4070-80-8
PHR2245-1G	SODIUM TRIPOLYPHOSPHATE HEXAHYDRATE	1G	15091-98-2
PHR1950-500MG	SPIRAMYCIN	500MG	8025-81-8
PHR1001-1G	SUCROSE	1G	57-50-1
PHR1498-1G	SULFACETAMIDE	1G	144-80-9
PHR1577-1G	SULFAMERAZINE	1G	127-79-7
PHR1809-150MG	TACROLIMUS	150MG	109581-93-3
PHR1810-500MG	TADALAFIL	500MG	171596-29-5
PHR1811-100MG	TEICOPLANIN	100MG	61036-62-2
PHR1855-500MG	TELMISARTAN	500MG	144701-48-4
PHR1592-1G	TENOFOVIR	1G	206184-49-8
PHR1957-500MG	TENOFOVIR DISOPROXIL FUMARATE	500MG	202138-50-9
PHR2027-250MG	TESTOSTERONE	250MG	58-22-0
PHR2026-500MG	TESTOSTERONE PROPIONATE	500MG	57-85-2
PHR1853-500MG	TETRAHYDROZOLINE HYDROCHLORIDE	500MG	522-48-5
PHR2205-1G	TETRASODIUM PYROPHOSPHATE DECAHYDRATE	1G	13472-36-1
PHR1837-100MG	THEOBROMINE	100MG	83-67-0
PHR1890-50MG	TIOTROPIUM BROMIDE MONOHYDRATE	50MG	411207-31-3
PHR1778-1G	TOLNAFTATE	1G	2398-96-1
PHR1891-500MG	TOLTERODINE TARTRATE	500MG	124937-52-6
PHR1600-1G	TOPIRAMATE	1G	97240-79-4
PHR2217-100MG	TOXOPYRIMIDINE	100MG	73-67-6
PHR1812-500MG	TRANEXAMIC ACID	500MG	1197-18-8
PHR2201-200MG	TRANS-RESVERATROL	200MG	501-36-0
PHR1622-3ML	TRAVOPROST	3ML	157283-68-6
PHR1495-500MG	TRIENTINE HYDROCHLORIDE	500MG	38260-01-4
PHR1749-1G	TRIMEPRAZINE TARTRATE	1G	4330-99-8
PHR1578-5ML	TROLAMINE	5ML	102-71-6
PHR1929-500MG	UBIDECARENONE	500MG	303-98-0
PHR1732-4X250MG	VANCOMYCIN HYDROCHLORIDE	4X250MG	1404-93-9
PHR1736-1G	VENLAFAXINE HYDROCHLORIDE	1G	99300-78-4
PHR1892-1G	VORICONAZOLE	1G	137234-62-9
PHR2039-500MG	XYLOMETAZOLINE HYDROCHLORIDE	500MG	1218-35-5
PHR2102-500MG	XYLOSE	500MG	58-86-6


Residual Solvents

- Residual solvents are defined as Organic Volatile Chemicals
- Used or produced in Drugs/Excipients manufacturing
- Trace amounts can impact Drug Safety & Efficacy
- Evaluation is quite crucial due to possible risk to human health
- Divided in 3 categories (Class-I, Class-II & Class-III)
- Strict guidelines to be followed by all Pharma companies



Cat. No.	Description	Pack Size	CAS No.
Residual Solvents			
PHR2108-3X1.2ML	1,2-DICHLOROETHANE	3X1.2ML	107-06-2
RTC000081-20ML	1,2-DICHLOROETHANE	20ML	107-06-2
PHR1549-3X1.2ML	1,4-DIOXANE	3X1.2ML	123-91-1
RTC000093-20ML	1,4-DIOXANE	20ML	123-91-1
PHR1209-3X1.2ML	1-BUTANOL	3X1.2ML	71-36-3
RTC000092-20ML	1-BUTANOL	20ML	71-36-3
PHR1217-3X1.2ML	1-PENTANOL	3X1.2ML	71-41-0
PHR1208-3X1.2ML	1-PROPANOL	3X1.2ML	71-23-8
PHR1563-3X1.2ML	2-BUTANOL	3X1.2ML	78-92-2
PHR1550-3X1.2ML	2-ETHOXYETHANOL	3X1.2ML	110-80-5
PHR1544-3X1.2ML	2-METHOXY ETHANOL	3X1.2ML	109-86-4
PHR2729-5X1.2ML	2-METHYL-1-PROPANOL	5X1.2ML	78-83-1
PHR1072-3X1.2ML	2-PROPANOL	3X1.2ML	67-63-0
RTC000079-20ML	2-PROPANOL	20ML	67-63-0
PHR2337-5X1.2ML	3-METHYL-1-BUTANOL	5X1.2ML	123-51-3
PHR2107-3X1.2ML	ACETALDEHYDE	3X1.2ML	75-07-0
RTC000075-20ML	ACETALDEHYDE	20ML	75-07-0
PHR1748-3X1.5ML	ACETIC ACID	3X1.5ML	64-19-7
PHR1219-3X1.2ML	ACETONE	3X1.2ML	67-64-1
RTC000078-20ML	ACETONE	20ML	67-64-1
PHR1551-3X1.2ML	ACETONITRILE	3X1.2ML	75-05-8
RTC000080-20ML	ACETONITRILE	20ML	75-05-8
PHR1212-3X1.2ML	ANISOLE	3X1.2ML	100-66-3
PHR1211-3X1.2ML	BUTYL ACETATE	3X1.2ML	123-86-4
PHR1665-3X1.2ML	CHLOROBENZENE	3X1.2ML	108-90-7
RTC000099-20ML	CHLOROBENZENE	20ML	108-90-7
PHR1552-3X1.2ML	CHLOROFORM	3X1.2ML	67-66-3
PHR1210-3X1.2ML	CUMENE	3X1.2ML	98-82-8
PHR1687-3X1.2ML	CYCLOHEXANE	3X1.2ML	110-82-7
RTC000089-20ML	CYCLOHEXANE	20ML	110-82-7
RTC000077-20ML	DIETHYL ETHER	20ML	60-29-7
PHR2183-3X1.2ML	DIISOPROPYL ETHER	3X1.2ML	108-20-3
RTC000085-20ML	DIISOPROPYL ETHER	20ML	108-20-3
RTC000103-20ML	DIMETHYL SULFOXIDE	20ML	67-68-5

Cat. No.	Description	Pack Size	CAS No.
PHR1553-3X1.2ML	DIMETHYLFORMAMIDE	3X1.2ML	68-12-2
RTC000098-20ML	DIMETHYLFORMAMIDE	20ML	68-12-2
PHR1373-5X1.2ML	ETHANOL	5X1.2ML	64-17-5
PHR1481-3X1.2ML	ETHYL ACETATE	3X1.2ML	141-78-6
RTC000087-20ML	ETHYL ACETATE	20ML	141-78-6
PHR1229-3X1.2ML	ETHYLBENZENE	3X1.2ML	100-41-4
RTC000100-20ML	ETHYLBENZENE	20ML	100-41-4
PHR1554-3X1.2ML	HEPTANE	3X1.2ML	142-82-5
RTC000091-20ML	HEPTANE	20ML	142-82-5
PHR1555-3X1.2ML	HEXANE	3X1.2ML	110-54-3
RTC000084-20ML	HEXANE	20ML	110-54-3
PHR1915-3X1.2ML	ISOOCTANE	3X1.2ML	540-84-1
PHR2109-3X1.2ML	ISOPROPYL ACETATE	3X1.2ML	108-21-4
RTC000090-20ML	ISOPROPYL ACETATE	20ML	108-21-4
PHR1547-3X1.2ML	MESITYL OXIDE	3X1.2ML	141-79-7
RTC000097-20ML	MESITYL OXIDE	20ML	141-79-7
PHR1372-3X1.5ML	METHANOL	3X1.5ML	67-56-1
RTC000076-20ML	METHANOL	20ML	67-56-1
RTC000086-20ML	METHYL ETHYL KETONE	20ML	78-93-3
PHR1556-3X1.2ML	METHYL ISOBUTYL KETONE	3X1.2ML	108-10-1
RTC000095-20ML	METHYL ISOBUTYL KETONE	20ML	108-10-1
PHR1557-3X1.2ML	METHYLENECHLORIDE	3X1.2ML	75-09-2
RTC000101-20ML	M-XYLENE	20ML	108-38-3
PHR2110-3X1.2ML	N,N-DIMETHYLACETAMIDE	3X1.2ML	127-19-5
RTC000104-20ML	N,N-DIMETHYLACETAMIDE	20ML	127-19-5
RTC000102-20ML	O-XYLENE	20ML	95-47-6
PHR1548-3X1.2ML	PROPYL ACETATE	3X1.2ML	109-60-4
RTC000094-20ML	PROPYL ACETATE	20ML	109-60-4
PHR1558-3X1.2ML	PYRIDINE	3X1.2ML	110-86-1
PHR2184-3X1.2ML	TERT-BUTANOL	3X1.2ML	75-65-0
RTC000082-20ML	TERT-BUTANOL	20ML	75-65-0
PHR1560-3X1.2ML	TERT-BUTYL METHYL ETHER	3X1.2ML	1634-04-4
RTC000083-20ML	TERT-BUTYL METHYL ETHER	20ML	1634-04-4
PHR1561-3X1.2ML	TETRAHYDROFURAN	3X1.2ML	109-99-9
RTC000088-20ML	TETRAHYDROFURAN	20ML	109-99-9
PHR1317-3X1.2ML	TOLUENE	3X1.2ML	108-88-3
RTC000096-20ML	TOLUENE	20ML	108-88-3
PHR1319-3X1.2ML	TRICHLOROETHYLENE	3X1.2ML	79-01-6



To know more contact:

Merck Life Science Private Limited

8th Floor, Godrej One, Pirojsha Nagar, Eastern Express Highway, Vikhroli (E), Mumbai – 400079

Website : www.sigmaaldrich.com