



QUALITY CONTROL RANGES

RAT/MOUSE FGF-21 ELISA CATALOG # EZRMFGF21-26K

STANDARD, LYOPHILIZED CATALOG # E8026-K

STANDARD	RAT/MOUSE FGF-21 VALUE
LOT NUMBER: 826-10A	12,644 pg/mL

The exact value of lyophilized standard is shown above.

Reconstitute with **0.5mL** distilled or deionized water.

QUALITY CONTROLS 1 AND 2, LYOPHILIZED CATALOG #E6026-K

Reconstitute each vial with **0.5mL** distilled or deionized water.

CONTROL	LOT NUMBER	EXPECTED RANGE	UNITS
QC 1	RMFGF-110A	155 - 322	pg/mL
QC 2	RMFGF-210A	1097 - 2279	pg/mL

E6026-K: RMFGF-110A.210A (26 March 2021)

EMD MILLIPORE • 1-800-645-5476 • 1-800-645-5439 (FAX) www.emdmillipore.com

The M logo is a trademark of Merck KGaA, Darmstadt, Germany.

© 2013 EMD Millipore Corporation, Billerica, MA 01821 USA.