

## 29231 Chloramphenicol Selective Supplement

Antibiotic Supplement recommended for selective isolation of yeasts and molds.

### Composition

Chloramphenicol 50 mg

(A vial is sufficient for 500 ml of medium)

### Directions

Rehydrate the contents of one vial of Cat. No. 29231 with 2 ml of ethanol and mix with 500 ml of sterile molten Dichloran Rose bengal Agar Base (Cat. No. 17147) or Aspergillus Differentiation Agar Base (Cat. No. 17121). Mix well and pour plates.

### Precautions and Disclaimer

This product is for R&D use only, not for drug, household, or other uses. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

