

3050 Spruce Street
Saint Louis, Missouri 63103 USA
Telephone 800-325-5832 • (314) 771-5765
Fax (314) 286-7828
email: techserv@sial.com
sigma-aldrich.com

ProductInformation

Tolnaftate

Product Number **T 6638** Storage Temperature 2-8 °C

Product Description

Molecular Formula: C₁₉H₁₇NOS Molecular Weight: 307.4 CAS Number: 2398-96-1 Melting Point: 110.5-111.5 °C¹

Synonyms: methyl(3-methylphenyl)carbamothioic acid *O*-2-naphthalenyl ester; m,*N*-dimethylthiocarbanilic

acid O-2-naphthyl ester1

Tolfanate is an antifungal compound which has activity against such species as *Epidermophyton*, *Microsporum*, *Trichophyton* species, and *Malassezia furfur*. ^{1,2} Its proposed mechanism of action includes the inhibition of squalene epoxidase. ³ The use of tolfanate to inhibit squalene epoxidase activity in *Candida albicans* (500 μ M) ⁴ and in *Trichophyton rubrum* (IC₅₀ = 51.5 nM) has been studied. ⁵

The activities of various antifungal compounds, including tolfanate, against dermatophytes from different species have been investigated. An *in vitro* investigation of tolfanate penetration into the human nail plate in the presence of N-acetyl-L-cysteine or 2-mercaptoethanol has been described. The effect of tolfanate on the biosynthesis of lipids in *Aspergillus niger* protoplasts has been studied.

Several LC methods for the analysis of tolfanate in formulations have been described. 9,10

Precautions and Disclaimer

For Laboratory Use Only. Not for drug, household or other uses.

Preparation Instructions

This product is soluble in chloroform (50 mg/ml), yielding a clear, colorless to very faint yellow solution. It is also soluble in acetone and sparingly soluble in ethanol and methanol.¹

References

- 1. The Merck Index, 12th ed., Entry# 9656.
- Martindale The Extra Pharmacopoeia, 31st ed., Reynolds, J. E. F., ed., Royal Pharmaceutical Society (London, UK: 1996), p. 416.
- Barrett-Bee, K., and Dixon, G., Ergosterol biosynthesis inhibition: a target for antifungal agents. Acta Biochim. Pol., 42(4), 465-479 (1995).
- Georgopapadakou, N. H., and Bertasso, A., Effects of squalene epoxidase inhibitors on Candida albicans. Antimicrob. Agents Chemother., 36(8), 1779-1781 (1992).
- Favre, B., and Ryder, N. S., Characterization of squalene epoxidase activity from the dermatophyte *Trichophyton rubrum* and its inhibition by terbinafine and other antimycotic agents. Antimicrob. Agents Chemother., 40(2), 443-447 (1996).
- Favre, B., et al., Comparison of *in vitro* activities of 17 antifungal drugs against a panel of 20 dermatophytes by using a microdilution assay. J. Clin. Microbiol., 41(10), 4817-4819 (2003).
- Kobayashi, Y., et al., Enhancing effect of N-acetyl-L-cysteine or 2-mercaptoethanol on the *in vitro* permeation of 5-fluorouracil or tolnaftate through the human nail plate. Chem. Pharm. Bull. (Tokyo), 46(11), 1797-1802 (1998).
- 8. Oh, K., et al., Effects of antimycotics on the biosynthesis of cellular macromolecules in *Aspergillus niger* protoplasts. Mycopathologia, **122(3)**, 135-141 (1993).
- 9. Thompson, R. D., and Carlson, M., Liquid chromatographic determination of tolnaftate in commercial products. J. Assoc. Off. Anal. Chem., **74(4)**, 603-607 (1991).
- 10. Dash, A. K., A liquid chromatographic method for the determination of tolnaftate in pharmaceutical formulations. J. Pharm. Biomed. Anal., 11(9),

Sigma brand products are sold through Sigma-Aldrich, Inc.

Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.

847-853 (1993).