



## MOUSE ANTI-ADENOVIRUS 1, 2, 5, 6 MONOCLONAL ANTIBODY

CATALOG NUMBER:	MAB8044
LOT NUMBER:	
QUANTITY:	100 µL
SPECIFICITY:	Adenovirus 1,2,5,6
IMMUNOGEN:	Adenovirus 1,2,5,6 hexon protein
CLONE NAME:	7/48-7a
APPLICATIONS:	IFA. Optimal working dilutions must be determined by end user.
FORMAT:	IgG1
PRESENTATION:	Ascites
STORAGE/HANDLING:	Maintain at -20°C in convenient aliquots for up to 12 months. Do not store in a diluted format.
REFERENCE:	Sharp, I.R.: (1996) <i>Rapid Detection and Serotyping of Adenovirus by Direct Immunofluorescence</i> Central Public Health Laboratory, London UK. Wood et al., (1997) <i>J. Med. Virology</i> <b>51</b> :198-201.

**Important Note:** During shipment, small volumes of product will occasionally become entrapped in the seal of the product vial. For products with volumes of 200 µL or less, we recommend gently tapping the vial on a hard surface or briefly centrifuging the vial in a tabletop centrifuge to dislodge any liquid in the container's cap.

FOR RESEARCH USE ONLY; NOT FOR USE IN DIAGNOSTIC  
PROCEDURES. NOT FOR HUMAN OR ANIMAL CONSUMPTION

Unless otherwise stated in our catalog or other company documentation accompanying the product(s), our products are intended for research use only and are not to be used for any other purpose, which includes but is not limited to, unauthorized commercial uses, in vitro diagnostic uses, ex vivo or in vivo therapeutic uses or any type of consumption or application to humans or animals.

©2002 - 2009: Millipore Corporation. All rights reserved. No part of these works may be reproduced in any form without permission in writing.

28820 Single Oak Drive • Temecula, CA 92590  
Technical Support: T: 1-800-MILLIPORE (1-800-645-5476) • F: 1-800-437-7502  
www.millipore.com