

55181 Brilliant Green + Novobiocin Selective Supplement

A supplement used for the selective enrichment of Salmonella.

Composition

(per vial, sufficient for 500 ml medium)

Brilliant Green	5 mg
Novobiocin, sodium salt	20 mg
Distilled water	4 ml

Directions

Mix the liquid with the powder by pressing down on the cap. Shake till total dissolution and aseptically add the vial contents to 500 ml of boiled Muller-Kauffmann Tetrathionate Broth base (Cat. No. 43052) cooled to 50°C.

Distribute into sterile tubes and dispense 20 ml in each.

Basal medium with supplement, kept in the refrigerator, it is effective for 2 months. Once iodine solution is added it is only effective for 48 hours.

Note: Don't heat the media once the supplements have been added.

Precautions and Disclaimer

This product is for R&D use only, not for drug, household, or other uses. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

