



114420


MQuant®
Nickel Test
Nickel-Test
Test Nickel
Test Níquel

Ni


Color card
Farbkarte
Carte colorimétrique
Tarjeta colorimétrica



1.



2.




7.69144.2005-XXXXXX 05/2019

Merck KGaA, 64271 Darmstadt, Germany,
Tel. +49(0)6151 72-2440
www.analytical-test-kits.com


EMD Millipore Corporation, 400 Summit Drive
Burlington MA 01803, USA, Tel. +1-978-715-4321

Sigma-Aldrich Canada Co. or Millipore (Canada) Ltd.
2149 Winston Park, Dr. Oakville, Ontario, L6H 6J8
Phone: +1 800-565-1400


0 0.02 0.04 0.07 0.10 0.15 0.2 0.3 0.4 0.5 mg/l Ni



20 ml



4 x Ni-1



+ Ni-1 ?



1 min



8 x Ni-2



pH 10 - 12



8 x Ni-3



3 min

