

Product Information

HORSE SERUM

Product Number **H 0146**

Product Description

Horse Serum is a pooled preparation obtained from a normal donor herd. Normal horse serum is often used as a blocking agent or negative control in immunoassays.

Reagent

This product is provided as a liquid containing preservative.

Precautions and Disclaimer

For R&D use only. Not for drug, household, or other uses. Consult the MSDS for information regarding hazards and safe handling practices.

Storage/Stability

Store at 2-8 °C for up to one month.

For extended storage, the solution may be frozen in working aliquots. Repeated freezing and thawing is not recommended. Storage in "frost-free" freezers is not recommended.

MAM,DAA,SM 08/20-1