



CHEMICON International, Inc. 28835 Single Oak Drive • Temecula, CA 92590

Phone: 909-676-8080 • 800-437-7500 • Fax: 909-676-9209

custserv@chemicon.com • techserv@chemicon.com • www.chemicon.com

MOUSE ANTI-HUMAN P SELECTIN MONOCLONAL ANTIBODY

CATALOG NUMBER:	MAB2154
LOT NUMBER:	21041929
QUANTITY:	100 µg
CONCENTRATION:	1 mg/mL
SPECIFICITY:	Human P Selectin
IMMUNOGEN:	IL-1β activated HUVEC (human)
ISOTYPE:	IgG ₁
CLONE:	P8G6
APPLICATIONS:	Immunohistochemistry: for use on acetone or paraformaldehyde fixed tissue Immunocytochemistry Immunoprecipitation ELISA Inhibits P Selectin dependent interactions Inhibits platelet binding to monocytes Optimal working dilutions must be determined by end user.
FORMAT:	Purified in 0.02M PBS, 0.25M NaCl, pH 7.6 with 0.1% NaN ₃
PRESENTATION:	Liquid
STORAGE/HANDLING:	Maintain at 2-8°C in undiluted aliquots for up to 6 months.

For research use only; not for use as a diagnostic.

Important Note: *During shipment, small volumes of antibody will occasionally become entrapped in the seal of the product vial. For antibodies with volumes of 200 µl or less, we recommend gently tapping the vial on a hard surface or briefly centrifuging the vial in a tabletop centrifuge to dislodge any liquid in the container's cap.*

Have you published a paper using a Chemicon primary antibody? You may be eligible for credit toward future purchases of Chemicon products! Contact our Technical Service department for more information.