



MOUSE ANTI-ENTEROVIRUS 71 MONOCLONAL ANTIBODY

CATALOG NUMBER:	MAB979
LOT NUMBER:	
QUANTITY:	100 µL
SPECIFICITY:	Reacts with Enterovirus 71. Neutralizes enterovirus 71 BrCr strain at a titer of less than 1:14.
ISOTYPE:	IgG ₁
APPLICATIONS:	Indirect immunofluorescence at 1:1,000. (Also cross-reacts with Coxsackie A16). Optimal working dilutions must be determined by end user.
FORMAT:	Ascites.
STORAGE/HANDLING:	Maintain at -20°C in undiluted aliquots for up to 12 months. Avoid repeated freeze/thaw cycles.
REFERENCES:	<i>Pediatr. Infect. Dis. J.</i> (1995) 14 :1095

Important Note: During shipment, small volumes of product will occasionally become entrapped in the seal of the product vial. For products with volumes of 200 µL or less, we recommend gently tapping the vial on a hard surface or briefly centrifuging the vial in a tabletop centrifuge to dislodge any liquid in the container's cap.

For research use only; not for use as a diagnostic.

Unless otherwise stated in our catalog or other company documentation accompanying the product(s), our products are intended for research use only and are not to be used for any other purpose, which includes but is not limited to, unauthorized commercial uses, in vitro diagnostic uses, ex vivo or in vivo therapeutic uses or any type of consumption or application to humans or animals.

©2007: Millipore Corporation. All rights reserved. No part of these works may be reproduced in any form without permission in writing.

USA & Canada • Phone: +1(800) 437-7500 • Fax: +1 (951) 676-9209
www.millipore.com

042605/mab979/ef