



## SIGMA QUALITY CONTROL TEST PROCEDURE

## Product Information

### Enzymatic Assay of $\alpha$ -1-ANTITRYPSIN Inhibition of Trypsin Activity

#### PRINCIPLE:



Abbreviation used:

BAEE = N $\alpha$ -Benzoyl-L-Arginine Ethyl Ester

$\alpha$ -1-Antitrypsin inhibits this reaction.

**CONDITIONS:** T = 25°C, pH = 7.6, A<sub>253nm</sub>, Light path = 1 cm

**METHOD:** Continuous Spectrophotometric Rate Determination

#### REAGENTS:

- A. 67 mM Sodium Phosphate Buffer, pH 7.6 at 25°C  
(Prepare 100 ml in deionized water using Sodium Phosphate, Dibasic, Anhydrous, Sigma Prod. No. S-0876. Adjust to pH 7.6 at 25°C with 1 M HCl.)
- B. 0.25 mM N $\alpha$ -Benzoyl-L-Arginine Ethyl Ester Solution (BAEE)  
(Prepare 50 ml in Reagent A using N $\alpha$ -Benzoyl-L-Arginine Ethyl Ester, Hydrochloride, Sigma Prod. No. B-4500.)
- C. 1 mM Hydrochloric Acid Solution (HCl)  
(Prepare 50 ml in deionized water using concentrated Hydrochloric Acid, Sigma Prod. No. H-7020.)
- D. Trypsin Enzyme Solution (Trypsin)  
(Immediately before use, prepare a solution containing 0.5 mg/ml of Trypsin, Sigma Prod. No. T-8003, in cold Reagent C.)
- E.  $\alpha$ -1-Antitrypsin Solution ( $\alpha$ -1-Antitryp)  
(Immediately before use, prepare a solution containing 2 mg solid/ml of  $\alpha$ -1-Antitrypsin in cold Reagent A for A-9024. A-6150 requires a concentration of 4 mgs solid/ml.)

**Enzymatic Assay of  $\alpha$ -1-ANTITRYPSIN  
Inhibition of Trypsin Activity**

**PROCEDURE:**

Step 1:

Pipette (in milliliters) the following reagents into suitable containers:

	<u>Uninh</u>	<u>Test 1</u>	<u>Test 2</u>	<u>Test 3</u>	<u>Test 4</u>	<u>Test 5</u>
Reagent E ( $\alpha$ -1-Antitryp)	---	0.10	0.20	0.30	0.40	0.50
Reagent D (Trypsin)	0.50	0.50	0.50	0.50	0.50	0.50
Reagent C (HCl)	9.50	9.40	9.30	9.20	9.10	9.00

Mix by inversion.

Step 2:

Pipette (in milliliters) the following reagents into suitable containers:

	<u>Uninh</u>	<u>Test 1</u>	<u>Test 2</u>	<u>Test 3</u>	<u>Test 4</u>	<u>Test 5</u>	<u>Blank</u>
Reagent B (BAEE)	3.00	3.00	3.00	3.00	3.00	3.00	3.00
Reagent C (HCl)	---	---	---	---	---	---	0.20

Mix by inversion and equilibrate to 25°C. Monitor the  $A_{253nm}$  until constant, using a suitably thermostatted spectrophotometer. Then add:

	<u>Uninh</u>	<u>Test 1</u>	<u>Test 2</u>	<u>Test 3</u>	<u>Test 4</u>	<u>Test 5</u>	<u>Blank</u>
Uninh (Step 1)	0.20	---	---	---	---	---	---
Test 1 (Step 1)	---	0.20	---	---	---	---	---
Test 2 (Step 1)	---	---	0.20	---	---	---	---
Test 3 (Step 1)	---	---	---	0.20	---	---	---
Test 4 (Step 1)	---	---	---	---	0.20	---	---
Test 5 (Step 1)	---	---	---	---	---	0.20	---

Immediately mix by inversion and record the increase in  $A_{253nm}$  for approximately 5 minutes. Obtain the  $\Delta A_{253nm}/\text{minute}$  using the maximum linear rate for the Tests, Blank, and Uninhibited Solution.

## Enzymatic Assay of $\alpha$ -1-ANTITRYPSIN Inhibition of Trypsin Activity

### CALCULATIONS:

Trypsin Activity:

$$\text{BAEE units/mg enzyme} = \frac{(A_{253\text{nm}}/\text{min Test} - A_{253\text{nm}}/\text{min Blank})(\text{df})}{(0.001)(\text{mg Trypsin/RM})}$$

df = Dilution factor

0.001 = The change in  $A_{253\text{nm}}$ /minute per unit of Trypsin as per the Unit Definition (One BAEE unit =  $\Delta A_{253}$  of 0.001 per minute with BAEE as substrate at pH 7.6 at 25°C. Reaction volume = 3.2 ml (1 cm light path))

RM = Reaction Mixture

Plot the Trypsin activity versus ml of  $\alpha$ -1-Antitrypsin. Determine from the plot the amount of  $\alpha$ -1-Antitrypsin needed for 100% inhibition.

The amount (in mg) of  $\alpha$ -1-Antitrypsin needed for 100% inhibition = (ml of  $\alpha$ -1-Antitrypsin for 100% inhibition)(conc. of  $\alpha$ -1-Antitrypsin (mg/ml))

$$\text{mg } \alpha\text{-1-Antitrypsin to inhibit 1 mg Trypsin} = \frac{\text{mg of } \alpha\text{-1-Antitrypsin needed for 100\% Inhibition}}{\text{mg Trypsin/Reaction Mix}}$$

### FINAL ASSAY CONCENTRATIONS:

In a 3.20 reaction mix, the final concentrations are 0.23 mM N $\alpha$ -benzoyl-L-arginine ethyl ester, 63 mM sodium phosphate, 0.06 M hydrochloric acid, 0.005 mg trypsin, 0.004 - 0.02 mg  $\alpha$ -1-antitrypsin.

### REFERENCE:

Bergmeyer, H.U., Gawehn, K. and Grassl, M. (1974) in *Methods of Enzymatic Analysis* (Bergmeyer, H.U. ed.) Volume I, 2nd ed., 516-517, Academic Press, New York, NY

Crawford, I.P. (1973) *Archives of Biochemistry and Biophysics* **156**, 215-222

Eriksson, S. (1965) *Acta Medica Scandinavica*, Suppementum 432, 6-12

**Enzymatic Assay of  $\alpha$ -1-ANTITRYPSIN  
Inhibition of Trypsin Activity**

**NOTES:**

1. This assay is based on the cited references.
2. Where Sigma Product or Stock numbers are specified, equivalent reagents may be substituted.

**Sigma warrants that the above procedure information is currently utilized at Sigma and that Sigma products conform to the information in Sigma publications. Purchaser must determine the suitability of the information and products for its particular use. Upon purchase of Sigma products, see reverse side of invoice or packing slip for additional terms and conditions of sale.**