



QUALITY CONTROL RANGES
MILLIPIXEL[®] MAP Human Apolipoprotein Magnetic Bead Panel Kit
Catalog # AP0MAG-62K

Standard Catalog # HAP0-8062
Lot# 862M-12K

Apolipoproteins	Standard Value		Apolipoproteins	Standard Value
Apo AI	48089 ng/mL		Apo AII	4066 ng/mL
Apo CIII	3878 ng/mL		Apo B	37064 ng/mL
Apo E	1825 ng/mL		Apo CII	1817 ng/mL

Control Catalog # HAP0-6062
Lot# AP0M-112 and AP0M-212

<u>Apolipoproteins</u>	<u>QC Level</u>	<u>Lot #</u>	<u>Expected Range</u>	<u>Units</u>
Apo-AI	Control I	AP0M-112	216 - 450	ng/mL
	Control II	AP0M-212	1075 - 2233	ng/mL
Apo-CIII	Control I	AP0M-112	18 - 37	ng/mL
	Control II	AP0M-212	82 - 170	ng/mL
Apo-E	Control I	AP0M-112	3.9 - 8.3	ng/mL
	Control II	AP0M-212	22 - 46	ng/mL
Apo-AII	Control I	AP0M-112	10 - 21	ng/mL
	Control II	AP0M-212	66 - 137	ng/mL
Apo-B	Control I	AP0M-112	113 - 234	ng/mL
	Control II	AP0M-212	613 - 1274	ng/mL
Apo-CII	Control I	AP0M-112	6.3 - 13.1	ng/mL
	Control II	AP0M-212	33 - 69	ng/mL

HAP0-6062: AP0M-112.212 (19 JUNE 12025)

SigmaAldrich.com/support • SigmaAldrich.com

© 2025 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved.

The vibrant M logo and MILLIPIXEL are trademarks of Merck KGaA, Darmstadt Germany or its affiliates.