

QUALITY CONTROL RANGES

**HUMAN S100B ELISA
CATALOG # EZHS100B-33K**

**STANDARD, LYOPHILIZED
CATALOG # E8033-K**

STANDARD	HUMAN S100B VALUE
LOT NUMBER: 833-16K	1891 pg/mL

The exact value of lyophilized standard is shown above.

Reconstitute with **0.5mL** distilled or deionized water.

QUALITY CONTROLS 1 AND 2, LYOPHILIZED CATALOG #E6033-K

Reconstitute each vial with **0.5mL** distilled or deionized water.

CONTROL	LOT NUMBER	EXPECTED RANGE	UNITS
QC 1	S100B-116	17 - 35	pg/mL
QC 2	S100B-216	182 - 378	pg/mL

E6033-K: S100B-116.216 (07 OCTOBER 2025)

• SigmaAldrich.com/support • SigmaAldrich.com

© 2025 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved.

The vibrant M logo and MILLIPLEX are trademarks of Merck KGaA, Darmstadt
Germany or its affiliates