



MOUSE ANTI-PARAINFLUENZA III MONOCLONAL ANTIBODY

CATALOG NUMBER:	MAB855-1
LOT NUMBER:	
QUANTITY:	100 µg
CONCENTRATION:	1 mg/mL
SPECIFICITY:	Specific for the HN protein, common to all Parainfluenza III strains. Exhibits no cross-reactivity to Parainfluenza I or II.
IMMUNOGEN:	C243
ISOTYPE:	IgG _{2b}
CLONE NAME:	240/12D
APPLICATIONS:	ELISA Indirect immunofluorescence (fresh frozen tissue sections). Optimal working dilutions must be determined by end user.
FORMAT:	Liquid purified immunoglobulin.
STORAGE/HANDLING:	Maintain at 2-8C for up to 6 months after date of receipt.

Important Note: *During shipment, small volumes of product will occasionally become entrapped in the seal of the product vial. For products with volumes of 200 µL or less, we recommend gently tapping the vial on a hard surface or briefly centrifuging the vial in a tabletop centrifuge to dislodge any liquid in the container's cap.*

**FOR RESEARCH USE ONLY; NOT FOR USE IN DIAGNOSTIC
PROCEDURES. NOT FOR HUMAN OR ANIMAL CONSUMPTION**

Unless otherwise stated in our catalog or other company documentation accompanying the product(s), our products are intended for research use only and are not to be used for any other purpose, which includes but is not limited to, unauthorized commercial uses, in vitro diagnostic uses, ex vivo or in vivo therapeutic uses or any type of consumption or application to humans or animals.

©2002 - 2011: Millipore Corporation. All rights reserved. No part of these works may be reproduced in any form without permission in writing.