

P9477 Perfringens S.F.P. Supplement

Recommended for the selective isolation of *Clostridium perfringens*.

Composition:

(per vial, sufficient for 500 ml medium)

Kanamycin sulphate 6mg

Polymyxin B 15'000IU

Store at 4°C.

Directions:

Rehydrate the contents of the vial aseptically with 2 ml of sterile filtered distilled water. Mix well and aseptically add it to 475 ml of sterile, molten Perfringens Agar Base (Fluka 39727) or TSC Agar (Fluka 93745), along with 25 ml of sterile Egg Yolk Emulsion (17148). Mix well and pour into sterile petri plates.