



SIGMA-ALDRICH

3050 Spruce Street  
Saint Louis, Missouri 63103 USA  
Telephone 800-325-5832 • (314) 771-5765  
Fax (314) 286-7828  
email: techserv@sial.com  
sigma-aldrich.com

## Product Information

### Magnesium sulfate heptahydrate ACS Reagent

Product Number **M 9397**  
Store at Room Temperature  
**23,039-1 is an exact replacement for M 9397**

#### Product Description

Molecular Formula:  $\text{MgSO}_4 \cdot 7\text{H}_2\text{O}$   
Molecular Weight: 246.5  
CAS Number: 10034-99-8

This product is designated as ACS Reagent grade and meets the specifications of the American Chemical Society (ACS) for reagent chemicals.

Magnesium sulfate is used in chemistry and molecular biology as a source of magnesium ions. Magnesium has a variety of biological roles in enzymology, cell membrane and wall structural integrity, muscle cell physiology, and nucleic acid structure.<sup>1,2</sup> Magnesium is an essential co-factor in many enzymes, including deoxyribonuclease (DNase), the restriction enzymes *EcoR* I and *EcoR* V, and Ribonuclease H.<sup>3,4</sup> Magnesium also stabilizes polymeric nucleic acids such as transfer RNA and ribozymes.<sup>5</sup>

A protocol for chromosome and nuclei isolation that incorporates  $\text{MgSO}_4$  has been published.<sup>6</sup> The use of  $\text{MgSO}_4$  for leukotoxin production from cultured *Pasteurella haemolytica* has been reported.<sup>7</sup>

#### Precautions and Disclaimer

For Laboratory Use Only. Not for drug, household or other uses.

#### Preparation Instructions

This product is soluble in water (100 mg/ml), yielding a clear, colorless solution.

#### References

1. Cowan, J. A., in *The Biological Chemistry of Magnesium*, Cowan, J. A., ed., VCH Publishers (New York: 1995), pp. 1-23.
2. *The Biological Chemistry of the Elements*, Frausto da Silva, J. J. R., and Williams, R. J. P., Clarendon Press (Oxford, UK: 1991), pp. 243-267.
3. Brooks, J. E., Properties and uses of restriction endonucleases. *Methods Enzymol.*, **152**, 113-129 (1987).
4. Black, C. B., and Cowan, J. A., in *The Biological Chemistry of Magnesium*, Cowan, J. A., ed., VCH Publishers (New York: 1995), pp. 137-157.
5. *Principles of Bioinorganic Chemistry*, Lippard, S. J., and Berg, J. M., University Science Books (Mill Valley, CA: 1994), pp. 192-196.
6. Trask, B., and van den Engh, G., Chromosome and nuclei isolation with the  $\text{MgSO}_4$  procedure. *Methods Cell. Biol.*, **33**, 363-367 (1990).
7. Sun, Y., and Clinkenbeard, K.D., Serum-free culture of *Pasteurella haemolytica* optimized for leukotoxin production. *Am. J. Vet. Res.*, **59(7)**, 851-855 (1998).

GCY/NSB 3/03

Sigma brand products are sold through Sigma-Aldrich, Inc.

Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.