

**RABBIT ANTI-Kv1.2
AFFINITY PURIFIED
POLYCLONAL ANTIBODY**

CATALOG NUMBER: AB5924

LOT NUMBER: 22040933

QUANTITY: 100 µg

CONCENTRATION 1.3 mg/mL

SPECIFICITY: Kv1.2

IMMUNOGEN: Synthetic peptide from the C-terminus of rat Kv1.2.

APPLICATIONS: Immunoblotting
Immunohistochemistry on 4% paraformaldehyde fixed tissue.
Immunocytochemistry
Optimal working dilutions must be determined by end user.

SPECIES REACTIVITIES: Human, rat and mouse. Other species have not yet been tested.

FORMAT: Affinity purified immunoglobulin.

PRESENTATION: Liquid. Buffer = 0.1M Phosphate buffer saline, pH 7.4. Contains no preservative.

STORAGE/HANDLING: Maintain at -20°C in undiluted aliquots for up to 6 months after date of receipt. Avoid repeated freeze/thaw cycles.

REFERENCE: Malin, S. and Nerbonne, J., *Journal of Neuroscience* (2001) **21**:8004-8014.

For research use only; not for use as a diagnostic.

Important Note: *During shipment, small volumes of product will occasionally become entrapped in the seal of the product vial. For products with volumes of 200 µL or less, we recommend gently tapping the vial on a hard surface or briefly centrifuging the vial in a tabletop centrifuge to dislodge any liquid in the container's cap.*