

Technical Data Sheet

Kanamycin Esculin Azide Agar

GranuCult® plus, For the isolation, differentiation, and enumeration of Enterococci in foodstuffs, water, and other materials according to MOSSEL *et al.* (1978).

1052220500

Kanamycin esculin azide agar is, unlike culture media containing bile which sometimes exhibits a fluctuating selectivity towards D-Streptococci, always highly selective for this group of bacteria.

Mode of Action

Kanamycin and azide largely inhibit the accompanying bacterial flora. D-Streptococci are, however, only slightly sensitive to these substances, so they can grow almost normally and hydrolyze the glucoside esculin to give glucose and esculetin. Esculetin forms an olive-green to black complex with iron (III) ions.

Typical Composition (g/L)

Kanamycin Esculin Azide Agar	
Peptones from casein	20.0
Yeast extract	5.0
Sodium chloride	5.0
Sodium citrate	1.0
Sodium azide	0.15
Kanamycin sulfate	002
Esculin	1.0
Ammonium iron (III) citrate	0.5
Agar-agar**	15.0

**Agar-agar is equivalent to other different terms of agar.

Preparation

Suspend 47.5 g/liter, autoclave (15 minutes at 121 °C), and pour plates.

Do not overheat.

pH: 7.1 ± 0.2 at 25 °C.

The plates are clear and brown-bluish.

Experimental Procedure and Evaluation

Inoculate by spreading the samples on the surface of the plates.

Incubation: up to 3 days at 35 °C or 42 °C aerobically. The higher temperature increases the selectivity of the medium.

Enterococci colonies are surrounded by a dark zone. Confirmatory tests for example catalase test, glucose utilization, and growth at 45 ± 1 °C, may be carried out.

Quality Control

Control strains	Recovery	Color change to olive green-black
<i>Enterococcus faecalis</i> ATCC 11700	≥ 70 %	+
<i>Enterococcus hirae</i> ATCC 8043 (WDCM 00089)	≥ 70 %	+
<i>Enterococcus durans</i> ATCC 6056	≥ 70 %	+
<i>Staphylococcus aureus</i> ATCC 6538-P (WDCM 00033)	No limit	-
<i>Bacillus cereus</i> ATCC 11778 (WDCM 00001)	≤ 0.01 %	-
<i>Escherichia coli</i> ATCC 11775 (WDCM 00090)	≤ 0.01 %	-

Please refer to the actual batch-related Certificate of Analysis.

A recovery rate of 50 % is equivalent to a productivity value of 0.5.



Enterococcus faecalis ATCC 29212



Streptococcus pyogenes ATCC 19615

Literature

BRANDL, E., ASPERGER, H., PFLEGER, F., u. IBEN, CH.: Zum Vorkommen von D-Streptokokken in Käse. - Arch. Lebensmittelhyg., 36; 18-22 (1985).

MOSSEL, D.A.A., BIJKER, P.G.H., a. EELDERING, J.: Streptokokken der Lancefield-Gruppe D in Lebensmitteln und Trinkwasser - Ihre Bedeutung, Erfassung und Bekämpfung. - Arch. f. Lebensmittelhyg., 29; 121-127 (1978).

Ordering Information

Product	Cat. No.	Pack size
Kanamycin Esculin Azide Agar	1052220500	500 g

Notice

We provide information and advice to our customers on application technologies and regulatory matters to the best of our knowledge and ability, but without obligation or liability. Existing laws and regulations are to be observed in all cases by our customers. This also applies in respect to any rights of third parties. Our information and advice do not relieve our customers of their own responsibility for checking the suitability of our products for the envisaged purpose.

The information in this document is subject to change without notice and should not be construed as a commitment by the manufacturing or selling entity, or an affiliate. We assume no responsibility for any errors that may appear in this document.

Technical Assistance

Visit the tech service page at SigmaAldrich.com/techservice.

Terms and Conditions of Sale

Warranty, use restrictions, and other conditions of sale may be found at SigmaAldrich.com/terms.

+Contact Information

For the location of the office nearest you, go to SigmaAldrich.com/offices.

The life science business of Merck KGaA, Darmstadt, Germany operates as MilliporeSigma in the U.S. and Canada.

Lit No. MS_DS13804EN

MilliporeSigma, Millipore, GranuCult and Sigma-Aldrich are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.
© 2024 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved.

Document 105222dat Rev 04/24

**MILLIPORE
SIGMA**